



Sep 05, 2025

Chemical transformation of plasmid into E. coli competent cells

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqp675vk5/v1>

Ainsley Lederer¹

¹University of Pittsburgh



Ainsley Lederer

University of Pittsburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqp675vk5/v1>

Protocol Citation: Ainsley Lederer 2025. Chemical transformation of plasmid into E. coli competent cells. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wgqp675vk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 04, 2025



Last Modified: September 05, 2025

Protocol Integer ID: 226466

Keywords: chemical transformation of plasmid, plasmid, supercoiled dna, chemical transformation, method for the transformation, heat shock method, dna

Abstract

From hundreds of transformations, I have optimized this method for the transformation of supercoiled DNA into *E. coli* using the heat shock method. This is meant as a general guideline, but it is what I have found to work best for my lab.

Guidelines

This is meant to be a general protocol for cells prepared using a modified version of the Hanahan method.

Materials

1. *E. coli* competent cells (50 uL) - found in -80°C
2. Plasmid for transformation - found in -20°C
3. SOC media - found in -20°C
4. 42°C water bath and 37°C shaker
5. LB agar plate with antibiotic

Before start

Understand that the total time of this procedure is 1+ day due to the overnight incubation step. Plan the experiment ahead of time!



- 1 Thaw the competent cells and plasmid on ice for ~10 minutes.
- 2 Add 10-100 ng of plasmid (typically 1-5 μ L) to the competent cells. *Gently* flick the tube to mix the cells and DNA. DO NOT vortex or pipette to mix.
- 3 Allow the cells to sit on ice for 30 minutes.
- 3.1 A few minutes before the next step, take an aliquot of SOC media and place in the 37°C water bath.
- 4 After 30 minutes, put cells in foam rack and place in 42°C water bath for EXACTLY 45 seconds.
Do not shake the cells in the bath, and do not exceed 45 seconds. This step is crucial.
- 5 Transfer the cells back to ice and allow them to sit for 2 minutes. After 2 minutes, open the caps of the tubes and add 250 μ L pre-warmed SOC media to the cells. Allow the tubes to sit with caps open for about a minute.
- 6 Tightly close the lids to the tubes. Take some tape, and tape the tubes horizontally to the 37°C shaker incubator. Allow them to shake at 37°C for 1 hour. This allows for the cells to recover from heat shock.
- 6.1 During this time, take your LB agar plate(s) and place upside down in the 37°C incubator.
- 7 After the incubation, bring the pre-warmed plates and cells back to the bench. If the plates have excess moisture, use an aspirator to remove the liquid by tilting plate to the side. DO NOT touch the agar with the aspirator.
- 8 Pipette the required volume of cells on the LB agar plate with the correct antibiotics for your cell strain. Use a sterile glass spreader to uniformly distribute the cells.
- 8.1 The required volume of cells is highly variable depending on cell strain, plasmid type, and plasmid amount. As a general guideline, I recommend the following:
 - BL21 (DE3) cells = 80-100 μ L
 - Rosetta or Rosetta2 (DE3) cells = 150 μ L
 - BL21 (DE3) with pBirA = 100 μ L
 - DH5 α = 20 μ L or just perform a streak using inoculation loop. These cells tend to grow a lawn.



- 9 Allow the plate to sit at room temperature for about 10 minutes in order to absorb the liquid. Then, place plates upside down in the 37°C incubator and leave overnight.