

Apr 30, 2018

Version 1

## Chemically competent *V. natriegens* cells V.1

DOI

[dx.doi.org/10.17504/protocols.io.pskdncw](https://dx.doi.org/10.17504/protocols.io.pskdncw)

Carlos H<sup>1</sup>

<sup>1</sup>iGEM Marburg 2018



Carlos H

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.pskdncw>

**Protocol Citation:** Carlos H 2018. Chemically competent *V. natriegens* cells. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.pskdncw>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

We are still developing and optimizing this protocol

**Created:** April 26, 2018



**Last Modified:** April 30, 2018

**Protocol Integer ID:** 11820

**Keywords:** vibrio natriegens cell, vibrio natriegen, cell

## Abstract

This protocol describes how to make chemically competent *Vibrio natriegens* cells. The protocol was described and published by Weinstock et al., 2016

## Guidelines

All steps are done at room temperature (RT).  
This protocol was published by Weinstock et al., 1016

## Materials

### MATERIALS

- ⊗ PIPES **P212121**
- ⊗ Potassium chloride **P212121**
- ⊗ Sodium Chloride **Fisher Scientific Catalog #S271**
- ⊗ Magnesium Chloride **Fisher Scientific Catalog #AC223210010**
- ⊗ Manganese chloride **Fisher Scientific Catalog #7773-01-5**
- ⊗ brain Heart Infusion Broth **Catalog #Oxoid CM1135-UK**
- ⊗ Calcium chloride, dihydrate **Bio Basic Inc. Catalog #CD0050.SIZE.500g**

## Troubleshooting



## Before start

Make sure you have all solutions there:

- **Brain heart infusion + V<sub>2</sub>-Salts**
- **V<sub>2</sub>-Salts**
- NaCl 204mM (58.44g/mol) → 11.92176g for 200ml
- 
- KCl 4,2mM (74.55g/mol) → 0,31311g for 200ml
- 
- MgCl<sub>2</sub> 23,14mM (203.3g/mol) → 7.704362g for 200ml
- 
- **MgCl<sub>2</sub> x 2H<sub>2</sub>O** 100mM (203.3g/mol) → 0.10165g for 200ml
- **CaCl<sub>2</sub> x 2H<sub>2</sub>O** 100mM (147.02g/mol) → 2.9404g for 200ml
- **Modified Inoue buffer**
- MnCl<sub>2</sub> x 2H<sub>2</sub>O 55mM (161.87g/mol) → 1.78057g for 200ml
- 
- CaCl<sub>2</sub> x 2H<sub>2</sub>O 15mM (147.02g/mol) → 0.44106g for 200ml
- 
- KCl 250mM (75.55g/mol) → 3.7275g for 200ml
- 
- PIPES 10mM (302.37g/mol) → 4ml (from 0.5M stock solution) for 200ml
- 
- **PIPES** 500mM (302.37g/mol) → 7.55925g in 50ml (adjust pH 6.7)
- **DMSO**



- 1 Inoculate 150ml BHI + V2-Salts in a baffled flask
- 2 Incubate shaking: OD = 0.4, 30°C, 200rpm  



🌡 30 °C

🕒 02:30:00 Competent *V. natriegens*
- 3 Split into three 50ml falcons
- 4 Centrifuge: 3000g, 5min, RT  

🕒 00:05:00 1. Centrifuge competent *V. natriegens*
- 5 Remove supernatant completely
- 6 Resuspend by gently inversion in 5ml 100mM MgCl<sub>2</sub>
- 7 Pool the cells in two 50ml falcons
- 8 Fill up to 30ml with 100mM MgCl<sub>2</sub>
- 9 Centrifuge: 3000g, 4min, RT  

🕒 00:04:00 2. Centrifuge competent *V. natriegens*
- 10 Remove supernatant completely
- 11 Resuspend the pellet by gently inversion in 5ml 100mM CaCl<sub>2</sub>
- 12 Pool the cells into one 50ml falcon
- 13 Fill up to 30ml with 100mM CaCl<sub>2</sub>



- 14 Incubate: 20min, RT  
 00:20:00 Incubation competent *V. natriegens*
- 15 Centrifuge: 3000g, 4min, RT  
 00:04:00 3. Centrifuge competent *V. natriegens*
- 16 Remove supernatant completely
- 17 Resuspend the pellet by gently inversion in 1.5ml modified Inoue buffer
- 18 Add DMSO to a volume concentration of 7% (=105µl)
- 19 Aliquot the cells into chilled tubes (50µl aliquots)
- 20 Freeze at -80°C