

Aug 17, 2020

# Chicken immunization with Keyhole limpet hemocynin (KLH)-gp120 fragment (254-274) conjugate raises anti-KLH antibodies in egg yolks.

DOI

[dx.doi.org/10.17504/protocols.io.bjvtkn6n](https://dx.doi.org/10.17504/protocols.io.bjvtkn6n)

Angel A Justiz-Vaillant<sup>1</sup>

<sup>1</sup>University of the West Indies St. Augustine

University of the West In...

angel.vaillant@sta.uwi.e...



Angel A Justiz-Vaillant

University of the West Indies St. Augustine

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bjvtkn6n](https://dx.doi.org/10.17504/protocols.io.bjvtkn6n)

**Protocol Citation:** Angel A Justiz-Vaillant 2020. Chicken immunization with Keyhole limpet hemocynin (KLH)-gp120 fragment (254-274) conjugate raises anti-KLH antibodies in egg yolks.. [protocols.io https://dx.doi.org/10.17504/protocols.io.bjvtkn6n](https://dx.doi.org/10.17504/protocols.io.bjvtkn6n)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 17, 2020

**Last Modified:** August 17, 2020

**Protocol Integer ID:** 40595



## Abstract

Chicken immunization with peptides is ineffective if only just the peptides are being inoculated. However, to make the immune response effective the fragment 254-274 of HIV-1 was conjugated with a carrier protein (KLH) that produced a critical immune response, assessed by ELISA, Immunoblot analysis and dot blot [1-4]. The Polson method (1990) can be used effectively to separate the IgY antibody from the egg yolk of immunized chickens [5].

## References

1. Schade R, Pfister C, Italtatsch R, Henklien P (1991) Polyclonal antibodies from chicken egg yolk-an alternative to the production of mammalian IgG type antibodies in rabbits. *ATLA* 19: 403-419.
2. Justiz Vaillant AA, Akpaka PE, Smikle M, McFarlane-Anderson N (2012) In vitro Inhibition of *Staphylococcus aureus* Isolates by Anti-Anti-Idiotypic Antibodies to Staphylococcal Protein (SpA). *J Vaccines Vaccin* 3: 127-131.
3. Collins AM, Robertson DM, Hosking CS, Flannery GR (1991) Oral immunization with xenogeneic antibodies stimulates the production of systemic and mucosal anti-idiotypic antibodies. *Immunology* 73: 388-393.
4. Ho DD, Kaplan JC, Rackauskas IE, Gurney ME (1988) Second conserved domain of gp120 is important for HIV infectivity and antibody neutralization. *Science* 239: 1021-1023.
5. Polson A (1990) Isolation of IgY from the yolks of eggs by a chloroform polyethylene glycol procedure. *Immunol Invest* 19: 253-258.

## Materials

### MATERIALS

⊗ 5 × 2ml No-Waste™ Freund's Incomplete Adjuvant (FIA) **G-Biosciences Catalog #786-099**

⊗ 2ml No-Waste™ Freund's Complete Adjuvant (FCA) **G-Biosciences Catalog #786-709**

⊗ Fragment 254-274 of Gp120 of HIV (peptide)

- 1 Two healthy layer chickens (brown Leghorn), aged approximately 6 months, are injected intramuscularly at multiple sites on the breast with the peptide-keyhole limpet hemocyanin (KLH) conjugate.
- 2 The chickens are immunized on day 0, with 0.2  $\mu\text{mol/ml}$  of the fragment 254-274 of HIV gp120-conjugated to KLH (immunogen) in 0.5 ml complete Freund's adjuvant (Sigma-Aldrich Co, St. Louis Missouri).
- 3 On days 15, 60, and 90 chickens are immunized with 0.2  $\mu\text{mol/ml}$  of the immunogen in 0.5 ml incomplete Freund's adjuvant (Sigma-Aldrich Co, St. Louis Missouri).
- 4 The eggs are collected post-immunization. The immunoglobulin Y is separated using the Polson method (1990).
- 5 The 96 well microtitre plate is coated overnight at 4°C with 1  $\mu\text{g}/\mu\text{l}$  per well of a mixture of Keyhole limpet hemocynin (KLH) in carbonate-bicarbonate buffer pH 9.6.
- 6 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 7 50  $\mu\text{l}$  of water soluble fraction from eggs is added and incubated for 1h at room temperature and the microplate is rewashed 4X with PBS-Tween.
- 8 Then 50  $\mu\text{l}$  of rabbit anti-chicken IgY-peroxidase conjugate diluted 1:30,000 in PBS-non-fat milk is added to each well and incubated for 1h at RT. The plate is washed 4X with PBS-Tween.
- 9 50  $\mu\text{l}$  of 4 mg/ml o-phenylenediamine solution (OPD) is added and the plate is incubated 15 minutes at RT in the dark.
- 10 The reaction is stopped with 50  $\mu\text{l}$  of 3M H<sub>2</sub>SO<sub>4</sub> solution.
- 11 The plate is visually assessed for the development of colour and read in a microplate reader at 492 nm.
- 12 A cut-off point should be calculated as the mean of the optical density of negative controls x 2 SD. The higher the OD value the higher will be the quantities of anti-KLH in hens.