

Aug 25, 2025

ChIP-seq library generation

 [Nature Cancer](#)

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbwz73gpk/v1>

hanjingnk¹

¹Tsinghua university



hanjingnk

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbwz73gpk/v1>

External link: <https://doi.org/10.1038/s43018-025-01065-3>

Protocol Citation: hanjingnk 2025. ChIP-seq library generation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgbwz73gpk/v1>

**Manuscript citation:**

Han, J., Lu, X., Guo, M. *et al.* Spatiotemporal control of SMARCA5 by a MAPK–RUNX1 axis distinguishes mutant KRAS-driven pancreatic malignancy from tissue regeneration. *Nat Cancer* (2025). <https://doi.org/10.1038/s43018-025-01065-3>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 23, 2025

Last Modified: August 25, 2025

Protocol Integer ID: 225309

Keywords: seq library generation chip

Abstract

ChIP-seq library generation

Guidelines

Protocol summary (as presented on this page):

1. Collect 10 million “806” cells and fix with 1% formaldehyde in DMEM medium at room temperature for 8 min; quench with glycine (final concentration 0.125 M) for 5 min at room temperature.
2. Wash cells with 1 × PBS and resuspend in cold cell lysis buffer (50 mM HEPES pH 7.5, 140 mM NaCl, 0.1% sodium deoxycholate, 1% Triton X-100, 1 mM EDTA, 0.1% SDS, 0.25% sarkosyl, 1 × complete protease inhibitor cocktail).
3. Sonicate cells with 10 sec on / 20 sec off cycles for a total of 15 min.
4. Centrifuge chromatin (10 min, 13,000 × g, 4 °C) and dilute supernatants with dilution buffer (50 mM HEPES pH 7.5, 140 mM NaCl, 1% Triton X-100, 1 mM EDTA, 1 × complete protease inhibitor cocktail) for immunoprecipitation.
5. Incubate Dynabeads (Invitrogen) with antibody in high salt buffer (20 mM Tris-HCl pH 7.9, 500 mM NaCl, 2 mM EDTA, 1% Triton X-100, 0.1% SDS, 1 × complete protease inhibitor cocktail) for 3 hr at 4 °C with rotation.
6. Add antibody/beads to chromatin and incubate overnight at 4 °C with rotation.
7. Wash beads 7 times with high salt buffer and twice with TE buffer (10 mM Tris-HCl pH 8.0, 1 mM EDTA).
8. Elute DNA/protein complexes twice for 15 min at 65 °C in 100 µL elution buffer (10 mM Tris-HCl pH 8.0, 1 mM EDTA, 200 mM NaCl, 1% SDS).
9. Incubate supernatants with RNase A and proteinase K and reverse crosslink at 65 °C for 6 h or overnight.
10. Purify DNA by phenol/chloroform extraction and ethanol precipitation.
11. Resuspend DNA in water and repair to 5' phosphorylated and 3' dA-tailed DNA using NEBNext® Ultra™ II End Repair/dA-Tailing Module (NEB, E7546S).
12. Ligate barcodes (BIOO Scientific, 514153) using Quick Ligation kit (NEB, M2200).
13. After size selection, amplify library with 16 cycles of PCR.
14. Purify PCR-amplified libraries using KAPA beads (Kapa Biosystems) and send for sequencing.

Materials

- 1% formaldehyde (in DMEM medium)
- Glycine (final concentration 0.125 M)
- 1 × PBS
- Cold cell lysis buffer (50 mM HEPES pH 7.5, 140 mM NaCl, 0.1% sodium deoxycholate, 1% Triton X-100, 1 mM EDTA, 0.1% SDS, 0.25% sarkosyl, 1 × complete protease inhibitor cocktail)
- Sonicator (capable of 10 sec on/20 sec off cycles)
- Centrifuge (able to reach 13,000 × g at 4 °C)
- Dilution buffer (50 mM HEPES pH 7.5, 140 mM NaCl, 1% Triton X-100, 1 mM EDTA, 1 × complete protease inhibitor cocktail)
- High salt buffer (20 mM Tris-HCl pH 7.9, 500 mM NaCl, 2 mM EDTA, 1% Triton X-100, 0.1% SDS, 1 × complete protease inhibitor cocktail)
- Dynabeads (Invitrogen)
- Antibodies (appropriate for target protein A or protein G IP)
- TE buffer (10 mM Tris-HCl pH 8.0, 1 mM EDTA)
- Elution buffer (10 mM Tris-HCl pH 8.0, 1 mM EDTA, 200 mM NaCl, 1% SDS)
- RNase A
- Proteinase K
- Phenol/chloroform for DNA purification
- Ethanol (for precipitation)
- Nuclease-free water (for resuspension)
- NEBNext® Ultra™ II End Repair/dA-Tailing Module (NEB, E7546S)
- Barcodes (BIOO Scientific, 514153)
- Quick Ligation kit (NEB, M2200)
- KAPA beads (Kapa Biosystems)
- PCR reagents and thermocycler (used for 16 cycles after size selection)
- Consumables: tubes, rotator, magnetic rack for beads

Before start

- Starting material: 10 million “806” cells (collect and proceed to fixation as described).
- Prepare fixation (1% formaldehyde in DMEM) and quench (glycine final 0.125 M) solutions prior to use.
- Prepare all buffers (cold cell lysis buffer, dilution buffer, high salt buffer, TE buffer, elution buffer) and have protease inhibitors on ice.

- 1 Collect 10 million "806" cells and fix with 1% formaldehyde in DMEM medium at room temperature for 8 min; quench with glycine (final concentration 0.125 M) for 5 min at room temperature.
- 2 Wash cells with 1 × PBS and resuspend in cold cell lysis buffer (50 mM HEPES pH 7.5, 140 mM NaCl, 0.1% sodium deoxycholate, 1% Triton X-100, 1 mM EDTA, 0.1% SDS, 0.25% sarkosyl, 1 × complete protease inhibitor cocktail).
- 3 Sonicate the cells with 10 sec on / 20 sec off cycles for a total of 15 min.
- 4 Centrifuge the chromatin (10 min, 13,000 × g, 4 °C) and dilute the supernatants with dilution buffer (50 mM HEPES pH 7.5, 140 mM NaCl, 1% Triton X-100, 1 mM EDTA, 1 × complete protease inhibitor cocktail) for immunoprecipitation.
- 5 Prepare Dynabeads (protein A or protein G) by incubating with antibody in high salt buffer (20 mM Tris-HCl pH 7.9, 500 mM NaCl, 2 mM EDTA, 1% Triton X-100, 0.1% SDS, 1 × complete protease inhibitor cocktail) for 3 hr at 4 °C with rotation.
- 6 Add the antibody-bound beads to the chromatin and incubate overnight at 4 °C with rotation.
- 7 Wash the beads 7 times with high salt buffer and then wash twice with TE buffer (10 mM Tris-HCl pH 8.0, 1 mM EDTA).
- 8 Elute the DNA/protein complexes twice for 15 min at 65 °C in 100 µL elution buffer (10 mM Tris-HCl pH 8.0, 1 mM EDTA, 200 mM NaCl, 1% SDS).
- 9 Incubate the combined supernatants with RNase A and proteinase K and reverse crosslink at 65 °C for 6 h or overnight.
- 10 Purify DNA by phenol/chloroform extraction and ethanol precipitation.
- 11 Resuspend DNA in nuclease-free water and perform end repair/dA-tailing to generate 5' phosphorylated and 3' dA-tailed DNA using NEBNext® Ultra™ II End Repair/dA-Tailing Module (NEB, E7546S).
- 12 Ligate barcodes (BIOO Scientific, 514153) using the Quick Ligation kit (NEB, M2200).
- 13 After size selection, amplify the library by PCR (16 cycles).



- 14 Purify the PCR-amplified libraries using KAPA beads (Kapa Biosystems) and submit the libraries for sequencing.