

Dec 04, 2024

ChIP-seq protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzkb62vx1/v1>

Renhe Luo¹, Michael Beer², Danwei Huangfu¹

¹Sloan Kettering Institute; ²Johns Hopkins University

IGVF



Michael Beer

Johns Hopkins University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzkb62vx1/v1>

Protocol Citation: Renhe Luo, Michael Beer, Danwei Huangfu 2024. ChIP-seq protocol. protocols.io
<https://dx.doi.org/10.17504/protocols.io.rm7vzkb62vx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 03, 2024

Last Modified: December 04, 2024

Protocol Integer ID: 113702

Keywords: seq protocol for esc, seq protocol chip, seq protocol, esc, seq, de differentiation, protocol, differentiation

Funders Acknowledgements:

NHGRI

Grant ID: HG012051

Abstract

ChIP-seq protocol for ESC-DE differentiation

ChIP-seq protocol

- 1 For each sample, around 30 million cells were crosslinked with 1% formaldehyde and quenched with 0.125M glycine.
- 2 Fixed cells were then lysed in 700uL SDS buffer (1% SDS, 10 mM EDTA, 50 mM Tris-HCl, pH 8), and incubated for 10min on ice.
- 3 Sonication was performed on a Branson Sonifier 150 set at 30% amplitude for 5min30s total on (10s on/off pulsing).
- 4 Clear supernatant was collected for antibody binding overnight, followed by Dynabeads (Thermo Fisher Scientific;10004D) incubation for 6 hours at 4°C.
- 5 Then the beads were pelleted and washed twice with low salt (0.1% SDS, 1% Triton X-100, 2 mM EDTA, 20 mM Tris-HCl, pH 8, 150 mM NaCl), high salt (0.1% SDS, 1% Triton X-100, 2 mM EDTA, 20 mM Tris-HCl, pH 8, 500 mM NaCl), and TE buffer (10 mM Tris-HCl, pH 8, 1 mM EDTA) respectively.
- 6 The DNA was eluted from the beads by incubating in elution buffer (1% SDS, 0.1 M NaHCO₃) at 65°C for 15min and decrosslinked with 5M NaCl at 65°C overnight.
- 7 A total of 10 µl 0.5 M EDTA, 20 µl 1 M Tris-HCl, pH 6.5, and 1 µl Proteinase K (20 mg ml⁻¹) were added to decrosslinked product and incubated for 1 h at 45°C.
- 8 DNA was isolated by using QIAquick PCR purification kit (Qiagen, 28104; QIAGEN).
- 9 Then the sequencing library was generated by using the NEBNext® Ultra II DNA Library Prep Kit (New England Biolabs; E7103S) and NEBNext® Multiplex Oligos for Illumina® (New England Biolabs Index Primers Set 1; NEB, E7335S).
- 10 Samples were pooled and submitted to MSKCC Integrated Genomics Operation core for quality control and sequencing on Illumina HiSeq 4000 platform.