

Apr 08, 2025

ChIPseq libraries preparation

 [Nature Cancer](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v97j21g3e/v1>

Laura Antonucci¹, Isidoro Cobo¹, Christopher K. Glass¹, Michael Karin¹

¹University of California, San Diego



Laura Antonucci

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v97j21g3e/v1>

External link: <https://doi.org/10.1038/s43018-025-01003-3>

Protocol Citation: Laura Antonucci, Isidoro Cobo, Christopher K. Glass, Michael Karin 2025. ChIPseq libraries preparation. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.x54v97j21g3e/v1>

Manuscript citation:

Antonucci, L., Li, N., Duran, A. *et al.* Self-amplifying NRF2–EZH2 epigenetic loop converts KRAS-initiated progenitors to invasive pancreatic cancer. *Nat Cancer* **6**, 1263–1282 (2025). <https://doi.org/10.1038/s43018-025-01003-3>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 08, 2025

Last Modified: April 08, 2025

Protocol Integer ID: 126322

Keywords: chipseq libraries preparation chip, h3k27me3 in ctrl, h3k27ac, h3k27me3, chromatin, specific antibody, using specific antibody, antibody, hiseq, cell, chip

Abstract

ChIP-seq was performed to analyze H3K27ac and H3K27me3 in CTRL and EZH2KO UN-KC6141 cells. Cells were crosslinked, lysed, and sonicated. Chromatin was immunoprecipitated using specific antibodies and Protein G Dynabeads. Libraries were prepared on beads using the NEBNext Ultra II kit, PCR-amplified, size-selected (200–500 bp), and sequenced on HiSeq 4000 or NextSeq 500.

- 1 ChIPseq experiments for H3K27ac and 3meH3K27 were performed as described elsewhere(1, 2). Briefly, CTRL and EZH2^{KO} UN-KC6141 cells were fixed with 1% formaldehyde (ThermoFisher, #28906) at room temperature for 10 minutes. Cells were washed once again with 0.01% Triton-X-100 in PBS, pelleted, snap frozen, and stored at -80°C.
Cells were thawed and lysed in a buffer containing 10 mM Tris/HCl pH 7.5, 100 mM NaCl, 1 mM EDTA, 0.5 mM EGTA, 0.1% deodeoxycholate, 0.5% sarkosyl, and protease inhibitor cocktail (SIGMA). Lysate was sonicated using a Covaris for 12 cycles with the following settings: time, 60 s; duty, 5.0; PIP, 140; cycles, 200; amplitude, 0.0; velocity, 0.0; dwell, 0.0. One percent of sonicated lysate was kept as ChIP input for analysis.
Immunoprecipitation mix consisting of Protein G Dynabeads (Invitrogen, #10003D) and 1 µg of H3K27ac (Active Motif, #39133) or 3meH3K27 (Abcam, ab195477). Beads were washed four times each with wash buffer I (20 mM tris-HCl, 150 mM NaCl, 0.1% SDS, 1% Triton X-100, and 2 mM EDTA), wash buffer II (10 mM tris-HCl, 250 mM LiCl, 1% IGEPAL CA-630, 0.7% Na-deoxycholate, and 1 mM EDTA), tris-EDTA (TE) 0.2% Triton X-100, and TE 50 mM NaCl and then resuspended 25µl of TT buffer (10 mM tris-HCl (pH 8.0) and 0.05% Tween 20, and sequencing libraries were prepared as following.
ChIP libraries were prepared while bound to Dynabeads using a NEBNext Ultra II Library preparation kit (NEB) using half reactions. DNA was polished, poly(A)-tailed, and ligated using dual UDI (IDT) or single (Bioo Scientific). Then crosslinks were reversed by adding 33.5 µl of a mix containing 18.4 µl nuclease free water, 4 µl 10% SDS, 3 µl 0.5 M EDTA, 1.6 µl 0.2M EGTA, 1 µl 10mg/ml Proteinase K (Biolabs, #P8107S), 1 µl 10mg/ml RNase A, 4.5 µl 5M NaCl by incubating at 55°C for 1h followed by 75°C for 30 min. Dynabeads were removed, and libraries were cleaned by adding 3 µl of SpeedBeads in 125 µl of 20% PEG 8000/1.5M NaCl, washed by adding 120 µl 80% etOH, air dried, and eluted in 12.5 µl of buffer containing 10 mM Tris/HCl pH 8.0 and 0.05% Tween-20. DNA was PCR-amplified for 15 cycles using NEBNext Ultra II PCR master mix using Solexa 1GA and Solexa 1GB primers. Libraries were size selected 200-500 bp by running in 10% TBE acrylamide gels (ThermoFisher, #EC62752BOX) and sequenced using a HiSeq 4000 or a NextSeq 500.

Protocol references

1. D. Gosselin *et al.*, An environment-dependent transcriptional network specifies human microglia identity. *Science* **356** (2017).
2. J. S. Seidman *et al.*, Niche-Specific Reprogramming of Epigenetic Landscapes Drives Myeloid Cell Diversity in Nonalcoholic Steatohepatitis. *Immunity* **52**, 1057-1074 e1057 (2020).