

Apr 11, 2024

Version 1

Chloroform-free DNA Extraction - Ammonium Acetate Precipitation Method V.1

DOI

dx.doi.org/10.17504/protocols.io.x54v9pekzg3e/v1

NERC Environmental Omics Facility (NEOF) Visitor Facility¹

¹University of Sheffield



NEOF - NERC Environmental Omics Facility VF

University of Sheffield

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9pekzg3e/v1>

Protocol Citation: NERC Environmental Omics Facility (NEOF) Visitor Facility 2024. Chloroform-free DNA Extraction - Ammonium Acetate Precipitation Method. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9pekzg3e/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We have been regularly using this method for decades.

Created: June 12, 2023

Last Modified: April 11, 2024

Protocol Integer ID: 83246

Keywords: dna extraction, effective dna extraction method, chloroform, ammonium acetate precipitation method, dna, extraction

Abstract

A chloroform-free, cost-effective DNA extraction method for a variety of sample types.

Materials

Digsol: To make 500ml

- 20ml 0.5M EDTA (pH 8.0)
- 3.425g NaCl
- 25ml 1M Tris-HCl (pH 8.0)
- 430ml ddH₂O


Autoclave then add 25ml of 20% SDS.


Low TE (Tris_{10mM}, EDTA_{0.1mM}): To make 500ml


- 5ml 1M Tris-HCl (pH 8.0)
- 100μl 0.5M EDTA (pH8.0)
- 495ml ddH₂O

Autoclave.

Protocol materials

 Proteinase K, 2mL Qiagen Catalog #19131

 1M DTT Merck MilliporeSigma (Sigma-Aldrich) Catalog #43816

 1M DTT Merck MilliporeSigma (Sigma-Aldrich) Catalog #43816

 Ammonium Acetate Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1542-500G

Troubleshooting

Safety warnings



Use GLP and wear protective equipment. Avoid contact with skin.

All chemicals can be disposed of down the sink with copious amounts of water to dilute the ethanol down to >20% of the waste.
















If inhaled: If unconscious, place in recovery position and seek medical advice. Keep the respiratory tract clear. If symptoms persist, call a physician.

In case of skin contact: Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. If symptoms persist, call a physician.










In case of eye contact: Remove contact lenses. Protect unharmed eye. Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed: If accidentally swallowed, obtain immediate medical attention. Rinse mouth with water. Never give anything by mouth to an unconscious person.



- 1 Add  250 μL Digsol buffer (see Materials for recipe) and  10 μL  Proteinase K, 2mL **Qiagen Catalog #19131** to a labelled 1.5ml tube. 1m
- 2 **Tissue:** Cut into small pieces ($<1\text{cm}^2$) with a sterile razor blade on a sterile glass plate before adding to the 1.5mL tube. 5m
Blood: Centrifuge blood sample at 13,000rpm for about 1 min (to pellet sample). Remove sample from ethanol with toothpick and blot onto tissue. When almost dry, transfer the toothpick into the 1.5mL tube and jiggle to dislodge the blood. Remove toothpick and place in disinfectant.
Swab: Dry of ethanol and place into the 1.5mL tube. Snap off the end of the swab so the lid may close.
Feather: Cut the calamus of 1-3 feathers into small pieces with a sterile razor blade on a sterile glass plate before adding to the tube. If feathers are very small and blood spots are present on the tips, feathers can be added whole.
Hair: Place hair(s) into the 1.5mL tube preferably with the largest root and add  5 μL  1M DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #43816**
Insect: Degut and add to 1.5mL tube. If the insect is very small, starve for 48hrs before freezing. Add  5 μL  1M DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #43816** . Crush with a pestle or add a metal lysing bead and place on TissueLyser (Qiagen) for 2min, 30/s.
- 3 Vortex and place samples in a rotating oven at  55 $^{\circ}\text{C}$ for  03:00:00 or  Overnight for maximum digestion. 3h
- 4 Add  300 μL 4M  Ammonium Acetate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1542-500G** . 1m
- 5 Vortex several times over a period of at least  00:15:00 at  Room temperature . 15m
- 6  13000 rpm, 00:10:00 10m



- 7 Aspirate the supernatant (clear liquid containing the DNA) into a new labelled 1.5ml tube. Discard tube containing the pelleted protein debris. 1m
- 8 Add  1 mL 100% ethanol and invert the tube gently several times to precipitate DNA. 1m
- 9  13000 rpm, 00:10:00 10m
- 10 Pour off ethanol in a smooth motion, taking care not to lose the DNA pellet. 1m
- 11 Add  500 μ L 70% ethanol and invert gently several times to clean the pellet. 1m
- 12  15000 rpm, 00:05:00 5m
- 13 Pour off ethanol in a smooth movement or using a pipette gently draw off the supernatant if fear of losing the pellet. Stand tubes upside-down on clean tissue until dry (approx. **30-60 minutes**). This can be sped up by using the heat of a lamp from above. 30m
- 14 Once fully dry add approx.  100 μ L Low TE (see Materials for recipe). Add less if a very tiny pellet or no pellet is observed. Flick sample to dislodge pellet. 1m
- 15 Place tubes in a thermomixer or shaking oven at  50 °C for  00:30:00 to dissolve the pellet. If the pellet has not completely dissolved add more Low TE. 30m
- 16 Store at  -20 °C (long term) or  4 °C (short term).



Protocol references

<https://onlinelibrary.wiley.com/doi/epdf/10.1046/j.0962-1083.2001.01355.x>

Richardson DS, Jury FL, Blaakmeer K, Komdeur J, Burke T (2001) Parentage assignment and extra-group paternity in a cooperative breeder: the Seychelles warbler (*Acrocephalus sechellensis*). *Molecular Ecology***10**, 2263-2273.

<https://onlinelibrary.wiley.com/doi/abs/10.1034/j.1600-048X.2000.310208.x>

The Evolution of Cooperative and Pair Breeding in Thornbills *Acanthiza* (Pardalotidae)

James A. Nicholls, Michael C. Double, David M. Rowell and Robert D. Magrath

Journal of Avian Biology

Vol. 31, No. 2 (Jun., 2000), pp. 165-176