


Mar 14, 2023

Version 2

## Choanoflagellate Ciliogenesis Live Imaging V.2

DOI

[dx.doi.org/10.17504/protocols.io.q26g7y9n3gwz/v2](https://dx.doi.org/10.17504/protocols.io.q26g7y9n3gwz/v2)



Maxwell C Coyle<sup>1</sup>

<sup>1</sup>University of California Berkeley

King Lab



Maxwell C Coyle

Harvard University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7y9n3gwz/v2>

**Protocol Citation:** Maxwell C Coyle 2023. Choanoflagellate Ciliogenesis Live Imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7y9n3gwz/v2> Version created by **Maxwell C Coyle**

**Manuscript citation:**

Coyle, M. C. *et al.* An RFX transcription factor regulated ciliogenesis in the progenitors of choanoflagellates and animals. *bioRxiv* 2022.11.11.515474 (2022) doi:10.1101/2022.11.11.515474

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**This protocol works well, although de-ciliation efficiency is in the 80-90% range, not 100. Also you may need to adjust the timing of the -20C step depending on the exact temp and heat exchange of your freezer. The goal is to go as long as possible before the solution freezes.**

**Created:** March 05, 2023

**Last Modified:** March 14, 2023

**Protocol Integer ID:** 78159

**Keywords:** flagella from choanoflagellate cell, ciliogenesis live imaging this protocol, choanoflagellate cell, ciliogenesis live imaging, portable into other choanoflagellate species, ciliary removal, other choanoflagellate species, ciliogenesis, cilia, flagella, cells for live imaging, species *salpingoeca rosetta*, cell

**Funders Acknowledgements:**

HHMI

## Disclaimer

### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **[protocols.io](#)** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **[protocols.io](#)**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

This protocol removes the cilia/flagella from choanoflagellate cells and sets up the cells for live imaging of ciliogenesis. It has been developed for the species *Salpingoeca rosetta*, but may be portable into other choanoflagellate species. Cells begin to re-generate their cilia/flagella right after removal. The idea of using a cold shock with glycerol for ciliary removal came from Brokaw et al 1960 (doi: 10.1016/0014-4827(60)90027-6).

## Materials

High Nutrient Media: 4% AKCGM3 + 4% AKSWC in AKSW - Artificial Known Sea Water (See Booth 2018 Molecular Biology of the Cell for sea water details).

AKSW without addeed supplements

Incubator

Tabletop centrifuge

75cm<sup>2</sup> vented flasks

Haemocytomer or automated cell counter

16% Paraformaldehyde

0.1 mg/ml Poly-D-lysine

Forceps

Surface corona treater

Fluorodishes

Circular (22 mm diameter) coverslips

70% ethanol

50% glycerol

-20C freezer for incubation

Widefield microscope

## Protocol materials

 Fluorodish **World Precision Instruments Catalog #FD35-100**

 Poly-D-lysine hydrobromide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P6407-5MG**

## Troubleshooting




## Concentrate cells and remove cilia

- 1 Grow choanoflagellate cells (*Salpingoeca rosetta* fed with *Echinicola pacifica*, ATCC PRA-390) in High Nutrient Media to a density of  $1\text{--}2 \times 10^6$  cells/ml. Grow at 22 °C, 60% humidity

We grow 30 ml of culture in 75 cm<sup>2</sup> vented flask. Typically, inoculating this flask with a choanoflagellate cell density of 8,000 cells/ml 48 hours before ciliogenesis works well.

- 2 Count cells by haemocytometer or automated cell counter\*. Shake culture flask vigorously to homogenize cell population and then mix 99 µl of cell culture with 1 µl of 16% paraformaldehyde to fix cells for counting. Typically 10 µl of fixed cells can be loaded into a haemocytometer or automated cell counting slide.

\*We use LUNA-FL Automated Fluorescence Cell Counter (Logos Biosystems L20001)

- 3 Aliquot and pellet  $6 \times 10^6$  cells  2000 x g, 00:10:00

- 4 Corona treat a fluorodish 5–10 seconds

 Fluorodish **World Precision Instruments Catalog #FD35-100**

### Equipment

BD-20AC Laboratory Corona Treater

NAME

Corona Treater

TYPE

Electro-Technic Products

BRAND

12051A

SKU

<https://www.electrotechnicproducts.com/bd-20ac-laboratory-corona-treater/><sup>LINK</sup>

- 5 Rinse fluorodish for 5 seconds with 1 ml of 0.1 mg/ml poly-D-lysine, followed by 3x washes with 1 ml water. Dry by air or by capillary action of kimwipe, being careful to minimize contact with surface of imaging dish.



Poly-D-lysine hydrobromide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P6407-5MG**

6 Rinse a coverslip (circular - 22mm diameter) in 70% EtOH followed by water and then lay on kimwipe to dry. Easiest to hold coverslip by forceps and dunk into 50 ml conical tubes with the ethanol or water.

1m

7 When cells are done pelleting, remove supernatant and resuspend cell pellet in 800  $\mu$ l of AKSW and transfer to 1.5 ml eppendorf tube.

1m

8 Add 200  $\mu$ l of 50% glycerol (final concentration: 10% glycerol) and mix by pipetting

1m

9 Add cells to a second Fluorodish (i.e. one not treated with poly-D-lysine) and incubate

-20 °C 7 mins

8m



Standard laboratory freezer is fine, but depending on exact temperature of your freezer or where in the freezer you place the cells, you may need to adjust the timing.

## Set up cells for live imaging of ciliogenesis

13m

10 Transfer cells to 1.5 ml eppendorf tube and pellet 4200 x g, 00:08:00

8m

11 Remove supernatant and resuspend cells in 25  $\mu$ l of AKSW

1m

12 Transfer cells to Fluorodish coated with poly-D-lysine and lay clean coverslip slowly on top using forceps

1m

13 Mount dish on microscope\* and find focus. Let cells settle for 1 minute.

2m

\*We use a Zeiss Axio Observer.Z1/7 widefield with a 100x NA 1.40 Plan-Apochromatic oil immersion objective and a Hamamatsu Orca-Flash 4.0 LT CMOS digital camera

14 Float coverslip off of cells by adding AKSW around the side of the coverslip drop by drop with a plastic transfer pipette. If you don't do this, the cells will eventually suffocate.

1m



## 15 Image!



On our system we use a short (5 ms exposure) with high light intensity (12.2 V) and a DIC condenser to get the best imaging of ciliogenesis. We use Zeiss Definite Focus and take a 10  $\mu\text{m}$  z-stack with 1  $\mu\text{m}$  between slices every 30 seconds for one hour.

## Protocol references

Brokaw, C. J. Decreased adenosine triphosphatase activity of flagella from a paralyzed mutant of *Chlamydomonas moewusii*. *Exp. Cell Res.***19**, 430–432 (1960)

Booth, D. S., Szmidt-Middleton, H. & King, N. Transfection of choanoflagellates illuminates their cell biology and the ancestry of animal septins. *Mol. Biol. Cell***29**, 3026–3038 (2018)