

Jun 20, 2017

🌐 Chrome Azural S (CAS) Plate Assay for Iron-Binding Compounds

DOI

dx.doi.org/10.17504/protocols.io.idvca66



Dr. Steven Wilhelm

The Aquatic Microbial E...



Steven W Wilhelm

The University of Tennessee, Knoxville

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.idvca66

Protocol Citation: Dr. Steven Wilhelm 2017. Chrome Azural S (CAS) Plate Assay for Iron-Binding Compounds. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.idvca66>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 10, 2017

Last Modified: March 21, 2018

Protocol Integer ID: 6293


Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.


Modified from Schwyn, B. & Neilands, J. B. Universal chemical assay for the detection and determination of siderophores. *Anal Biochem*, 160:147-156 (1987).


Materials

STEP MATERIALS


 Piperazine anhydrous

 Iron(III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #44944**

 CTAB (Hexadecyltrimethylamm onium bromide) **BBI Biotech Catalog #CB0108-100g**


 Piperazine anhydrous


 Iron(III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #44944**

 CTAB (Hexadecyltrimethylamm onium bromide) **BBI Biotech Catalog #CB0108-100g**


Protocol materials


 Iron(III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #44944**


 CTAB (Hexadecyltrimethylamm onium bromide) **BBI Biotech Catalog #CB0108-100g**

 Piperazine anhydrous


 Iron(III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #44944**

 CTAB (Hexadecyltrimethylamm onium bromide) **BBI Biotech Catalog #CB0108-100g**

 Piperazine anhydrous


 Piperazine anhydrous


 Iron(III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #44944**


 CTAB (Hexadecyltrimethylamm onium bromide) **BBI Biotech Catalog #CB0108-100g**


Assay Solution Preparation


- 1 Take 6 mL 10 mM CTAB (HDTMA) stock and mix with 40 mL Milli-Q H₂O

 CTAB (Hexadecyltrimethylamm onium bromide) **BBI Biotech Catalog #CB0108-100g**
- 2 Mix 1.5 mL FeCl₃-HCl stock (1 mM FeCl₃ dissolved in 10 mM HCl) with 7.5 mL 2 mM CAS dye stock

 Iron(III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #44944**
- 3 Slowly add Fe-CAS mixture to CTAB solution while stirring
- 4 Add 6.5 mL of 12 N HCl slowly to 25 mL Milli-Q H₂O
- 5 Add 4.3 g anhydrous piperazine to the acid solution


 4.3 g

 Piperazine anhydrous
- 6 Mix the piperazine acid solution slowly into the Fe-CAS*CTAB solution
- 7 Bring CAS solution to 100 mL final volume


 100 mL
- 8 Add enough 5-sulfosalicyclic acid to get 4 mM final concentration

Assay Part I

- 9 Add 60.5 mg of 2 mM CAS dye to 50 mL Milli-Q H₂O
- 10 Add 10 mL FeCl₃ solution

 10 mL
- 11 Dissolve 73 mg of 10 mM CTAB (HDTMA) in 40 mL Milli-Q H₂O



- 12 Mix the two solutions together
- 13 Bring solution to 1 L with Milli-Q H₂O
 1 L
- 14 Autoclave at 121°C for 20 min

Assay Part II

- 15 Make medium depending on the type of microorganism that you are trying to grow.

Note

For example, add ESAW-FE with casamino acids, peptone and agar for marine microorganisms.

- 16 Autoclave medium at 121°C for 20 min
- 17 Once both part I and part II are cooled (55°C), mix the two solutions together and pour plates.
- 18 Incubate plates according to the microorganism. Screen plates once colonies are grown: yellow halos = siderophore and blue = no siderophore