


Feb 02, 2023

Version 1

# 🌐 Chronic Recoverable Neuropixels in Mice V.1

DOI

[dx.doi.org/10.17504/protocols.io.e6nvwo87lmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvwo87lmk/v1)



Emily A Aery Jones<sup>1</sup>

<sup>1</sup>Stanford University



Emily A Aery Jones

Stanford University

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwo87lmk/v1>

**Collection Citation:** Emily A Aery Jones 2023. Chronic Recoverable Neuropixels in Mice. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e6nvwo87lmk/v1>

**License:** This is an open access collection distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this collection and it's working**

**Created:** February 02, 2023

**Last Modified:** February 02, 2023

**Collection Integer ID:** 76314

**Keywords:** electrophysiology, Neuropixels, silicon probe, hippocampus, entorhinal cortex, spatial navigation, freely moving recording, electrode, chronic recoverable neuropixels in mice, chronic recoverable neuropixel, neuropixels 2021 course lecture, other chronic recoverable design, probes for future use, neuropixel, probe during surgery, mice, beta probe, same probe, implantation into any brain region, probe, probe adjustable before gluing,

**Funders Acknowledgements:**

**Simons Collaboration on the Global Brain**

Grant ID: 542987SPI

**Stanford School of Medicine Dean's Postdoctoral Fellowship**

**A.P. Giannini Foundation Postdoctoral Fellowship**

## Abstract

This protocol collection explains how to build a low-cost, lightweight system to implant 1 Neuropixels 1.0 probe or 2 Neuropixels 2.0 probes into mice, record during freely moving behavior, then recover the probes for future use.

### Other chronic recoverable designs:

- single 1.0 probe in mice: Juavinett et al, 2019 (<https://elifesciences.org/articles/47188>)
- multiple 1.0 probes in rats: Luo & Bondy et al, 2020 (<https://elifesciences.org/articles/59716>)
- single 2.0 probe in mice and dual or moveable 2.0 probes in rats: van Daal, Aydin, & Michon et al, 2021 (<https://www.nature.com/articles/s41596-021-00539-9>)
- overview of unpublished designs: Neuropixels 2021 course lecture 3.7 by Yoh Isogai ([https://www.youtube.com/watch?v=7cZqYliGvBQ&list=PLfhWmWntvjI64ti\\_a-MzHlwqwEU0ZIALb&index=18](https://www.youtube.com/watch?v=7cZqYliGvBQ&list=PLfhWmWntvjI64ti_a-MzHlwqwEU0ZIALb&index=18))

### Advantages of this design:


- Lightweight (entire assembly, including headbar, dental cement, and tape, weighs <3g for single 1.0 probe or <4g for dual 2.0 probes)
- Allows headfixed or freely moving recordings
- Quickly attaches to a headstage holder, which provides LEDs for tracking and a counterweight to encourage running
- Position and angle of internal mount on the probe adjustable before gluing, allowing implantation into any brain region, but custom adjustment is not required, allowing the same probe to be re-inserted into a variety of sites
- Unlike completely enclosed designs, the shank is uncovered during insertion for better visualization, yet doesn't require delicately surrounding the shank with a glue column
- Dual 2.0 probes: probes are independently insertable, so can set custom positions, depths, and angles for each

### Disadvantages of this design:


- Not suitable for larger animals (lightweight design likely can't withstand larger forces, flex cable remains exposed, moisture known to wick up shank into PCB in rats)
- External components are assembled around the probe during surgery rather than during assembly, so surgeries take slightly longer
- Less elegant than completely enclosed designs and requires a larger skull surface area for gluing
- Dual 2.0 probes: how close together the probes can get is limited by some of the components

### Assembly preview:

Single 1.0 probe



Dual 2.0 beta probes



## Troubleshooting



## Files

### Protocol



NAME

**Assembly: Chronic recoverable Neuropixels in mice**

VERSION 7

CREATED BY

**Emily A Aery Jones**  
Stanford University[OPEN](#) →

### Protocol



NAME

**Implant Surgery: Chronic recoverable Neuropixels in mice**

VERSION 5

CREATED BY

**Emily A Aery Jones**  
Stanford University[OPEN](#) →

### Protocol



NAME

**Explant Surgery: Chronic recoverable Neuropixels in mice**

VERSION 3

CREATED BY

**Emily A Aery Jones**  
Stanford University[OPEN](#) →

### Protocol



NAME

**Building a SpikeGLX Rig with camera: Chronic recoverable Neuropixels in mice**

VERSION 2



CREATED BY



Emily A Aery Jones  
Stanford University

OPEN →

## Protocol



NAME

Freely moving recording: Chronic recoverable Neuropixels in mice

VERSION 5

CREATED BY



Emily A Aery Jones  
Stanford University

OPEN →