



Mar 19, 2025

Circular Dichroism Spectroscopy of Clusterin

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr91b2vmk/v1>

Andreas Bracher¹, Patricia Yuste-Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry, Martinsried, Germany



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr91b2vmk/v1>

Protocol Citation: Andreas Bracher, Patricia Yuste-Checa, F Ulrich Hartl 2025. Circular Dichroism Spectroscopy of Clusterin. protocols.io <https://dx.doi.org/10.17504/protocols.io.6qpvr91b2vmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 11, 2025

Last Modified: March 19, 2025

Protocol Integer ID: 124382

Keywords: ASAPCRN, circular dichroism spectroscopy of clusterin, curves of the extracellular chaperone clusterin, extracellular chaperone clusterin, circular dichroism spectroscopy, clusterin, circular dichroism, spectra

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to record circular dichroism (CD) spectra and CD melting curves of the extracellular chaperone Clusterin.

Materials

Buffers and reagents:

- CD buffer: [M] 50 millimolar (mM) K-phosphate  7.0 .

Note

- Prepare 1 M K-phosphate pH 7.0 stock solution by titrating 1 M K_2HPO_4 with 1 M KH_2PO_4 . Dilute stock solution 1:20 with deionized water, sterile-filter and de-gas in vacuum.
- Buffers for CD experiments should not absorb strongly in the far UV range. Chloride salts should be avoided.
- Clusterin variant (see protocol "Clusterin purification from HEK293E cells", [dx.doi.org/10.17504/protocols.io.bvvkn64w](https://doi.org/10.17504/protocols.io.bvvkn64w)).
- Amersham NAP-5 columns (Cytiva, 17085301)

Note

Columns can be re-used when stored at 4 °C.

- De-ionized water
- Methanol

Safety information

Methanol is toxic. During handling, wear gloves and safety goggles.

Instruments

- Nanodrop microvolume spectrophotometer
- Jasco J-715 CD spectrometer equipped with a Peltier-thermostat
- 0.1 mm light path UV quartz cuvette

Preparation

- 1 Equilibrate NAP-5 column at  Room temperature with  10 mL volume of CD buffer according to the manufacturer's protocol using gravity flow.
- 2 Prepare six 1.5 ml plastic vials in a microtube rack for fraction collection.
- 3 Thaw Clusterin variant stock solution aliquot.
- 4 Apply  25 μL Clusterin stock (approximately  200 micromolar (μM)) to NAP-5 column. Allow liquid to completely enter gel bed.
- 5 Apply  475 μL CD buffer. Allow liquid to completely enter gel bed.
- 6 Place first fraction collection vial under the column outlet.
- 7 Add  100 μL CD buffer. Allow liquid to completely enter gel bed. Collect fraction. Move to next vial. Repeat five times. 
- 8 Pool fractions 2-4. Determine absorbance at 280 nm (A_{280}) of pool using Nanodrop, in triplicate. From the A_{280} average value, calculate Clusterin variant concentration using absorbance coefficients calculated with ProtParam (<https://web.expasy.org/protparam/>) and the respective sequence of the processed (starting at residue 23) Clusterin variant.
- 9 Prepare  300 μL solution of  5 micromolar (μM) Clusterin variant for CD melting curve measurement. Dilute with CD buffer. Store  On ice .
- 10 Prepare  300 μL solution of  0.1 μL Clusterin variant (~  2 micromolar (μM)) for CD spectrum measurement. Dilute with CD buffer. Store  On ice .

CD Spectrometer Start-up



- 11 Open nitrogen valve to flush Jasco J-715 CD spectrometer with gaseous nitrogen.
- 12 Start cooling water flow at low rate.
- 13 Switch on detector power source (Power 1).
- 14 Switch on Xenon lamp (Power 2).
- 15 Start control computer and open Spectra Manager program.
- 16 Start external Peltier thermostat (Peltier accessory).

Note

NOTE: For CD spectrometer shut-down, stop the items in steps 11-16 in reverse sequence.

CD Spectrum Measurement

- 17 Start "Spectrum measurement". Select Peltier thermostat under "Accessories" and set the temperature to  20 °C in "Control" under "Accessory".
- 18 Set the "Parameters" like in the following:

Spectrum Measurement - Parameter

Parameters | Data Mode | Data File | Option

Sensitivity: Standard(100mdeg)

Start: 250 nm

End: 195 nm

Data Pitch: 0.1 nm

Scanning Mode: Continuous

Scanning Speed: 50 nm/min

Response: 1 sec

Band Width: 1.0 nm

Slit Width: um

Accumulation: 4

OK Cancel Open... Save

- 19 Wash the UV cuvette with 3 times water and 3 times methanol. Dry with pressurized nitrogen.
- 20 Pipette ~ 300 μ L CD buffer into the cuvette. Do not trap bubbles. Insert into spectrometer cuvette holder. Let cuvette equilibrate for 00:05:00 .
- 21 Record a spectrum of CD buffer alone with above parameters.



5m





22 Wash the cuvette as above (step 19).



23 Pipette ~  300 μL Clusterin variant ( 0.1 μL) CD buffer into the cuvette. Do not trap bubbles. Insert into spectrometer cuvette holder. Let cuvette equilibrate for  00:05:00 .

5m



24 Record a CD spectrum of the Clusterin variant with above parameters.

Note

NOTE: The additional High Tension (HT) signal is a measure for sample absorbance at the wavelength. Ideally it should be below 500; above 800 the CD signal gets too noisy for quantitative measurement.

25 Subtract the CD buffer spectrum from the CD spectrum of the Clusterin variant in Spectra Manager.

26 In Spectra Manager I, under "CD Options" and "Optical Constants", convert the units into molar ellipticities using the molar concentration and the residue number.

27 In Spectra Manager II, under "CD Options" and "CDPro", analyze the secondary structure composition using the SMP56 dataset as a reference and the CONTIN algorithm.

Note

The CONTIN algorithm gives the best fit for Clusterin.




Fig. 1 Comparison of the CD spectrum of wildtype Clusterin (green) with the CONTIN algorithm fitted curve (blue) from CDPro analysis. Residuals are shown in brown.

CD Melting Curve Measurement

- 28 Start "Variable temperature measurement" in Spectra Manager. Select Peltier thermostat under "Accessories" and set the temperature to 4 °C in "Control" under "Accessory".
- 29 Set the "Parameters" like in the following:

Variable Temperature Measurement - Parameter

Parameters | Data Mode | Data File | Option

Wavelength: 222.0 nm

Start: 4.0 C

End: 90.0 C

Data Pitch: 0.1 C

Delay Time: 10 sec

Temperature slope: 60 C/hour

Options

☒ Return to Start Temperature

☐ Use Measure Cell :

☐ Reverse Temperature Scan

Reverse Hold Time 60 sec

Sensitivity: Standard(100mdeg)

Response: 16 sec

Band Width: 1.0 nm

Slit Width: um

OK Cancel Open... Save

- 30 Wash the UV cuvette with 3 times water and 3 times methanol. Dry with pressurized nitrogen.
- 31 Pipette ~ 300 μL Clusterin variant (5 micromolar (μM)) CD buffer into the cuvette. Do not trap bubbles. Insert into spectrometer cuvette holder. Let cuvette equilibrate for 00:05:00 .



5m



**Note**

The higher Clusterin concentration compared to the CD spectrum measurement increases the sensitivity. The absorbance (HT signal) increase at 222 nm wavelength is negligible.

- 32 Record a CD melting curve of the Clusterin variant with above parameters.

Note

A steeper temperature slope will keep the sample further away from equilibrium, resulting in more pronounced overshoot. Heat denaturation of wildtype Clusterin at 5 μ M concentration is not reversible. Therefore, the reverse temperature scan was omitted.

- 33 Under "Processing", "Common Options" and "Denatured Protein", analyze the melting curve.






Fig. 2 Comparison of the CD melting curve of wildtype Clusterin (blue) with the fitted curve (green). The fitting range suggested by the program was used.