

Jun 02, 2020

Citrate-Phosphate Buffer

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.bfydjps6

Neillier Junior¹

¹University of Manitoba



Neillier Junior

University of Manitoba, Universidade Federal de Viçosa



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bfydjps6

Protocol Citation: Neillier Junior 2020. Citrate-Phosphate Buffer. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bfydjps6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 04, 2020

Last Modified: June 02, 2020

Protocol Integer ID: 36581

Keywords: pH, Henderson-Hasselback, Biochemistry, Molecular Biology,



Abstract

A buffer solution has the function of resisting changes in pH even when adding powerful acids or bases. However, in the physiological environment the buffered system also provides cofactors for enzymatic reactions, critical salts and even essential nutrients for cells and tissues. Therefore, when trying to reproduce biological conditions in vitro, we must make the appropriate choice of the buffer. After all, it will provide the appropriate medium in which reactions will occur.

Materials

- Deionized Water
- pH Meter (sensitive)
- Citric Acid
- Dibasic Sodium Phosphate (dihydrate and heptahydrate)

Safety warnings

! Wear personal protective equipment: gloves, lab coat and mask.

Before start

Organize your workspace.

Make sure all solutions and equipment are available.

Citrate-Phosphate Buffer


1

Note

pH range:  2.6 to  7.0

(a) 0.1 M Citric acid; 19.21 g L⁻¹ (M.W. 192.1 g mol⁻¹)

(b) 0.2 M Dibasic sodium phosphate; 35.6 g L⁻¹ (dihydrate; M.W. 178.0 g mol⁻¹) or 53.6 g L⁻¹ (heptahydrate; M.W. 268.0 g mol⁻¹)

Mix citric acid and sodium phosphate solutions in the proportions indicated and adjust the final volume to  100 mL with deionized water.

mL of Citric acid	44.6	39.8	35.9	32.3	29.4	26.7	24.3	22.2	19.7	16.9	13.6	6.5
mL of Sodium phosphate	5.4	10.2	14.1	17.7	20.6	23.3	25.7	27.8	30.3	33.1	36.4	43.6
pH	2.6	3.0	3.4	3.8	4.2	4.6	5.0	5.4	5.8	6.2	6.6	7.0

2 Adjust the final pH using a sensitive pH meter