

Jun 13, 2018

Version 1

## Cloning guides to lentiCRISPR v2 V.1

DOI

[dx.doi.org/10.17504/protocols.io.qx3dxqn](https://dx.doi.org/10.17504/protocols.io.qx3dxqn)

Amit Weiner<sup>1</sup>

<sup>1</sup>University of Toronto



Amit Weiner

University of Toronto

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.qx3dxqn>

**Protocol Citation:** Amit Weiner 2018. Cloning guides to lentiCRISPR v2. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.qx3dxqn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** June 13, 2018



**Last Modified:** June 13, 2018

**Protocol Integer ID:** 13019

**Keywords:** cloning guide rna sequences to the lenticrispr vector, cloning guide rna sequence, guides to lenticrispr v2, cloning guide, lenticrispr v2, lenticrispr vector, rna

## Abstract

I have been struggling with cloning guide RNA sequences to the lentiCRISPR vector (<https://www.addgene.org/52961/>) for months now and have finally figured out a method that consistently works for me, so I felt like I should share this information with whomever is struggling as much as I did!

## Troubleshooting

### Before start

Make Stbl3 chemically competent cells.



## Vector preparation:

- 1 Digest the lentiCRISPRv2 vector:  
3µg vector  
2µl 10X Tango buffer (Thermo Fisher)  
1µl 20mM DTT  
1-1.5µl Esp3I (Thermo Fisher)  
Water to 20µl  
  
In a thermocycler:  
37°C for 4 hours, inactivate at 65°C for 20 mins, keep at 4°C.  
**SKIP** the alkaline phosphatase step.
- 2 Use a gel purification/PCR cleanup kit (Qiagen) **WITHOUT** running on a gel.  
Optional: run ~200ng of the purified vector on a gel to verify digestion.

## Guide insert preparation:

- 3 Anneal and phosphorylate gRNA oligos pair:  
1µl of each oligo (100µM stock)  
1µl 10X T4 ligation buffer (not PNK buffer)  
0.5µl T4 PNK  
Water to 10µl  
  
In a thermocycler: 37°C for 30 mins, 95°C for 5 mins, ramp down to 25°C at 0.1°C/sec (or 5-6°C/min).  
Optional: keep at 4°C.
- 4 Serially dilute the annealed oligos to 1:500

## Ligation:

- 5 Out of ligation at a vector:insert molar ratios of 1:5, 1:10, 1:20, I found that 1:5 works best.  
  
50ng vector  
1.5µl 10X T4 ligation buffer (NEB)  
1µl T4 ligase (NEB)  
2µl diluted oligos (1:500)  
Water to 15µl



Incubate at RT for 1-2 hours.

## Transformation:

- 6 Transform 5µl of the ligation reaction to 50µl Stbl3 chemically competent cells.

If you incubate the Stbl3 cells at 30°C the colonies will be VERY small so look for them carefully. Incubating them at 37°C didn't result in LTR recombination in my hands.

## Colony PCR:

- 7 Replica plate colonies on a new LB-Amp plate prior to inserting the tip to 10µl water.  
Mix well/ vortex  
Boil at 98°C for 10 mins
- 8 Using primers upstream and downstream of the guide insert sequence, perform colony PCR:  
F primer: gca tat acg ata caa ggc tgt tag aga ga  
R primer: gag cca gta cac gac atc act t

		µl
	10X buffer	2
	10mM dNTPs	0.4

	F p r i m e r ( 1 0 u M )	0.3
	R p r i m e r ( 1 0 u M )	0.3
	B o i l e d c o l o n y	0.5
	T a q p o l	0.1
	d d H 2 O	16.4
	T o t a l:	20

\* Can be scaled up or down.

PCR program:

	1	94°C	5 min
--	---	------	-------



2	94°C	30 sec
3	54°C	1 min
4	72°C	1 min/Kb
5	Go to step 2	X30
6	72°C	5 min
7	4°C	Hold

- 9 Run on a gel -  
product of positive colonies: 500 bp  
product of negative colonies: 2500 bp
- 10 Send ~2 colonies per guide for sequencing with U6\_fwd commercial primer.