

May 22, 2023

Cloning individual AAV capsid variants

Forked from a private protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj87ebgk5/v1>

Sripriya Ravindra Kumar¹, Timothy F. Shay¹, Xinhong Chen¹, David Brown¹, Tatyana Dobрева¹, Qin Huang¹, Xiaozhe Ding¹, Yicheng Luo¹, Pétur . Einarsson¹, Alon Greenbaum¹, Min J. Jang¹, Benjamin E. Deverman¹, Viviana Gradinaru¹

¹California Institute of Technology



Miguel Chuapoco

California Institute of Technology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj87ebgk5/v1>

External link: <https://rdcu.be/c9Vzw>

Protocol Citation: Sripriya Ravindra Kumar, Timothy F. Shay, Xinhong Chen, David Brown, Tatyana Dobрева, Qin Huang, Xiaozhe Ding, Yicheng Luo, Pétur . Einarsson, Alon Greenbaum, Min J. Jang, Benjamin E. Deverman, Viviana Gradinaru 2023. Cloning individual AAV capsid variants. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj87ebgk5/v1>

**Manuscript citation:**

Ravindra Kumar, S., Miles, T.F., Chen, X.*et al.* Multiplexed Cre-dependent selection yields systemic AAVs for targeting distinct brain cell types. *Nat Methods* **17**, 541–550 (2020). <https://doi.org/10.1038/s41592-020-0799-7>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 08, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 81615

Keywords: ASAPCRN, individual aav capsid variants for single variant characterization, individual aav capsid variants this protocol, cloning individual aav, capsid variant, individual aav, aav, capsid, single variant characterization

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020495

Abstract

This protocol describes the procedure to clone individual AAV capsid variants for single variant characterization.

Protocol materials

✂ MscI - 1,250 units **New England Biolabs Catalog #R0534L**

✂ Agel - 1,500 units **New England Biolabs Catalog #R0552L**

✂ pUCmini-iCAP-PHP.B plasmid **addgene Catalog #103002**

✂ Q5 Hot Start High-Fidelity 2X Master Mix - 500 rxns **New England Biolabs Catalog #M0494L**

✂ NEB Stable Competent E.coli (High Efficiency) - 20×0.05 ml **New England Biolabs Catalog #C3040H**



- 1 Digest  pUCmini-iCAP-PHP.B plasmid **addgene Catalog #103002** with
 MscI - 1,250 units **New England Biolabs Catalog #R0534L** and
 AgeI - 1,500 units **New England Biolabs Catalog #R0552L**
- 2 Design 100 bp primers for the desired plasmid variant that cover the insertion region with ~20 bp overlap of the backbone
- 3 To synthesize the double-stranded DNA fragment, anneal the primers by PCR with 20 cycles of  98 °C for  00:00:10 ,  60 °C for  00:00:30 and  72 °C for  00:00:10 using
 Q5 Hot Start High-Fidelity 2X Master Mix - 500 rxns **New England Biolabs Catalog #M0494L**
- 4 Assemble the variant-specific fragment into the digested backbone by Gibson assembly
- 5 Transform the assembled plasmid into
 NEB Stable Competent E.coli (High Efficiency) - 20×0.05 ml **New England Biolabs Catalog #C3040H**
- 6 Select colonies on carbenicillin/ampicillin-LB-agar plates.

50s