Oct 24, 2020

# Cloning of *Bacillus mycoides*

DOI

dx.doi.org/10.17504/protocols.io.bkuwkwxe

Andreea S<sup>1</sup>

<sup>1</sup>University of Groningen

iGEM Groningen 2020

a.stan.6





DOI: dx.doi.org/10.17504/protocols.io.bkuwkwxe

Protocol Citation: Andreea S 2020. Cloning of Bacillus mycoides . protocols.io <u>https://dx.doi.org/10.17504/protocols.io.bkuwkwxe</u>


**License:** This is an open access protocol distributed under the terms of the **<u>Creative Commons Attribution License</u>**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Other Protocol has been developed based on literature, but hasn't been tested yet.

Created: September 04, 2020

Last Modified: October 24, 2020

Protocol Integer ID: 41590



## Abstract

**Introducing NLP14a in the genome of B. mycoides:** In order to create plasmid pYCR-gamyNLP, the backbone of Pycr will be digested with the PCR product of gamy\_Fw and gamy\_Rv to produce the pYCR\_gamy (CRISPR vector containing the gRNA). The NLP14a sequence will be ordered from Twist Biosciences And PCR-amplified with the primers NLP\_Fw and NLP\_Rv. The resulting PCR product, as well as the pYCR\_gamy vector will be digested with Sfil. The digested mixtures will be ligated using T4 ligase.

**Kill switch cloning:** The suggested kill switch mechanism is based on Trp auxotrophic strains as well as a toxinantitoxin mechanism. The toxin will be introduced in the genome using crispr in place of trpE. The antitoxin will be maintained as a cytoplasmic plasmid (pAD-YqcF) and expressed only in the presence of solanine. Unfortunately, to the best of our knowledge, no operator that binds solanine of B. mycoides has been described in literature. If it is provebn that no such operator is present, we suggest on oly keeping the tryptophan dependence strategy.

The plasmid pAD-KPS12-Pman will used as a starting point for introducing the antitoxin gene (ygcF) in *B. mycoides.* The vector will be cut with Xbal and Sphl and ligated with the PCR amplified ygcF (using ygcF\_Fw and ygcF\_Rv as primers). In order to make the antitoxin expression inducible by solanine we plan on replacing the mannose inducible promoter with a promoter induced by solanine. To accomplish this, we will use the primers Psol\_Fw and Psol\_Rv to amplify the solanine promoter. The pAD-ygcFvector will be digested with EcoRI and Xbal, and ligated with the solanine promoter PCR product that will be digested as well with the up mentioned restriction enzymes.

In order to create plasmid pYCR-gtrpEYqcG plasmid, the backbone of Pycr will be digested with Bsal and ligated with the annealing product of gtrpe\_Fw and gtrpe\_Rv to produce the pYCR\_gtrp (CRISPR vector containing the gRNA). Synthetic dna for *yqcG* will be ordered from Twist Biosciences and PCR-amplified with the primers Yqcg\_Fw and Ygcg\_Rv. The resulting PCR product, as well as the pYCR\_gtrpE vector will be digested with Sfil. The digested mixtures will be ligated using T4 ligase.

#### Making electro competent cells

- 1 Pick 1 colony of B. mycoides M2E\_15 and inoculate it in BHIS
- 2 When OD600nm reached 0.85 add 2% glycine and 2% threonine in order to weaken the cell wall.
- 3 Grow the cells overnight at **\$5** 200 rpm, 30°C
- 4 Dilute the overnight culture 50 times in LBSP medium until the OD600nm reaches 0.65
- 5 Collect the cells by 🚯 4000 x g, 4°C, 00:10:00 . Discard supernatant
- Wash cell pellet with pre-chilled electroporation buffer (10% glycerol, 0.25 M sorbitol, 0.25M trehalose) x4
- 7 Suspended in 1 ml electroporation buffer (10% glycerol, 0.25 M sorbitol, 0.25M trehalose).
- 8 Snap-freeze the electrocompetent cells in liquid nitrogen and sore at 📲 -80 °C

### **PCR - Phusion cloning**

- 9 For the PCR mix pipette:
  - $45 \text{ }_{\mu\text{L}}$  primers mixture Fw + Rv ( 5 uM) (see table "Primers")
  - 👗 10 µL HF buffer
  - 🗕 1 μL dNTP (10mM)
  - DNA template (0.1 ng)
  - Δ 0.25 μL phusion polymerase
  - MIIiQ to 🛛 🕹 50 µL
- 10 PCR reaction is performed with the following protocol:

Temperature (*C)	Time (mm: ss)
98	05:00
98	10:00
Variable	00:30
 72	00:30
repeat above steps	30x
72	10:00
11	Infinit e

11 Load  $\underline{4} 5 \mu \underline{L}$  of sample plus  $\underline{4} 1 \mu \underline{L}$  of staining solution 6x. Check the length of the band (bp).

Purify sample : PCR cleanup/ gel extraction according to manufacturer.

## Cloning

- 12 Digest the polymerized insert and the desired vector with 2U of apropriate restriction enzyme. See table "Primers" in order to choose the enzyme.
- 13 Digestion mixture (  $\underline{A}$  20  $\mu$ L ):

-Sfil 2 U

- Plasmid ( 👗 50 ng )
- -Insert ( 👗 20 ng )
- $\Delta_{2 \mu L}$  CutSmart 10x (or appropriate buffer)
- -MiliQ to 🛛 🕹 20 µL
- 14 Incubate 🕑 01:00:00 at the temperature recommended by the manufacturer of the restriction enzyme

15	Inactivate the restriction enzyme by incubationg 20 min	8	65 °C	
	**Sfil can't be heat inactivated so in order to inactivate it a PCR clean up is necesary			

<u>Ligate</u> the digested insert and vector.
Ligation mixture :
-10 ul digestion product
-1 ul T4
-2 ul T4 ligase buffer
-7 ul MQ

17 Incubate 🕥 02:00:00

**B** Room temperature

## Transformation (electroporation)

18 Thaw on ice  $4 100 \,\mu\text{L}$  of electrocompetent cells

- 19 Add  $\underline{\underline{A}}_{2 \mu g}$  of plasmid and the aliquot of electrocompetent cells to ice-cold electroporation cuvettes
- 20 Electroporate 25 uF, 10 kV/cm, 200 Ohms
- Add 1 ml of BHIS and incubate (5 100 rpm, 30°C, 05:00:00 for recovery
- Plate on LB + Agar (1.5%) + 100 μg/ml spectinomycin.For crispr cloning, to activate cas9 expression, add 0.2% mannose.

#### **Mutant selection**

- 23 Incubate at 🖁 30 °C 🚫 Overnight
- Randomly pick potential mutants and test by colony PCR (for primers use the ones that flank the region of insertion sites).

25 After selecting colonies that show the expected band size, purify the desired plasmid (miniprep) and transform into B. mycoides.