


Sep 25, 2019

🌐 Cloning of L0 parts into pUAP-ye for Loop typellS

DOI

dx.doi.org/10.17504/protocols.io.4wdgxa6



Eftychis Frangedakis¹, Susana Sauret-Gueto²

¹University of Cambridge; ²Plant Sciences, University of Cambridge, OpenPlant

OpenPlant Project



Eftychis Frangedakis

University of Cambridge, Plant Sciences

OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4wdgxa6

Protocol Citation: Eftychis Frangedakis, Susana Sauret-Gueto 2019. Cloning of L0 parts into pUAP-ye for Loop typellS. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4wdgxa6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 29, 2019

Last Modified: September 25, 2019

Protocol Integer ID: 25253

Keywords: Loop assembly Type IIS Level 0 marchantia



- 1 Design primers as described in the Figure below to add the necessary overhangs for cloning into pUAP-pe and also the overhangs that correspond to the common syntax.
- 2 DNA parts are PCR amplified from the source DNA part (e.g. plasmid DNA or genomic DNA), using a high fidelity DNA polymerase such as Phusion. Use 10 ng of DNA if the template is plasmid DNA and 100ng if the template is genomic DNA. The cycling conditions are: Denaturation at 98 °C for 30 s. 35 cycles of: denaturation at 98 °C for 10 s, annealing at primer annealing temperature for 30 s, and extension at 72 °C for 15 s/kb. Final extension at 72 °C for 10 min.
- 3 Run the PCR products on a 1.5% (w/v) agarose gel.
- 4 Gel extract the band that corresponds to the size of the amplified DNA part using a kit such as QIAquick Gel Extraction Kit.
- 5 Determine DNA concentration with spectrophotometry (Nanodrop).
- 6 Prepare aliquots of the DNA part at a concentration of 15 nM and of the pUAP-pe acceptor vector at a concentration of 7.5 nM.
- 7 Set up a Type IIS assembly reaction into a 0.2 mL tube according to Table to clone the amplified DNA part into pUAP-pe.

Sterile water	5
pUAP-pe	1
DNA part	1



10x Tang o buffer (Ther mo Fishe r)	1
1 mg/m L bovin e seru m albu min (NEB)	0.5
T4 DNA ligase (5 U/ μL) (Ther mo Fishe r)	0.25
10mM ATP (SIG MA)	1
Sapl (Lgul) (5 U/ μL) (Ther mo Fishe r)	0.25
Final volum e	10

- 8 Place samples on a thermocycler and use the following program:
Assembly: 26 cycles of 37 °C for 3 min and 16 °C for 4 min.
Termination and enzyme denaturation: 50 °C for 5 min and 80 °C for 10 min.
- 9 Transform 20 μL of chemically competent E. coli cells (transformation efficiency of 1×10^7 transformants/μg plasmid DNA) using 2 μL of the assembly reaction and then spread



on LB agar plates containing 25 µg/mL chloramphenicol.

- 10 Incubate overnight at 37 °C.
- 11 Select white colonies for sequencing.
- 12 Confirm the presence of the correct insert with Sanger sequencing using the primers UAP_F and UAP_R and any additional DNA part specific primers.
- 13

A**B****C**