



Feb 06, 2024

🌐 Cloning, Protein Expression, and Purification of 20S CPs and Assembly Intermediates



In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldmd59l5b/v1>


Frank Adolf¹

¹Department of Molecular Machines and Signaling, Max Planck Institute of Biochemistry, 82152 Martinsried, Germany



Frank Adolf

Department of Molecular Machines and Signaling, Max Planck I...



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldmd59l5b/v1>

Protocol Citation: Frank Adolf 2024. Cloning, Protein Expression, and Purification of 20S CPs and Assembly Intermediates. protocols.io <https://dx.doi.org/10.17504/protocols.io.n92ldmd59l5b/v1>

Manuscript citation:

<https://www.biorxiv.org/content/10.1101/2024.01.27.577538v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 05, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94683

Keywords: ASAPCRN, proteasome, core particle, 20S proteasome, chaperone, molecular machine, multiprotein complex, POMP, PAC1, PAC2, PAC3, PAC4, propeptide, protease, protein expression, methods for cloning, protocol details methods for cloning, purification, cloning, assembly intermediate, protein, cp

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000282

Abstract

This protocol details methods for cloning, expression, and purification of 20S CPs and assembly intermediates for biochemical and structural analysis.

Guidelines

Please familiarise yourself with the laboratory safety rules and guidelines and follow these while performing the experiment. Please wear appropriate PE while performing the experiment.

Materials

X-tremeGENETM HP DNA Transfection Reagent - ROCHE

cOmplete PROTEASE INHIBITOR COCKTAIL - ROCHE

Strep-Tactin[®] Sepharose[®] resin - IBA

Safety warnings

- ⚠ Please refer to the Safety Data Sheets (SDS) for health and environmental hazards. Liquid nitrogen (LN2) and other cryogenics can cause severe damage to the skin and eyes. Always wear personal protective equipment when handling these cryogenics.



Cloning of baculovirus transfer vectors

1w



- 1 Baculovirus transfer vectors were assembled utilizing a combination of the biGBac and MultiBac systems

Citation

Weissmann F, Petzold G, VanderLinden R, Huis In 't Veld PJ, Brown NG, Lampert F, Westermann S, Stark H, Schulman BA, Peters JM (2016)

. biGBac enables rapid gene assembly for the expression of large multisubunit protein complexes.

<https://doi.org/10.1073/pnas.1604935113>

LINK

Citation

Vijayachandran LS, Viola C, Garzoni F, Trowitzsch S, Bieniossek C, Chaillet M, Schaffitzel C, Busso D, Romier C, Poterszman A, Richmond TJ, Berger I (2011)

. Robots, pipelines, polypeptides: enabling multiprotein expression in prokaryotic and eukaryotic cells.

<https://doi.org/10.1016/j.jsb.2011.03.007>

LINK

The MultiBac vector pACEBac1 was used as both library vector and acceptor vector for step 1 biGBac assembly reactions with the following primers:

ACEBac-1EC-BBA.rev

ATTTAAATCTTTAGACCATAGAGCGTTCTCGCGAATCGATACTAGTGTTTAAACTCGCTACCT
TAGGACC

ACEBac-2EC-BBA.fwd

ATTTAAATAAACCTAATGATGCCTGATGTTTCCTAGGGTATACCCATCTAATTGGAACCAGAT
AAGTGAAATC



ACEBac-3EC-BBA.fwd

ATTTAAATAAACGGTTCACATAGCTTAGTTTCCTAGGGTATACCCATCTAATTGGAACCAGAT
AAGTGAAATC

ACEBac-4EC-BBA.fwd

ATTTAAATAAACACTGACATTGACTTGTTTCCTAGGGTATACCCATCTAATTGGAACCAGAT
AAGTGAAATC

ACEBac-5EC-BBA.fwd

ATTTAAATAAATCTATATCTCAATCGGGGTTCTAGGGTATACCCATCTAATTGGAACCAGAT
AAGTGAAATC

- 2 Clone all 20S CP subunits and assembly chaperones into pACEBac1 using Gibson assembly

1d



For affinity purification add C-terminal TEV-cleavable twin strep tags on β 2 (PSMB7) and β 7 (PSMB4)

- 2.1 Screen for positive clones and sequence verify by Sanger sequencing

1d



- 3 Preparation of bigBac assembly inserts:

6h

Amplify expression cassette from library vectors (step 2) by PCR with the following primers:



Cas1-ACEBac.fwd

AACGCTCTATGGTCTAAAGATTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas2-ACEBac.fwd

AAACTGGATACTATTGCACGTTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas3-ACEBac.fwd

AAACCTAATGATGCCTGATGTTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas4-ACEBac.fwd

AAACGGTTCACATAGCTTAGTTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas5-ACEBac.fwd

AAACACTGACATTGACTTGTTTAAATCGACCTACTCCGGAATATTAATAGATCATGG



Cas1-ACEBac.rev

AAACGTGCAATAGTATCCAGTTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

Cas2-ACEBac.rev

AAACATCAGGCATCATTAGGTTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

Cas3-ACEBac.rev

AAACTAAGCTATGTGAACCGTTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

Cas4-ACEBac.rev

AAACCAAGTCAATGTCAGTGTTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

Cas5-ACEBac.rev

AACCCCGATTGAGATATAGATTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

3.1 Purify all inserts by agarose gel extraction

2h



4 Preparation of biGBac step 1 assembly acceptor vectors: linearize the acceptor vector pACEBac1 by PCR with the following primers:

6h



for cloning of 5 expression cassettes:

ACEBac-1EC-BBA.rev and ACEBac-5EC-BBA.fwd

for cloning of 4 expression cassettes

ACEBac-1EC-BBA.rev and ACEBac-4EC-BBA.fwd

for cloning of 3 expression cassettes

ACEBac1-1EC-BBA.rev and ACEBac-3EC-BBA.fwd

4.1 Purify all linearized vectors by agarose gel extraction

2h



5 Clone multi expression cassette Baculo transfer vectors listed below by biGBac step 1 assembly with expression cassette inserts from step 3 and acceptor vectors from step 4

1d



pACEBac1-POMP-PSMG1-PSMG2-PSMG3-PSMG4

pACEBac1-PSMA1-PSMA2-PSMA3-PSMA4,

pACEBac1-PSMA5-PSMA6-PSMA7,

pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4,

pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4-TEV-2xSTII,

pACEBac1-PSMB5-PSMB6-PSMB7,

pACEBac1-PSMB5-PSMB6-PSMB7-TEV-2xSTII.

- 5.1 Screen for positive clones and sequence verify by Sanger sequencing

1d



- 6 Assemble finale multi expression cassette baculo transfer vectors by multibac assembly

2d



pACEBac1-PSMA1-PSMA2-PSMA3-PSMA4-PSMA5-PSMA6-PSMA7

pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4-TEV-2xSTII-PSMB5-PSMB6-PSMB7

pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4-PSMB5-PSMB6-PSMB7-TEV-2xSTII

- 6.1 Digest acceptor vectors listed below with I-CeuI, CIP and purify by agarose gel extraction

4h



pACEBac1-PSMA5-PSMA6-PSMA7

pACEBac1-PSMB5-PSMB6-PSMB7

pACEBac1-PSMB5-PSMB6-PSMB7-TEV-2xSTII

- 6.2 Digest donor vector listed below with I-CeuI and BstXI, and purify inserts by agarose gel extraction

4h



pACEBac1-PSMA1-PSMA2-PSMA3-PSMA4

pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4-TEV-2xSTII

pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4

- 6.3 Ligate the following vector/insert pairs listed below and transform in DH5alpha *E. coli*

1d



pACEBac1-PSMA5-PSMA6-PSMA7 and EC1-4 PSMA1-PSMA2-PSMA3-PSMA4

pACEBac1-PSMB5-PSMB6-PSMB7 and EC1-4 PSMB1-PSMB2-PSMB3-PSMB4-TEV-2xSTII

pACEBac1-PSMB5-PSMB6-PSMB7 TEV-2xSTII and EC1-4 PSMB1-PSMB2-PSMB3-PSMB4

- 6.4 Screen for positive clones and sequence verify by Sanger sequencing

1d



Baculo virus amplification and insect cell expression

1w

- 7 Culture Sf9 insect cells (Thermo Fisher Scientific) for virus amplification in serum-free Ex-cell 420 medium (Sigma-Aldrich)



Culture Trichoplusia ni (Thermo Fisher Scientific) in protein free ESF 921 insect cell culture media (Expression Systems LLC)

All steps were carried out according to standard protocols



Citation

Fitzgerald DJ, Berger P, Schaffitzel C, Yamada K, Richmond TJ, Berger I (2006)
. Protein complex expression by using multigene baculoviral vectors.

<https://doi.org/>

LINK

Citation

Bieniossek C, Richmond TJ, Berger I (2008)
. MultiBac: multigene baculovirus-based eukaryotic protein complex production.

<https://doi.org/10.1002/0471140864.ps0520s51>

LINK

- 8 For bacmid preparation transform baculo transfer vector from the section above into EMBACY *E. coli* and plate on LB-agar plates with ampicillin, kanamycin, tetracyclin, gentamycin, IPTG, and XGal
- 8.1 Prepare  Overnight cultures from white colonies in LB with ampicillin, kanamycin, and gentamycin
- 8.2 Prepare bacmids by alkaline lysis and subsequent isopropanol and  70 % (v/v) EtOH precipitation
- 8.3 Air dry bacmid pellets, resuspend DNA in milliQ water, and store at  4 °C until usage
- 9 P1 virus production

1d



1d



2h



1h



2d 12h





- 9.1 For P1 virus production seed $0.7 - 0.8 \times 10^6$ cells/well in  3 mL medium in a six well plate and leave for  00:15:00 min at  27 °C  15m
- 9.2 Prepare bacmids by diluting  1 µg bacmid DNA in  20 µL milliQ water, and add  200 µL of medium  5m
- 9.3 For each bacmid mix  100 µL medium with  15 µL Extreme Gene HP DNA transfection reagent (Roche)  5m
- 9.4 Add  112 µL of Medium Extreme Gene Mix to each bacmid DNA  5m
- 9.5 Add  156 µL of the DNA-Extreme Gene Mix to each well in the six-well plate  5m
- 9.6 Wrap the plates with Parafilm and incubate for  60:00:00 h at  27 °C  2d 12h
- 9.7 Transfer medium containing P1 virus from each well into a  15 mL tube and store at  4 °C until further usage   10m
- 10 P2 and P3 virus amplification   2d
- 10.1 Seed $0.8 - 1.0 \times 10^6$ Sf9 cells/ml into an appropriate conical flask and infect cells with 0.5 - 1.0% P1 or P2, respectively  1h
- 10.2 Incubate cell suspensions at  27 °C at 80 rpm for  48:00:00 h  2d
- 10.3 Harvest P2 and P3 by centrifugation, transfer medium containing virus into a  50 mL tube or  250 mL bottle and store at  4 °C until further usage   15m



11 Protein Expression in High FiveTM insect cells

2d 12h



11.1 Seed 2.0×10^6 High FiveTM cells/ml into an appropriate conical flask and infect with 0.5 - 1.0 % P3 virus

30m



11.2 Incubate cell suspensions at $27\text{ }^{\circ}\text{C}$ at 80 rpm for 60:00:00 h

2d 12h



11.3 Harvest High FiveTM cells by centrifugation, resuspend cell pellets in PBS, transfer to a 50 mL tube, harvest cells by centrifugation, discard supernatant, snap freeze cell pellets in LN₂, and store until further usage at $-80\text{ }^{\circ}\text{C}$

30m



Purification of 20S CPs and 20S CP assembly intermediates

8h

12 Purification of mature 20S CPs together with their assembly intermediates

8h



12.1 Thaw cell pellets from section 2 - step 11, resuspend in buffer A 25 mM HEPES pH 7.5 (KOH), 150 millimolar (mM) NaCl, 1 millimolar (mM) DTT, and add 1 tablet cOmplete protease inhibitor mix (Merck) per 10 mL cell pellet resuspended in 40 mL buffer A in a 50 ml tube

30m



12.2 Lyse cells by sonication on ice

10m



12.3 Centrifuge lysate at 22k at $4\text{ }^{\circ}\text{C}$ for 01:00:00 h

1h



12.4 In parallel wash Strep-Tactin[®] Sepharose[®] resin (IBA) 3 times with buffer A

45m



12.5 Incubate supernatants from step 12.3 with resin for 01:00:00 h at $4\text{ }^{\circ}\text{C}$

1h



12.6 Wash resin 3 times with buffer A by centrifugation at 2000 x g for 00:15:00 min at

4 °C

15m



12.7 Load resin on a gravity flow column and elute proteins with 2.5 micromolar (μ M) *d*-Desthiobiotin in buffer A

15m

12.8 Pool fractions of interest, concentrate in a spin protein concentrator with a MWCO of 100 kDa to max 5-8 mg/mL

1h

12.9 In parallel equilibrate a Superose 6 10/300 GL column (Cytiva) with buffer A

1h

12.10 Apply concentrated eluate from step 12.8 onto the SEC column at a flow rate of 1 ml/min, record absorption at 280 nm, and fractionate elution a 300 μ L

1h

12.11 Analyze SEC fractions by SDS-PAGE, pool fractions of interest, snap freeze in LN2, and store until further usage at -80 °C for biochemical assays

2h



For structural analysis use fractions directly after SEC

Citations

Step 7

Fitzgerald DJ, Berger P, Schaffitzel C, Yamada K, Richmond TJ, Berger I. Protein complex expression by using multigene baculoviral vectors.

<https://doi.org/>

Step 7

Bieniossek C, Richmond TJ, Berger I. MultiBac: multigene baculovirus-based eukaryotic protein complex production.

<https://doi.org/10.1002/0471140864.ps0520s51>