


Dec 31, 2019 Version 2

Cloning shRNA Oligos into pLKO.1 V.2

DOI

dx.doi.org/10.17504/protocols.io.bawmifc6



Addgene The Nonprofit Plasmid Repository¹

¹Addgene



Addgene The Nonprofit Plasmid Repository

Addgene

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bawmifc6>

External link: <http://www.addgene.org/tools/protocols/plko/>

Protocol Citation: Addgene The Nonprofit Plasmid Repository 2019. Cloning shRNA Oligos into pLKO.1. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bawmifc6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 31, 2019

Last Modified: December 31, 2019

Protocol Integer ID: 31405

Keywords: cloning shrna oligo, designing shrna oligo, trc cloning vector, protocol

Abstract

This is the protocol accompanying the "pLKO.1 – TRC Cloning Vector". For information about the PLKO.1-TRC cloning vector and tips on designing shRNA oligos for pLKO.1 see Addgene's website:

<http://www.addgene.org/tools/protocols/plko/>

Materials

STEP MATERIALS

☒ NEBuffer 3 - 5.0 ml **New England Biolabs Catalog #B7003S**

☒ NEBuffer 1 - 5.0 ml **New England Biolabs Catalog #B7001S**

☒ Agel - 300 units **New England Biolabs Catalog #R0552S**

☒ EcoRI - 10,000 units **New England Biolabs Catalog #R0101S**

☒ EcoRI - 10,000 units **New England Biolabs Catalog #R0101S**

☒ T4 DNA Ligase Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0202S**

☒ T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**



Protocol materials

☒ NEBuffer 3 - 5.0 ml **New England Biolabs Catalog #B7003S**

☒ NEBuffer 1 - 5.0 ml **New England Biolabs Catalog #B7001S**

☒ Agel - 300 units **New England Biolabs Catalog #R0552S**

☒ EcoRI - 10,000 units **New England Biolabs Catalog #R0101S**

☒ EcoRI - 10,000 units **New England Biolabs Catalog #R0101S**

☒ T4 DNA Ligase Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0202S**

☒ T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

☒ NEBuffer 3 - 5.0 ml **New England Biolabs Catalog #B7003S**

☒ NEBuffer 1 - 5.0 ml **New England Biolabs Catalog #B7001S**

☒ Agel - 300 units **New England Biolabs Catalog #R0552S**

☒ EcoRI - 10,000 units **New England Biolabs Catalog #R0101S**

☒ EcoRI - 10,000 units **New England Biolabs Catalog #R0101S**









☒ T4 DNA Ligase Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0202S**

☒ T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

Troubleshooting



Annealing Oligos

- 1 Resuspend oligos in ddH₂O to a concentration of 20 μ M.
- 2 Add 5 μ L Forward oligo
 5 μ L
- 3 Add 5 μ L Reverse oligo
 5 μ L
- 4 Add 5 μ L 10x NEB buffer 2
 5 μ L
 NEBuffer 3 - 5.0 ml **New England Biolabs Catalog #B7003S**
- 5 Add 35 μ L ddH₂O
 35 μ L
- 6 Incubate for 4 minutes at 95°C in a PCR machine or in a beaker of boiling water.
 00:04:00
- 7 Incubate the sample at 70°C for 10 minutes in a PCR machine.
 00:10:00
- 8 Slowly cool to room temperature over the period of several hours.
 03:00:00

Note

This will take a few hours, but it is important for the cooling to occur slowly for the oligos to anneal.

Note

If using a beaker of water, remove the beaker from the flame, and allow the water to cool to room temperature.

Digesting pLKO.1 TRC Cloning Vector

- 9 Mix: 6 μ g pLKO.1 TRC-cloning vector (maxiprep or miniprep DNA)



 6 µg

10 with 5 µL 10x NEB buffer 1

 5 µL

 NEBuffer 1 - 5.0 ml **New England Biolabs Catalog #B7001S**

11 with 1 µL Agel


 1 µL

 Agel - 300 units **New England Biolabs Catalog #R0552S**

12 bring to 50 µL ddH₂O

 50 µL

13 Incubate at 37°C for 2 hours.


 02:00:00

14 Purify with Qiaquick gel extraction kit, eluting in 30 µL of ddH₂O.

15 Digest eluate with EcoRI by mixing: 30 µL pLKO.1 TRC-cloning vector digested with Agel


16 with 5 µL 10x NEB buffer for EcoRI

 5 µL

 EcoRI - 10,000 units **New England Biolabs Catalog #R0101S**

17 with 1 µL EcoRI


 1 µL

 EcoRI - 10,000 units **New England Biolabs Catalog #R0101S**

18 and 14 µL ddH₂O

 14 µL

19 Incubate at 37°C for 2 hours.

 02:00:00

20 Run digested DNA on 0.8% low melting point agarose gel until you can distinctly see 2 bands, one 7kb and one 1.9kb.


**Note**

When visualizing DNA fragments to be used for ligation, use only long-wavelength UV light. Short wavelength UV light will increase the chance of damaging the DNA.

- 21 Cut out the 7kb band and place in a sterile microcentrifuge tube.
- 22 Purify the DNA using a Qiaquick gel extraction kit. Elute in 30 μL of ddH₂O.
- 23 Measure the DNA concentration.

Ligating and Transforming into Bacteria

- 24 Use your ligation method of choice. For a standard T4 ligation, mix: 2 μL annealed oligo from "Annealing Oligos" section above.
- 25 With 20 ng digested pLKO.1 TRC-cloning vector from the "Digesting pLKO.1 TRC Cloning Vector" section above.

 2 μL


 20 ng

Note

If you were unable to measure the DNA concentration, use 1 μL

- 26 With 2 μL 10x NEB T4 DNA ligase buffer

 2 μL

 T4 DNA Ligase Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0202S**

- 27 With 1 μL NEB T4 DNA ligase


 1 μL

 T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

- 28 Bring up to 20ul with ddH₂O

- 29 Incubate at 16°C for 4-20 hours.



 04:00:00

- 30 Transform 2 μL of ligation mix into 25 μL competent DH5 alpha cells, following manufacturer's protocol.
- 31 Plate on LB agar plates containing 100 $\mu\text{g/mL}$ ampicillin or carbenicillin (an ampicillin analog).