



Mar 27, 2025

Clusterin Phospholipid Particles Flotation Assay Using an Optiprep Step Gradient

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vz6mq4gx1/v1>

Patricia Yuste-Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vz6mq4gx1/v1>

Protocol Citation: Patricia Yuste-Checa, F Ulrich Hartl 2025. Clusterin Phospholipid Particles Flotation Assay Using an Optiprep Step Gradient. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vz6mq4gx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: March 20, 2025

Last Modified: March 27, 2025

Protocol Integer ID: 124798

Keywords: ASAPCRN, clusterin phospholipid particles flotation assay, separate clusterin phospholipid particle, free clusterin, clusterin, optiprep step gradient this protocol, optiprep step density gradient, optiprep step gradient, binders by flotation, flotation, dmpc

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol is based on Yeh FL et al. Neuron (2016) and details how to efficiently separate Clusterin phospholipid particles (Clusterin:DMPC) from free Clusterin together with binders by flotation using an Optiprep step density gradient.

Materials

Reagents

- PBS
- 2% Na-deoxycholate
- Trichloroacetic acid
- acetone
- OptiPrep Density Gradient Medium (Merck, D1556)
- 3.5 mL Thickwall Polycarbonate Tube (349622, Beckman Coulter)
- NuPAGE™ LDS Sample Buffer (4X) buffer (Thermo Fisher Scientific)

Equipments

SW55 Ti rotor (Beckman Coulter)

Bioruptor sonication bath (Diagenode)

 OptiPrep™ Density Gradient Medium **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1556**

 3.5 mL, Open-Top Thickwall Polycarbonate Tube, 13 × 51mm - 25Pk **Beckman Coulter Catalog #349622**



Rhodanese Aggregation Assay

- 1 Perform the rhodanese aggregation reaction* in the presence of Clusterin or purified Clu:DMPC nanodiscs** at molar ratio D-Rhod:Clu 1:3.

Note

* See protocol 'Prevention of rhodanese aggregation assay'.

** See 'Formation and isolation of Clu-phospholipid particles'.
[dx.doi.org/10.17504/protocols.io.bp2l6x59zlqe/v1](https://doi.org/10.17504/protocols.io.bp2l6x59zlqe/v1)

- 2 Centrifuge the rhodanese aggregation reactions in the presence of Clusterin or Clu:DMPC at  20000 x g, 4°C, 00:15:00 to remove big aggregates.

15m



Note

NOTE: Here we use D-Rhodanese, but binding of other Clusterin substrates to Clusterin:DMPC nanodiscs can be tested and isolated using this assay.

Optiprep Gradient

- 3 Mix the supernatant with OptiPrep Density Gradient Medium (Merck, D1556, 60% iodixanol) to a final concentration of 36% iodixanol in PBS and final volume of  1 mL in a 3.5 mL Thickwall Polycarbonate Tube (349622, Beckman Coulter).
- 4 Carefully place  1 mL of 24% iodixanol diluted in PBS on top of the 36% iodixanol layer.
- 5 Carefully place  1 mL of PBS as top layer.
- 6 Centrifuge the gradient at  54000 rpm, 4°C, 03:00:00 using a SW55 Ti rotor (Beckman Coulter).



3h





7 After centrifugation, carefully collect 6 fractions of  0.5 mL .

TCA Precipitation for SDS-PAGE Analysis

8 Dilute each fraction with  0.5 mL PBS.

9 Add  40 μ L of 2% Na-deoxycholate and vortex. 

10 Incubate  On ice for  00:15:00 .

15m



11 Add  100 μ L of 100% Trichloroacetic acid and vortex. 

12 Incubate  On ice for  01:00:00 .

1h



13 Centrifuge the samples at  20000 x g, 4°C, 00:30:00 .

30m



14 Add  500 μ L of ice-cold acetone. 

15 Vortex or sonicate in a Bioruptor sonication bath (Diagenode) (2 cycles of  00:00:30 on –  00:00:30 off) for faster resuspension.

16 Centrifuge at  200000 x g, 4°C, 00:10:00 .

10m



17 Air dry the pellets.

18 Add  30 μ L  300 millimolar (mM) Tris  8.8 . 



19 Incubate  On ice for  00:05:00 .

5m



20 Add  30 μ L NuPAGE™ LDS Sample Buffer (4X) buffer (Thermo Fisher Scientific) containing  100 millimolar (mM) DTT and boil the samples for  00:10:00 .

10m



21 Separate the proteins by sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE) followed by immunoblotting.

Note

NOTE: Isolated Clusterin phospholipid particles bound to a substrate can be also used for further experiments like cellular uptake, see Yeh FL *et al. Neuron* (2016).

Protocol references

Yeh FL, Wang Y, Tom I, Gonzalez LC, Sheng M. TREM2 Binds to Apolipoproteins, Including APOE and CLU/APOJ, and Thereby Facilitates Uptake of Amyloid-Beta by Microglia. *Neuron*. 2016 Jul 20;91(2):328-40. doi: 10.1016/j.neuron.2016.06.015. PMID: 27477018.