


Mar 27, 2024

CMT-93 Cell Culture Protocol

DOI

dx.doi.org/10.17504/protocols.io.bp2l62k1dgqe/v1



Laura Gómez¹

¹CBM Severo Ochoa



Laura Gómez

CBM Severo Ochoa

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l62k1dgqe/v1>

External link: <https://www.atcc.org/products/ccl-223>

Protocol Citation: Laura Gómez 2024. CMT-93 Cell Culture Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l62k1dgqe/v1>

Manuscript citation:



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: March 18, 2024

Last Modified: March 27, 2024

Protocol Integer ID: 96872

Keywords: CMT-93, cell line, epithelial cell line, polyploid carcinoma, cell culture, epithelial morphology, exhibiting epithelial morphology, cell line, rectum, cell

Abstract

CMT-93 is a cell line exhibiting epithelial morphology that was isolated from the rectum of a mouse with polyploid carcinoma.

Image Attribution

ATCC - <https://www.atcc.org/products/ccl-223>

Guidelines

Working with cell cultures requires a laminar flow cabinet. It has to be radiated with UV light, cleaned with any highly effective terminal disinfectant (such as Tego® 2000 or Suredis®) and 70% ethanol. All material introduced into the cabinet must also be sprayed with ethanol.

Once the work is finished, we must clean the cabinet with the detergent, then with 70% etOH and turn on the UV light for 30 min.



Materials

Plasticware:

p60 cell culture plates

p100 cell culture plates

Cell culture flasks, 75 cm², treated for cell attachment.

15 and 50 mL centrifuge tubes

Cryovials

To prepare the complete medium:

DMEM 1X

Fetal bovine serum heat inactivated (FBS)

Glutamine 200 mM

To subculture the cells:

Trypsin-EDTA

PBS 1X

Troubleshooting

Safety warnings







Every reagent must be sterile in order to avoid contaminations.

Before start







Clean and prepare the laminar flow cabinet, turn on the water bath and warm up the culture media.



Preparation of complete growth medium (DMEM+)

- 1 Add  445 mL 1X DMEM,  50 mL FBS and  5 mL glutamine [M] 200 millimolar (mM) to a sterile 500 mL bottle and homogenize
- 2 Label the bottle with name, group, phone number, date and additions.
- 3 Close with parafilm and store at  4 °C .




Cell thawing procedure

- 4 Remove one vial of cell stock from the liquid nitrogen tank with gloves and forceps. Transfer them to the cell culture laboratory in an appropriate container or a box with ice.
- 5 Thaw the vial by gently shaking it in a  37 °C water bath. Thawing should be rapid (approximately 2 min).
- 6 Remove the vial from the water bath as soon as the contents are thawed, and decontaminate by dipping in or spraying with 70% ethanol
- 7 Transfer the contents of the vial to a centrifuge tube containing  9 mL of complete culture medium and  1200 rpm, Room temperature, 00:05:00 5m
- 8 Resuspend with  10 mL DMEM+ and distribute on 2 P60 plates.
- 9 Incubate cultures at  37 °C , 5% CO₂  Overnight 5m

Subculturing procedure

- 10 Remove and discard culture medium.



- 11 Rinse with PBS 1X solution and discard
- 12 Briefly rinse the cell layer with 0.25% (w/v) Trypsin-0.53mM EDTA solution to remove all traces of serum which contains trypsin inhibitor. Discard.
- 13 Add  2 mL of Trypsin-EDTA solution to flask and incubate  00:10:00 at  37 °C to facilitate detachment from the plate.



10m



Note

To avoid clumping, **do not agitate** the cells by hitting or shaking the flask while waiting for the cells to detach.

- 14 Observe cells under an inverted microscope until cell layer is dispersed.

Note

If the cells are not detached already, incubate  00:05:00 more at  37 °C .

- 15 Add  5 mL of DMEM+ and aspirate cells by **gently** pipetting. Pour the existing volume down the walls of the flask in order to drag and collect as many cells as possible.
- 16 Collect the cell suspension in a centrifuge tube and  1200 rpm, Room temperature, 00:05:00 .
- 17 Discard the supernatant into a beaker with 70% EtOH or 10% bleach.
- 18 Resuspend in medium according to the dilution to be made.





5m

Note

A subcultivation ratio of 1:4 to 1:10 is recommended



19 Add  1 mL of the cell suspension to new culture vessels containing  14 mL DMEM+.


20 Incubate cultures at  37 °C , 5% CO₂  Overnight

5m


Cryopreservation and storage procedure

1d



21 Repeat steps of *Subculturing procedure* until the "Resuspend in medium according to the dilution to be made" step.

22 Resuspend in medium taking into account that for every p100 we can storage up to 2 cryovials of cells, containing  1 mL .

23 Prepare the cryovials with  50 µL DMSO.

24 Add  950 µL of the cell suspension to every cryovial.

25 Label the cryovials with cell line, passage, date and lab number or phone number.

26 Store the cryovials in a slow freezing container at  -80 °C for  24:00:00 .

1d

27 Transfer the cryovials to the liquid nitrogen tank.

Protocol references

<https://www.atcc.org/products/ccl-223#detailed-product-information>