



Mar 12, 2025

Co-crystallization of CLK3 with Ligand

 [Science Advances](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk8zj5l5r/v1>

Nicolai D. Raig¹, Andreas Kramer¹, Stefan Knapp¹

¹Structural Genomics Consortium (SGC)



Nicolai D. Raig

Structural Genomics Consortium (SGC)

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk8zj5l5r/v1>

External link: <https://doi.org/10.1126/sciadv.adt2050>

Protocol Citation: Nicolai D. Raig, Andreas Kramer, Stefan Knapp 2025. Co-crystallization of CLK3 with Ligand. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk8zj5l5r/v1>

Manuscript citation:

Nicolai D Raig, Katherine J Surridge, Marta Sanz-Murillo, Verena Dederer, Andreas Krämer, Martin P Schwalm, Nicholas M Lattal, Lewis Elson, Deep Chatterjee, Sebastian Mathea, Thomas Hanke, Andres E Leschziner, Samara L Reck-Peterson, Stefan Knapp (2025) Type II kinase inhibitors that target Parkinson's disease-associated LRRK2. Science Advances doi: [10.1126/sciadv.adt2050](https://doi.org/10.1126/sciadv.adt2050)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 15, 2024

Last Modified: March 12, 2025

Protocol Integer ID: 109897

Keywords: clk3, ligand

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: 000519

Abstract

This is a step by step protocol, how to set up co-crystallization plates of CLK3 and compound, like it has been done for the PDB submission 9EZ3



Co-Crystallization of CLK3

1

1. Incubate CLK3 protein kinase domain (~11 mg/mL) with compound of choice. 0.5-1 mM of compound if soluble
2. Incubate the mixture on ice for 30 minutes.
3. Prepare the crystallization firescreen around this condition: 15% PEG 6000, 0.1 M Bicine buffer pH 8.0 using 3 lenses 96 well plate
4. Pipette 200-400 nL drops of the protein-reservoir mixture onto a sitting-drop vapor diffusion plate with robot or by hand (use 24 well plate), use protein:reservoir solution ratio 2:1 / 1:1 / 1:2
5. seal plate and incubate at 4°C and incubate at 4°C
6. Monitor crystal growth over time.
7. Before fishing the crystals use Approx. 20-25% Ethylene Glycol as cryo protection
8. Flash Freeze the crystal in Liquid N₂

Protocol references

<https://www.rcsb.org/structure/9EZ3>