

Jul 16, 2025

## Co-cultivation Protocol for *Porosira glacialis* Diatom and Bacterial Isolates

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5ow14g1b/v1>



William Arnli<sup>1</sup>, Nerea Johanna Aalto<sup>2,3</sup>, Stina Krsmanovic<sup>2,3</sup>, Hans Christopher Bernstein<sup>2,3,4</sup>

<sup>1</sup>University of Tromsø;

<sup>2</sup>The Norwegian College of Fishery Science, Faculty of Biosciences, Fisheries and Economics, UiT - The Arctic University of Norway, Tromsø, Norway;

<sup>3</sup>Microalgae & Microbiomes Research Group (M2RG), UiT - The Arctic University of Norway, Tromsø, Norway;

<sup>4</sup>The Arctic Centre for Sustainable Energy (ARC), UiT - The Arctic University of Norway, Tromsø, Norway



**William Arnli**

University of Tromsø

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5ow14g1b/v1>

**Protocol Citation:** William Arnli, Nerea Johanna Aalto, Stina Krsmanovic, Hans Christopher Bernstein 2025. Co-cultivation Protocol for *Porosira glacialis* Diatom and Bacterial Isolates . protocols.io  
<https://dx.doi.org/10.17504/protocols.io.8epv5ow14g1b/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 21, 2025

**Last Modified:** July 16, 2025

**Protocol Integer ID:** 220677

**Keywords:** Phycosphere , Growth dynamics, Diatom cultivation , Synthetic Consortia, Characterization, Diatom-Bacteria Growth Dynamics, IAA, Indole-3-Acetic Acid, chlorophyll a extraction, Cell quantification , Cold-adapted diatom, *Porosira glacialis* , *Sulfitobacter marinus*, Growth analysis, Bacterial associates, Synthetic consortia, Axenic cultivation, Paired cultivation, Chlorophyll a extraction, Marine growth media, Utermöhl settling method, Microalgae, Orbital shaker, Research protocol, Sustainable energy, Cocultivation , Co-cultivation, adapted marine diatom *porosira*, *porosira glacialis* diatom, marine diatom *porosira*, various algal strain, bacteria growth dynamic, bacterial isolate, uit201 diatom, microbiome, bacterial associate, microalgae, growth curves from measured cell density, paired cultivation

## Disclaimer


This protocol was developed as part of a master's thesis affiliated with **[Microalgae & Microbiomes - The M2 Research Group](#)** at **[Norges fiskerihøgskole](#)**, University of Tromsø, Norway

## Abstract

The protocol provides a platform for controlled experimentation on axenic and paired cultivation of the cold-adapted marine diatom *Porosira glacialis* UiT201 diatom with its bacterial associates in a synthetic consortia under variable environmental conditions, allowing the characterization of diatom-bacteria growth dynamics. Treatments include UiT201-*Sulfitobacter marinus* co-culture and UiT201 treated with Indole-3-Acetic Acid (IAA); however, this protocol can be adapted to various algal strains, bacterial isolates, and conditions to study diatom-bacteria growth dynamics. Expected results include growth curves from measured cell density and chlorophyll a extraction.

This protocol was developed as part of a master's thesis affiliated with **[Microalgae & Microbiomes - The M2 Research Group](#)** at **[Norges fiskerihøgskole](#)**, University of Tromsø, Norway

## Guidelines



Suggested experimental workflow

## Materials

### MEADIA

FMAP (1 L)

🧪 15 g Marine Broth BD Difco™ Dehydrated Culture Media: Marine Broth 2216 **Thermo Fisher Scientific**

🧪 5 g Peptone

Peptone from casein (Tryptone) **Merck MilliporeSigma (Sigma-Aldrich)**

🧪 300 mL Filtered seawater

🧪 700 mL MilliQ water



## Protocol materials

☒ Sodium metasilicate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #NA**

☒ 1H-Indole-3-acetic acid, 98% **Thermo Scientific**

☒ 96-100% ethanol

☒ Guillard's (F/2) Marine Water Enrichment Solution 50 x **Merck MilliporeSigma (Sigma-Aldrich)**

☒ Peptone from casein (Tryptone) **Merck MilliporeSigma (Sigma-Aldrich)**

☒ BD Difco™ Dehydrated Culture Media: Marine Broth 2216 **Thermo Fisher Scientific**

## Safety warnings

- ⚠ Refer to the Safety Data Sheets (SDS) for information on health and environmental hazards. Be sure to follow your institution's Health and Safety (HMS) guidelines, wear appropriate personal protective equipment (PPE), and adhere to all established safety protocols to minimize risks

## Before start

This protocol outlines the axenic and paired cultivation of *Porosira glacialis* UiT201 and *Sulfitobacter marinus*, including IAA treatment; however, it can easily be adapted to other algal strains, bacteria isolates, and variable conditions to investigate diatom-bacteria growth dynamics. FMAP media was utilized as an all-purpose marine growth media in the creation of this protocol (see Materials for media composition), however it is highly recommended to establish optimal growth media for each bacterial strain. Results are assessed through cell density and chlorophyll a extraction.

Please note that your experimental design may vary. For this protocol, samples are measured and Indole-3-acetic acid (IAA) is supplemented every other day. Chlorophyll a and cell counts are recorded continuously throughout the experiment, which lasts until two repeated measurements are taken during the stationary phase. A suggested workflow is provided in the Guidelines section.

Experimental overview

| Replicate          | A1                                                     | A2            | A3            | B1                                   | B2                                   | B3                                   | C1                     | C2                     | C3                     |
|--------------------|--------------------------------------------------------|---------------|---------------|--------------------------------------|--------------------------------------|--------------------------------------|------------------------|------------------------|------------------------|
| Treatment          | <i>UIT201</i>                                          | <i>UIT201</i> | <i>UIT201</i> | <i>UIT201</i> +<br><i>S. marinus</i> | <i>UIT201</i> +<br><i>S. marinus</i> | <i>UIT201</i> +<br><i>S. marinus</i> | <i>UIT201</i><br>+ IAA | <i>UIT201</i><br>+ IAA | <i>UIT201</i><br>+ IAA |
| Tot. Volume        | 500 ml                                                 |               |               |                                      |                                      |                                      |                        |                        |                        |
| Media              | F/2 8 ml/l (50x stock)+ Silicate 8 ml/l (3,5g/l stock) |               |               |                                      |                                      |                                      |                        |                        |                        |
| Treatment<br>Conc. | -                                                      |               |               | 50:1 bacteria:diatom ratio           |                                      |                                      | 50 nM                  |                        |                        |
| Treatment<br>added | -                                                      |               |               | Once at T=0                          |                                      |                                      | Every 2 days           |                        |                        |

Sampling and measurement scheme

|             |              |
|-------------|--------------|
| Chl a       | Every 2 days |
| Algae count | Every 2 days |

## Cultivation set-up

10m

- 1
- Position a large orbital shaker in a 

🌡️ 7 °C

 cold-room and arrange flexible 200 × 50W RGB LED strips along the shaker to provide a sufficient cultivation light source

5m

Equipment

PSU-20i Orbital Shaking Platform

NAME

Orbital Shaking Platform

TYPE


Grant bio

BRAND

NA

SKU

<https://www.grantinstruments.com/products/psu-20i-orbital-shaking-platform><sup>LINK</sup>



PSU-20i Orbital Shaking Platform with RGB light fixture

## Day 1 (T= -2) - Treatment preparations

3h 20m

- 2 Inoculate bacteria isolates in 5 mL FMAP media in 15 mL culture tubes and shake Overnight

10m



- 3 Prepare a stock solution of 1M 0.25 millimolar (mM) 1H-Indole-3-acetic acid, 98% Thermo Scientific dissolved in Milli-Q water. Wrap the stock solution in tinfoil to prevent degradation, and place in a refrigerator for short term storage. For long term storage, consider storing at -20 °C .

10m


- 4 The Utermöhl settling method (Edler, L., and Elbrächter, M., 2010) is utilized to count the stock algae culture: Sample  4 mL of stock diatom culture and fix  2 mL diatoms in duplicates with  7  $\mu\text{L}$  Lugol's iodine solution in  2 mL cell culture wells for 1-2 hours in a  4  $^{\circ}\text{C}$  refrigerator.

2h



- 5 Once the samples have settled, observe it under a microscope. Diatom cells are identified and counted in several randomly selected fields of view to ensure a representative sample. If the diatom concentration is too large for an accurate count ( $>100$ ), determine a dilution factor and proceed to dilute the culture. For each duplicate sample, count the diatom cells first horizontally, then vertically along an axis. For the horizontal count, exclude cells on the bottom edge of the microscopic view, and for the vertical count, exclude cells on the right edge of the microscopic view. Repeat the process for the second duplicate sample. This results in four values which are averaged to calculate the concentration of cells per unit volume of the original water sample.

1h



*Counting diatom cells using the Utermöhl method.* Cells are counted horizontally and vertically, excluding those on the bottom and right edges (outlined in red boxes). This process is repeated for the second duplicate sample.

- 6 Aiming for a starting culture of 1-2 mill cells/L, dilute with Milli-Q to 1:3 of the starting algae concentration.

#### Note

This step decreases lag-phase, ensuring that a quick transition to log phase. Skipping this step could result in an extended lag-phase.



## Day 2 (T= -1) - Treatment preparations

15m

- 7 Move the overnight bacteria inoculation tubes to a 7-degree cold room to acclimate. Leave Overnight for additional inoculation.

5m



- 8 Prepare a stock of [M] 3.5 g/L Sodium metasilicate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #NA** in Milli-Q

10m

## Day 3 (T=0) - Initiating cultivation

4h 35m

- 9 For determining stock diatom culture sample, utilize the Utermöhl settling method: 4 mL of stock diatom culture and fix 2 mL diatoms in duplicates with 7 µL Lugol's iodine solution in 2 mL cell culture wells for 1-2 hours.

2h

- 10 Prepare fresh media by adding [M] 8 mg/mL of Guillard's (F/2) Marine Water Enrichment Solution 50 x **Merck MilliporeSigma (Sigma-Aldrich)** and [M] 8 mg/mL of the stock sodium metasilicate to 5 L of filtered seawater\*

5m

\* 5 L is sufficient for 9 (3 treatments \* 3 biological replicates) replicates

- 11 Filter the media through 0.22 µm sterile filters and fill 9 500 mL tissue culture flasks

30m

### Note

Follow and do not exceed the recommended capacity of the filters

| Equipment                                                                                                                                                             |       |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| <b>Sterivex™ Pressure Filter</b>                                                                                                                                      | NAME  |
| 0.22 µm filter                                                                                                                                                        | TYPE  |
| Millipore®                                                                                                                                                            | BRAND |
| 00-00-0                                                                                                                                                               | SKU   |
| <a href="https://www.sigmaaldrich.com/NO/en/product/mm/svgp010#product-documentation">https://www.sigmaaldrich.com/NO/en/product/mm/svgp010#product-documentation</a> | LINK  |

| Equipment                                                                                                                                                                                                                                                                   |       |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| <b>Masterflex L/S</b>                                                                                                                                                                                                                                                       | NAME  |
| Peristaltic Pump                                                                                                                                                                                                                                                            | TYPE  |
| Masterflex                                                                                                                                                                                                                                                                  | BRAND |
| NA                                                                                                                                                                                                                                                                          | SKU   |
| <a href="https://www.avantorsciences.com/us/en/product/39208424/masterflex-ls-digital-pump-system-with-gore-sta-pure-pfl-tubing-avantor">https://www.avantorsciences.com/us/en/product/39208424/masterflex-ls-digital-pump-system-with-gore-sta-pure-pfl-tubing-avantor</a> | LINK  |

- 12

After 1-2 hours of fixation, utilize the Utermöhl method in step 5 to count the stock algae culture.

1h
- 13

Aiming for a starting culture of 1-2 mill cells/L, transfer diatoms to the 

🧪 500 mL

 tissue culture flasks

10m



- 14 Determine bacteria isolate concentrations by measuring OD600. Aiming for a 50:1 bacteria:diatom ratio, proceed to pipette the appropriate amount to each of the bacteria treatment culture flasks 30m
- 15 From the [M] 0.25 millimolar (mM) 5m  
⚗ 1H-Indole-3-acetic acid, 98% Thermo Scientific stock spike 🧪 100  $\mu$ L of IAA to each of the IAA treatment culture flasks to achieve a final concentration of [M] 50 nanomolar (nM) . Immediately proceed to step 16.
- 16 Turn over the culture flasks 10 times to homogenize the cultures and place the flasks on the orbital shaker in the 🌡 7 °C cold-room. Set the orbital shaker to 🌀 50 rpm, 7°C 15m  
with 300 reciprocal motion for 12 seconds and vibration set to 2.

#### Note

Measure light intensity ensuring consistent light intensity across all flasks. It is recommended to alternate the positioning of the bottles.



Culture bottle placement on orbital shaker in the 🌡️ 7 °C cold-room



Orbital shaker settings



## Sampling for chlorophyll extraction and cell count

30m

- 17 Retrieve and extract  25 mL samples from the each of the culture bottles in  50 mL falcon tubes. To minimize incubation room downtime, immediately return the culture flasks to the orbital shaker.

### Note

For each measurement  4 mL is needed for cell counting and  20 mL for chlorophyll extraction. Extracting  25 mL of sample allows for some margin. As the diatom culture enters log phase, the samples need to be diluted for measurements and less extracted sample volume is needed.

## Diatom cell count

2h

- 18 From each of the  25 mL samples, the utilize Utermöhl settling method is used to create duplicates count the cells (steps 4 and 5): pipette 2\*  2 mL of culture to  2 mL cell culture wells and fix with  7  $\mu$ L Lugol's iodine solution for 1-2 hours. Proceed to count the cells according to the previously described method.

2h

## Chlorophyll extraction

15h

- 19 Assemble a vacuum pump with a fitted filtration system. From the remaining  21 mL of each sample culture, create 5 technical replicates by depositing 5 \*  4 mL on  1.2  $\mu$ m glass microfiber filter papers. Proceed to create a vacuum and filter the culture samples.

1h

## Equipment

Grade GF/C Glass Microfiber Filter Papers

NAME

Microfiber Filter Papers

TYPE

Whatman

BRAND

NA

SKU

<https://www.tischscientific.com/whatman/filter-papers-glass-microfiber-filters-1822-025>

LINK

1.2um, 25mm Dia., Binder-Free

SPECIFICATIONS

## Equipment

Chemical Duty Pump

NAME

Vacuum pump

TYPE

Merck Millipore

BRAND

NA

SKU

<https://www.sigmaaldrich.com/NO/en/product/mm/wp6122050>

LINK

220 V/50 Hz

SPECIFICATIONS

## Equipment

1225 Sampling Manifold

NAME

Millipore

BRAND

NA

SKU

<https://www.sigmaaldrich.com/NO/en/product/mm/xx2702550><sup>LINK</sup>

15 mL process volume

SPECIFICATIONS

- 20 Transfer the glass microfiber filters to separate test tubes and add  5 mL   96-100% ethanol . For extraction, immediately cover the samples in tinfoil and place them in a  7 °C cold-room  Overnight . 
- 21 After allowing overnight extraction, proceed to measure chlorophyll a on a fluorometer with excitation wavelength set to 430-450 nm (typically around 435 nm) and emission wavelength to 675-685 nm. A chlorophyll a factor of 0.00417 is used to convert absorbance readings into the concentration of chlorophyll a (Chla µg/L) in a sample.  

## Equipment

TD-700 Laboratory Fluorometer

NAME

Fluorometer

TYPE

Turner Designs

BRAND

00-00-0

SKU

[Discontinued](#)

LINK



## Protocol references

Edler, L., and Elbrächter, M. (2010). "The utermöhl method for quantitative phytoplankton analysis," in *Microscopic and Molecular Methods for Quantitative Phytoplankton Analysis*, Vol. 110, eds B. Karlson, C. Cusack, and E. Bresnan (Paris: Intergovernmental Oceanographic Commission Manual and guides), 13–20