

Mar 02, 2024

Co-extraction of RNA and DNA from animal tissue

DOI

dx.doi.org/10.17504/protocols.io.ewov1qrkpgr2/v1



Dominik Buchner¹

¹University of Duisburg-Essen, Aquatic Ecosystem Research



Dominik Buchner

University of Duisburg-Essen, Aquatic Ecosystem Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1qrkpgr2/v1>

Protocol Citation: Dominik Buchner 2024. Co-extraction of RNA and DNA from animal tissue. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1qrkpgr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 02, 2024

Last Modified: March 02, 2024

Protocol Integer ID: 96045

Keywords: dna from animal tissue sample, dna from animal tissue, animal tissue sample, extraction, rna in the flow, rna, interested in the rna, dna, dna spin, animal tissue, cell debris, remaining protein, tissue, separate tube

Abstract

This protocol describes how to co-extract RNA and DNA from animal tissue samples. Samples are homogenized and simultaneously lyzed by bead-beating. Cell debris is pelleted by centrifugation, the DNA is then subsequently bound to a silica column, while the RNA passes the membrane. The RNA in the flow-through is then precipitated with 70% ethanol and bound to a second silica column. Both, DNA and RNA are washed with different wash buffers to remove remaining proteins and other contaminants and finally eluted in separate tubes. If the user is just interested in the RNA, the DNA spin-column can just be discarded.

Guidelines


Follow general lab etiquette. Wear gloves to prevent contamination of samples. Clean the workspace before starting and after finishing with 80% EtOH.


Materials

Materials required:

Below all materials needed for the protocol are listed. Vendors and part numbers are listed but interchangeable depending on the supply situation.

Chemicals:


Guanidinium thiocyanate  Guanidinium thiocyanate **Fisher Scientific Catalog #10503345**

Tris ultrapure 99.9%  Tris ultrapure 99.9% **Diagonal Catalog #A1086.1000**

Hydrochloric acid fuming 37%

 Hydrochloric acid fuming 37% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #1003171011**

Pre-filter columns  Pre Filter Columns - 850 µl **Biopolymer Isolation Technologies Catalog #MC-01P-100**


Guanidinium chloride  Guanidine hydrochloride **Fisher Scientific Catalog #10543325**


Ethanol absolute  Ethanol absolute 99.8% **Fisher Scientific Catalog #11994041**


Antifoam solution (optional):  Silicon-Antischaumemulsion 30 **Carl Roth Catalog #0734.1**

Labware:

2 mL screwcap tubes  2 mL screwcap tube **Sarstedt Catalog #72.693**




2 mm zirconia beads  Zirconia Beads 2 mm dia **BioSpec Products Catalog #11079124zx**


0.1 mm glass beads  Glass Beads 0.1 mm dia **BioSpec Products Catalog #11079101**




EconoSpin mini spin column  EconoSpin mini spin column with lid **Epoch Life Science Catalog #1920-050**

Stock solutions:

 1 L Tris stock solution **[M] 1 Molarity (M)** **pH 7.5**





- Add  121.1 g Tris ultrapure 99.9% to a beaker
- Adjust volume to  800 mL with ddH₂O
- Adjust pH to **pH 7.5** with HCl
- Adjust volume to  1 L with ddH₂O



 1 L sodium chloride stock solution **[M] 5 Molarity (M)**






- Add  292.2 g sodium chloride to a beaker
- Adjust volume to  1 L with ddH₂O
- Sterilize by filtering and store at  Room temperature

 1 L Tris stock solution **[M] 1 Molarity (M)** **pH 8.5**













- Add  121.1 g Tris ultrapure 99.9% to a beaker
- Adjust volume to  800 mL with ddH₂O
- Adjust pH to  8.5 with HCl
- Adjust volume to  1 L with ddH₂O



 1 L DNA wash buffer 2 stock solution  50 millimolar (mM) Tris  7.5








- Add  50 mL of  1 Molarity (M) Tris stock solution  7.5 to a beaker
- Adjust volume to  1 L with ddH₂O
- Sterilize by filtering and store at  Room temperature




Working solutions:








 1 L GITC lysis buffer ( 4 Molarity (M) Guanidinium thiocyanate ,  10 millimolar (mM) Tris)  7.5

- Add  472.6 g guanidinium thiocyanate to a beaker
- Add  10 mL of  1 Molarity (M) Tris stock solution  7.5
- Adjust volume to  1 L with ddH₂O
- Stir until the GITC is completely dissolved (heating will speed this up)
- Sterilize by filtering and store at  Room temperature




 1 L RNA wash buffer 1 ( 900 millimolar (mM) Guanidinium thiocyanate ,  10 millimolar (mM) Tris ,  20 % (v/v) Ethanol absolute)  7.5








- Add  106.3 g guanidinium thiocyanate to a beaker
- Add  10 mL of  1 Molarity (M) Tris stock solution  7.5
- Add  200 mL Ethanol absolute
- Adjust volume to  1 L with ddH₂O
- Sterilize by filtering and store at  Room temperature




 1 L RNA wash buffer 2 ( 100 millimolar (mM) sodium chloride ,  10 millimolar (mM) Tris ,  80 % (v/v) ethanol absolute)  7.5




- Add  20 mL of  5 Molarity (M) sodium chloride stock solution
- Add  10 mL of  1 Molarity (M) Tris stock solution  7.5
- Adjust volume to  200 mL with ddH₂O
- Adjust volume to  1 L with ethanol absolute




- Sterilize by filtering and store at  Room temperature






 1 L DNA wash buffer 1 ( 2.5 Molarity (M) Guanidinium chloride ,  10 millimolar (mM) Tris ,  57 % (v/v) Ethanol absolute)  7.5

- Add  238.9 g guanidinium chloride to a beaker
- Add  10 mL of  1 Molarity (M) Tris stock solution  7.5
- Adjust volume to  430 mL with ddH₂O to dissolve the GuHCl
- Adjust volume to  1 L with Ethanol absolute
- Sterilize by filtering and store at  Room temperature

 1 L DNA wash buffer 2 ( 10 millimolar (mM) Tris ,  80 % (v/v) ethanol absolute)  7.5

- Add  200 mL DNA wash buffer 2 stock solution to a beaker
- Adjust volume to  1 L with Ethanol absolute
- Sterilize by filtering and store at  Room temperature

 1 L elution buffer  10 millimolar (mM) Tris  8.5

- Add  10 mL of  1 Molarity (M) Tris stock solution  8.5 to a beaker
- Adjust the volume to  1 L with ddH₂O
- Sterilize by filtering and store at  Room temperature

Troubleshooting

Safety warnings

- ⚠ Buffers containing guanidine produce highly reactive compounds when mixed with bleach. Don't mix the extraction waste with bleach or solutions that contain bleach.
Reagents are potentially damaging to the environment. Dispose waste as mandated.


Before start

Make sure all buffers are prepared before starting.




Sample preparation and lysis

5m


- 1 For each sample prepare one 2 mL screwcap tube pre-filled with approximately  400 mg of 2 mm zirconia beads and 0.1 mm glass beads.


Note

Generally, we just add a small spoon of each type of beads to the tube. As long as the tissue is fully homogenized after bead-beating, the amount of beads is sufficient.


- 2 Add up to  30 mg of animal tissue to the prepared tube.


Note

For samples with a high RNA content, less starting material might lead to better results. For most sample types  10 mg of starting material will yield a sufficient amount of DNA and RNA for downstream analysis.

- 3 Add  1000 μ L GITC lysis buffer to the sample tube.

Note

For complete inactivation and destruction RNAses  10 μ L of 2-Mercaptoethanol can be added in addition. We usually don't because then the samples have to be handled under a fume hood until all lysate has been handled and discarded appropriately.

If you experience a lot of foam formation after bead-beating consider adding  30 Parts per Million (PPM) silicone antifoam to the lysis buffer when preparing it. See materials for a recommendation.

- 4 Immediately bead beat for  00:05:00 at maximum speed.

5m



Note

Depending on the bead beater used in this step the time might have to be adjusted. We recommend to bead beat the sample until the material is completely homogenized.


Lysate clearing

10s

- 5  Room temperature, 00:10:00 , at maximum speed

10m

DNA binding

- 6 Transfer  700 μL of the cleared lysate from step 5 to a silica spin column to bind the DNA in the lysate. **Keep the flow-through. Mark the spin column as the DNA column.**



Note

The protocol will work with all kinds of silica spin columns. See the materials section for what we use.

If you are only interested in RNA: If only RNA is of interest the DNA spin column can be discarded at this point in the protocol.

RNA precipitation and binding

15s

- 7 Add  700 μL 70% Ethanol to the flow-through from step 6 to adjust the binding conditions for RNA to bind to the silica column.
- 8 Vortex the samples to mix the lysate with the ethanol. Do not centrifuge.
- 9 Load the mixture on a second spin column. **Mark this column as the RNA spin column.**
 11000 x g, Room temperature, 00:00:15 and discard the flow-through.



15s

**Note**

Two loading steps will be necessary to pass the complete volume through the spin column.




Washing steps

15s




- 10 Add  700 μ L RNA wash buffer 1 to the **RNA spin column**,
 11000 x g, Room temperature, 00:00:15 and discard the flow-through.

Note

For less experienced users: If you are concerned about needing too much time to process both fractions at the same time and risk RNA degradation it is fine to first finish the RNA extraction until safe storage and then finish the DNA fraction.

- 11 Add  500 μ L RNA wash buffer 2 to the **RNA spin column**, add
 500 μ L DNA wash buffer 1 to the **DNA spin column**,
 11000 x g, Room temperature, 00:00:15 and discard the flow-through.





15s

- 12 Add  500 μ L RNA wash buffer 2 to the **RNA spin column**, add
 500 μ L DNA wash buffer 2 to the **DNA spin column**,
 11000 x g, Room temperature, 00:00:15 and discard the flow-through.

15s

Column drying and elution



4m

- 13  11.000 x g, Room temperature, 00:01:00 to dry the silica membrane of the spin columns. Transfer the spin column to a fresh 1.5 mL microcentrifuge tube.
- 14 Add  100 μ L elution buffer directly to the silica membrane. Incubate the column for
 00:03:00 at  Room temperature

1m

3m

**Note**

The volume of the elution buffer can be adjusted in this step if a higher concentration or higher volume is required for downstream analysis. Usually, every volume in the range from  30 μL to  200 μL is fine.

15



11.000 x g, Room temperature, 00:01:00 , store

the eluted RNA at



-80 °C

1m


and the eluted DNA at



-20 °C



Expected result



The described protocol was tested with different kinds of invertebrate samples, we expect it to work with all animal tissue.

Top row: Plecoptera sample and two Trichoptera samples.

Lower row: Three Gammarus samples.

28S/18S bands are clearly visible and should have a clear band. Genomic DNA is free from RNA contamination. There is some DNA contamination in the RNA extracts. If DNA-free RNA is needed for downstream analysis consider treating the RNA samples with DNase and cleaning them up with an RNA cleanup protocol afterward (see [RNA Cleanup with magnetic beads](#)).

