

Feb 10, 2023

## Co-extraction of RNA and DNA from plant tissue

DOI

[dx.doi.org/10.17504/protocols.io.n2bvj8qxngk5/v1](https://dx.doi.org/10.17504/protocols.io.n2bvj8qxngk5/v1)



Dominik Buchner<sup>1</sup>

<sup>1</sup>University of Duisburg-Essen, Aquatic Ecosystem Research



Dominik Buchner

University of Duisburg-Essen, Aquatic Ecosystem Research

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj8qxngk5/v1>

**Protocol Citation:** Dominik Buchner 2023. Co-extraction of RNA and DNA from plant tissue. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.n2bvj8qxngk5/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 08, 2023

**Last Modified:** February 10, 2023

**Protocol Integer ID:** 76687

**Keywords:** dna from plant tissue sample, dna from plant tissue, plant tissue sample, rna, plant tissue, rna in the flow, extraction, interested in the rna, dna, dna spin, cell debris, cell, tissue

## Abstract

This protocol describes how to co-extract RNA and DNA from plant tissue samples. Samples are homogenized and simultaneously lysed by bead-beating. Cell debris is then caught with a pre-filter column, the DNA is then subsequently bound to a silica column, while the RNA passes the membrane. The RNA in the flow-through is then precipitated with 100% ethanol and bound to a second silica column. Both, DNA and RNA are washed with different wash buffers to remove remaining proteins and other contaminants and finally eluted in separate tubes. If the user is just interested in the RNA, the DNA spin-column can just be discarded.

## Guidelines


Follow general lab etiquette. Wear gloves to prevent contamination of samples. Clean the workspace before starting and after finishing with 80% EtOH.

## Materials

### Materials required:

Below all materials needed for the protocol are listed. Vendors and part numbers are listed but interchangeable depending on the supply situation.

### Chemicals:


Guanidinium thiocyanate  Guanidinium thiocyanate **Fisher Scientific Catalog #10503345**


Tris ultrapure 99.9%  Tris ultrapure 99.9% **Diagonal Catalog #A1086.1000**

Hydrochloric acid fuming 37%


 Hydrochloric acid fuming 37% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #1003171011**


Pre-filter columns  Pre Filter Columns - 850 µl **Biopolymer Isolation Technologies Catalog #MC-01P-100**


Guanidinium chloride  Guanidine hydrochloride **Fisher Scientific Catalog #10543325**


Ethanol absolute  Ethanol absolute 99.8% **Fisher Scientific Catalog #11994041**

### Labware:

2 mL screwcap tubes  2 mL screwcap tube **Sarstedt Catalog #72.693**





2 mm zirconia beads  Zirconia Beads 2 mm dia **BioSpec Products Catalog #11079124zx**

0.1 mm glass beads  Glass Beads 0.1 mm dia **BioSpec Products Catalog #11079101**




EconoSpin mini spin column  EconoSpin mini spin column with lid **Epoch Life Science Catalog #1920-050**

### Stock solutions:

 1 L Tris stock solution [M] 1 Molarity (M)  7.5

- Add  121.1 g Tris ultrapure 99.9% to a beaker
- Adjust volume to  800 mL with ddH<sub>2</sub>O
- Adjust pH to  7.5 with HCl
- Adjust volume to  1 L with ddH<sub>2</sub>O

 1 L sodium chloride stock solution [M] 5 Molarity (M)

- Add  292.2 g sodium chloride to a beaker
- Adjust volume to  1 L with ddH<sub>2</sub>O
- Sterilize by filtering and store at  Room temperature

 1 L Tris stock solution [M] 1 Molarity (M)  8.5

- Add  121.1 g Tris ultrapure 99.9% to a beaker



- Adjust volume to 800 mL with ddH<sub>2</sub>O
- Adjust pH to 8.5 with HCl
- Adjust volume to 1 L with ddH<sub>2</sub>O

1 L DNA wash buffer 2 stock solution 50 millimolar (mM) Tris 7.5

- Add 50 mL of 1 Molarity (M) Tris stock solution 7.5 to a beaker
- Adjust volume to 1 L with ddH<sub>2</sub>O
- Sterilize by filtering and store at Room temperature

### Working solutions:

1 L GITC lysis buffer ( 4 Molarity (M) Guanidinium thiocyanate , 10 millimolar (mM) Tris ) 7.5

- Add 472.6 g guanidinium thiocyanate to a beaker
- Add 10 mL of 1 Molarity (M) Tris stock solution 7.5
- Adjust volume to 1 L with ddH<sub>2</sub>O
- Stir until the GITC is completely dissolved (heating will speed this up)
- Sterilize by filtering and store at Room temperature

1 L RNA wash buffer 1 ( 900 millimolar (mM) Guanidinium thiocyanate , 10 millimolar (mM) Tris , 20 % (v/v) Ethanol absolute ) 7.5

- Add 106.3 g guanidinium thiocyanate to a beaker
- Add 10 mL of 1 Molarity (M) Tris stock solution 7.5
- Add 200 mL Ethanol absolute
- Adjust volume to 1 L with ddH<sub>2</sub>O
- Sterilize by filtering and store at Room temperature

1 L RNA wash buffer 2 ( 100 millimolar (mM) sodium chloride , 10 millimolar (mM) Tris , 80 % (v/v) ethanol absolute ) 7.5

- Add 20 mL of 5 Molarity (M) sodium chloride stock solution
- Add 10 mL of 1 Molarity (M) Tris stock solution 7.5
- Adjust volume to 200 mL with ddH<sub>2</sub>O
- Adjust volume to 1 L with ethanol absolute
- Sterilize by filtering and store at Room temperature



🧪 1 L DNA wash buffer 1 ( [M] 2.5 Molarity (M) Guanidinium chloride , [M] 10 millimolar (mM) Tris , [M] 57 % (v/v) Ethanol absolute ) Ⓟ 7.5

- Add 🧪 238.9 g guanidinium chloride to a beaker
- Add 🧪 10 mL of [M] 1 Molarity (M) Tris stock solution Ⓟ 7.5
- Adjust volume to 🧪 430 mL with ddH<sub>2</sub>O to dissolve the GuHCl
- Adjust volume to 🧪 1 L with Ethanol absolute
- Sterilize by filtering and store at 🌡 Room temperature

🧪 1 L DNA wash buffer 2 ( [M] 10 millimolar (mM) Tris , [M] 80 % (v/v) ethanol absolute ) Ⓟ 7.5

- Add 🧪 200 mL DNA wash buffer 2 stock solution to a beaker
- Adjust volume to 🧪 1 L with Ethanol absolute
- Sterilize by filtering and store at 🌡 Room temperature

🧪 1 L elution buffer [M] 10 millimolar (mM) Tris Ⓟ 8.5

- Add 🧪 10 mL of [M] 1 Molarity (M) Tris stock solution Ⓟ 8.5 to a beaker
- Adjust the volume to 🧪 1 L with ddH<sub>2</sub>O
- Sterilize by filtering and store at 🌡 Room temperature

## Troubleshooting

## Safety warnings

- ❗ Buffers containing guanidine produce highly reactive compounds when mixed with bleach. Don't mix the extraction waste with bleach or solutions that contain bleach.  
Reagents are potentially damaging to the environment. Dispose waste as mandated.



## Before start

Make sure all buffers are prepared before starting.





## Sample preparation and lysis

5m

- 1 For each sample prepare one 2 mL screwcap tube pre-filled with approximately  400 mg of 2 mm zirconia beads and 0.1 mm glass beads.
- 2 Add up to  200 mg of plant tissue to the prepared tube.

### Note

For samples with a high RNA content less starting material might lead to better results. For most sample types  50 mg of starting material will yield a sufficient amount of DNA and RNA for downstream analysis.

- 3 Add  800 µL GITC lysis buffer to the sample tube.

### Note

For complete inactivation and destruction RNAses of 2-Mercaptoethanol can be added in addition. We usually don't because then the samples have to be handled under a fume hood until all lysate has been handled and discarded appropriately.

- 4 Immediately bead beat for  00:05:00 at maximum speed.

10s


### Note

Depending on the bead beater used in this step the time might have to be adjusted. We'd recommend to bead beat the sample until the material is completely homogenized.


## Lysate clearing and pre-filtering


10s



5  Room temperature, 00:00:10 , at maximum speed


10s

6 Transfer  700  $\mu$ L of the crude lysate to a pre-filter column.

7  Room temperature, 00:10:00 , at maximum speed

10m

## DNA binding


8 Transfer  700  $\mu$ L of the flowthrough from step 7 to a silica spin column to bind the DNA in the lysate. **Keep the flow-through. Mark the spin column as the DNA column.**

### Note


The protocol will work with all kinds of silica spin columns. See materials section for what we use.

## RNA precipitation and binding

15s

9 Add  350  $\mu$ L Ethanol absolute to the flow-through from step 8 to adjust the binding conditions to bind RNA to the silica column.

10 Vortex the samples to mix the lysate with the ethanol. Do not centrifuge.

11 Load the mixture on a second spin column. **Mark this column as the RNA spin column.**  
 11000 x g, Room temperature, 00:00:15 and discard the flow-through.

15s









### Note

Two loading steps will be necessary to pass the complete volume through the spin column.










## Washing steps

15s

- 12 Add  700  $\mu$ L RNA wash buffer 1 to the **RNA spin column**,  
 11000 x g, Room temperature, 00:00:15 and discard the flow-through. 15s
- 13 Add  500  $\mu$ L RNA wash buffer 2 to the **RNA spin column**, add  
 500  $\mu$ L DNA wash buffer 1 to the **DNA spin column**,  
 11000 rpm, Room temperature, 00:00:15 and discard the flow-through. 15s
- 14 Add  500  $\mu$ L RNA wash buffer 2 to the **RNA spin column**, add  
 500  $\mu$ L DNA wash buffer 2 to the **DNA spin column**,  
 11000 rpm, Room temperature, 00:00:15 and discard the flow-through. 15s


## Column drying and elution

4m

- 15  11.000 rpm, Room temperature, 00:01:00 to dry the silica membrane of the spin columns. Transfer the spin column to a fresh 1.5 mL microcentrifuge tube. 1m
- 16 Add  100  $\mu$ L elution buffer directly to the silica membrane. Incubate the column for  
 00:03:00 at  Room temperature 3m
- 17  11.000 rpm, Room temperature, 00:01:00 , store the eluted RNA at  -80 °C and  
the eluted DNA at  -20 °C 1m



## Expected result



Expected result of the described protocol. Extraction was carried out in 4 replicates, left part of the gel picture shows the DNA fraction of the sample, while the right part shows the RNA fraction.