

Jan 11, 2026

Co-Immunoprecipitation (co-IP) Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk8z6xl5r/v1>

Ali Ghoochani^{1,2,3,4}, Monther Abu-Remaileh^{1,2,3,4}

¹Department of Chemical Engineering, Stanford University, Stanford, CA 94305, USA;

²Department of Genetics, Stanford University, Stanford, CA 94305, USA;

³The Institute for Chemistry, Engineering and Medicine for Human Health (Sarafan ChEM-H), Stanford University, Stanford, CA 94305, USA;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA

asap



Ali Ghoochani

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk8z6xl5r/v1>

Protocol Citation: Ali Ghoochani, Monther Abu-Remaileh 2026. Co-Immunoprecipitation (co-IP) Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk8z6xl5r/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 13, 2024

Last Modified: January 11, 2026

Protocol Integer ID: 109773

Keywords: protein interaction, studying protein, protocol, cellular model, ip, protein

Abstract

Co-immunoprecipitation (co-IP) is a widely utilized technique for studying protein-protein interactions, particularly in cellular models.



Materials

- 1x PBS (Phosphate-Buffered Saline)
- 1:10 diluted extraction buffer:
- 50 mM HEPES, pH 7.4
- 40 mM NaCl
- 2 mM EDTA
- 1% Triton X-100
- 1.5 mM sodium orthovanadate (NaVO₄)
- 30 mM sodium fluoride (NaF)
- 10 mM sodium pyrophosphate (Na₄P₂O₇)
- 10 mM sodium β-glycerophosphate
- Protease inhibitor cocktail
- Anti-Flag magnetic beads (Pierce™ Anti-DYKDDDDK Magnetic Agarose (Anti-FLAG)(Thermo Fisher Scientific)
- SDS-PAGE materials (Thermo Fisher Scientific)
- Loading dye (5x and 1x)
- Centrifuge
- Shaker
- Incubator at 4°C
- Immunoblotting equipment

Extraction Buffer (1x):

- **50 mM HEPES, pH 7.4**
- **40 mM NaCl**
- **2 mM EDTA**
- **1% Triton X-100**
- **1.5 mM Sodium Orthovanadate (NaVO₄)**
- **30 mM Sodium Fluoride (NaF)**
- **10 mM Sodium Pyrophosphate (Na₄P₂O₇)**
- **10 mM Sodium β-Glycerophosphate**



Cell Culture

3d

- 1 Seed 15×10^6 SH-SY5Y cells expressing the Flag-fusion protein (here: Flag-tagged SLC45A1) and control cells in 15 cm dishes
- 1.1 Incubate the cells for 72:00:00 at 37 °C in a humidified atmosphere with 5% CO₂, or until the cells reach confluence

3d

Cell Lysis

1h

- 2 Wash the cells gently with 10 mL 1x PBS to remove any residual media.
- 2.1 Prepare 1 mL of **1:10 diluted extraction buffer** supplemented with a protease inhibitor cocktail.
Extraction buffer: Please see "Materials"
- 2.2 Add the extraction buffer and lyse the cells. Use a cell lifter to collect the cells, then transfer the lysate to a new tube.



It is highly recommended to use low-protein-binding tubes for all steps.

- 2.3 Incubate lysate at 4 °C with vortexing for 01:00:00
- 2.4 Centrifuge the lysate at maximum speed for 00:10:00 at 4 °C
- 2.5 Carefully collect the supernatant to new tube
- 2.6 Take 100 µL of the supernatant and add loading dye. Mix well and set aside for later use as the input sample.

1h

10m

Immunoprecipitation

1h



- 3 Transfer the remaining  900 μL of supernatant to a new tube containing  25 μL of **pre-washed Anti-Flag magnetic beads** (Beaded were washed with **1:10 diluted extraction buffer**).
- 3.1 Incubate the mixture at  4 $^{\circ}\text{C}$ with gentle rotation overnight. Alternatively, a shorter incubation of 3 hours can also be used if desired
- 3.2 After incubation, wash the beads **seven times** with **1 ml** of **1:10 diluted extraction buffer**. Use a magnet to facilitate the separation of the magnetic beads during the washing process

Important: Change tubes every two washes to ensure effective washing.
- 3.3 After the final wash, add  50 μL of 1x loading dye to the beads. Mix thoroughly.
- 3.4 Alternatively, if non-denatured protein samples are needed, incubate the beads with 200 $\mu\text{g/mL}$ Flag-peptide (prepared in 1:10 diluted extraction buffer  01:00:00 to  02:00:00 at  4 $^{\circ}\text{C}$

3h

It is recommended to optimize the Flag-peptide concentration and incubation time for your specific Flag-fusion protein