



Mar 21, 2025

Co-immunoprecipitation ER-phagy receptors

 [Nature Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbrkd1lpk/v1>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbrkd1lpk/v1>

External link: <https://doi.org/10.1038/s41556-025-01712-y>

Protocol Citation: Elias Adriaenssens 2025. Co-immunoprecipitation ER-phagy receptors. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgbrkd1lpk/v1>

Manuscript citation:

Adriaenssens, E., Schaar, S., Cook, A.S.I. *et al.* Reconstitution of BNIP3/NIX-mitophagy initiation reveals hierarchical flexibility of the autophagy machinery. *Nat Cell Biol* **27**, 1272–1287 (2025). <https://doi.org/10.1038/s41556-025-01712-y>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 22, 2025

Last Modified: March 21, 2025

Protocol Integer ID: 119233

Keywords: ASAPCRN, phagy receptor, phagy receptors this protocol

Funders Acknowledgements:

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Abstract

This protocol details the co-immunoprecipitation of ER-phagy receptors.

Materials

 pHAGE-TEX264-GFP **addgene Catalog #201925** or

 pHAGE-FAM134C-GFP **addgene Catalog #201927**

 Earle's Balanced Salts **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Earle's Balanced Salts**

 Pierce™ Detergent Compatible Bradford Assay Kit **Thermo Fisher Catalog #23246**

 ChromoTek GFP-Trap® Agarose **Proteintech Catalog # gta-20**

 NuPAGE™ 4-12% Bis-Tris Protein Gels, 1.0 mm, 12-well **Thermo Fisher Catalog #NP0322BOX**

Lysis buffer

	A	B
	KCl	150 mM
	MgCl ₂	2.5 mM
	Tris-HCl pH 7.4	20 mM
	NP-40	0.5%

Washing buffer

	A	B
	KCl	150 mM
	MgCl ₂	2.5 mM
	Tris-HCl pH 7.4	20 mM



Procedure

6h 35m

- 1 Wild-type HeLa cells are transiently transfected with pHAGE_TEX264-GFP (RRID:Addgene_201925) or pHAGE_FAM134C-GFP (RRID:Addgene_201927) using Lipofectamine 3000 (Thermo Fisher).
- 2 After seeding, treat the cells for  06:00:00 with Earle's balanced salt medium starvation medium to induce ER-phagy. 
- 3 After the treatment, collect the cells by trypsinization and wash the cell pellet with PBS once before cells are lysed in lysis buffer. 

Lysis buffer:

	A	B
	KCl	150 mM
	MgCl ₂	2.5 mM
	Tris-HCl pH 7.4	20 mM
	NP-40	0.5%

- 4 Lyse the samples for  00:20:00  On ice before cell lysates are cleared by centrifugation at  20000 x g, 4°C, 00:10:00 . 

- 5 Then, determine the protein concentrations of the cleared protein lysates with the Pierce Detergent Compatible Bradford Assay Kit.
- 6 Incubate the  6 mg of cell lysate  Overnight with  20 µL of GFP-Trap agarose beads. 
- 7 In the morning, wash the samples three times in washing buffer before the beads are resuspended in protein loading dye, supplement with  100 millimolar (mM) DTT, and boil for  00:05:00 at  95 °C . 


Washing buffer:



	A	B
	KCl	150 mM
	MgCl ₂	2.5 mM
	Tris-HCl pH 7.4	20 mM

- 8 Load the samples on 4-12% SDS-PAGE gels and analyze by western blotting as described above.