

Jul 25, 2023

Co-immunoprecipitation using GFP-trap

DOI

dx.doi.org/10.17504/protocols.io.6qpvr36o2vmk/v1

Dan Dou¹, Erika Holzbaur¹

¹University of Pennsylvania



Dan Dou

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr36o2vmk/v1>

Protocol Citation: Dan Dou, Erika Holzbaur 2023. Co-immunoprecipitation using GFP-trap. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.6qpvr36o2vmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 20, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 85310

Keywords: ASAPCRN, hek293 cell, trap bead, using gfp, cell

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000350

Abstract

Here, we describe performing co-immunoprecipitation experiments in HEK293 cells using GFP-trap beads

Troubleshooting

1 See "Protocol: HEK293 cell culture for co-immunoprecipitation experiments" for preceding culture and transfection. For the present study, HEK cells were transfected with EGFP-tagged bait proteins. Prey proteins were either endogenous or transfected with HA- or SNAP-tags.

2 Prepare wash buffer:

	A	B
	Tris-HCl pH 7.5	10 mM
	NaCl	150 mM
	EDTA	0.5 mM
	Triton-X (optional)	0.4%

Note

Triton-X should be included for prey proteins where appreciable non-specific binding is observed in the EGFP vector condition.

3 Prepare lysis buffer. Lysis buffer composition differed for experiments needing lambda phosphatase treatment

3.1 Lysis buffer for non-lambda phosphatase experiment:

	A	B
	Tris-HCl pH 7.5	10 mM
	NaCl	150 mM
	EDTA	0.5 mM
	NP-40	0.5%
	PMSF	1 mM
	TAME	0.01 mg/mL

	A	B
	Leupeptin	0.01 mg/mL
	Pepstatin A	0.001 mg/mL

3.2 Lysis buffer compatible with lambda phosphatase experiment:

	A	B
	1x NEBuffer for Protein MettaloPhosphatases (New England BioLabs)	1x
	NP-40	0.5%
	Leupeptin	0.01 mg/mL
	ddH ₂ O	To desired volume

- 4 24 hours after transfection, wash HEK cells twice in ice-cold PBS and lyse in appropriate lysis buffer. Use 600 μ L lysis buffer per experimental condition (pooled across the three 10cm dishes).
- 5 Clarify lysates at 10 x g at 4 degrees C for 10 minutes
- 6 Wash 25 μ L GFP-trap beads per experimental condition in 500 μ L wash buffer, in low protein binding Eppendorf tubes under rotating agitation.

Note

For large protein complexes, GFP-Trap Magnetic Particles M-270 should be used instead of GFP-Trap Magnetic Agarose beads

- 7 Equilibrate beads in lysis buffer for 5 min at 4 degrees C under rotating agitation



- 8 For lambda phosphatase experiments: add 60 μL of 10 mM MnCl_2 and 24 μL lambda phosphatase (2,400 units; New England BioLabs) for a final reaction volume of 600 μL per experimental condition.

Note

For conditions that are lambda phosphatase-negative, add 24 μL ddH₂O instead

- 8.1 Incubate at 30 degrees C for 30 minutes
- 9 Incubate beads with clarified lysate (or lambda phosphatase-treated lysate) for 1 hour at 4 degrees C under rotating agitation
- 10 Wash beads three times for 5 min in wash buffer at 4 degrees C under rotating agitation
- 11 Resuspend beads in 60 μL denaturing buffer, and boil for 10 minutes to release bound protein