

Aug 29, 2024

Version 1

# Collection and shipment of specimen for Visium Spatial Transcriptomics (vST/Visium ST) V.1

DOI

[dx.doi.org/10.17504/protocols.io.36wgqn193gk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqn193gk5/v1)

Laura Robinson<sup>1</sup>, Susan Sheehan<sup>1</sup>, Gaven Garland<sup>1</sup>, Ron Korstanje<sup>1</sup>

<sup>1</sup>The Jackson Laboratory, Bar Harbor, ME, USA

Cellular Senescence Net...



Harshpreet Chandok

The Jackson Laboratory

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.36wgqn193gk5/v1>

**Protocol Citation:** Laura Robinson, Susan Sheehan, Gaven Garland, Ron Korstanje 2024. Collection and shipment of specimen for Visium Spatial Transcriptomics (vST/Visium ST). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgqn193gk5/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 29, 2024

**Last Modified:** August 29, 2024

**Protocol Integer ID:** 106590

**Keywords:** shipment procedures for visium spatial transcriptomic, visium spatial transcriptomic, spatial transcriptomic, visium st, visium, details on specimen collection, shipment of specimen, specimen collection, vst, sennet consortium, sen project in the sennet consortium

**Funders Acknowledgements:**

**National Institute on Aging (NIA) JAX-Sen Senescence Tissue Mapping Center**

**Grant ID:** U54 AG079753

## Abstract

These samples are part of the JAX-Sen project in the SenNet Consortium. Here we provide details on specimen collection and shipment procedures for Visium Spatial Transcriptomics (vST/Visium ST).

## Troubleshooting

## Reagents and Materials:

- 1
  - 10% NBF fixative
  - Tweezers
  - Appropriate container for fixing

## Quality Key Points:

- 2
  - The tissue specimen should be always kept at 4 degrees Celsius and RNase-free.
  - It is crucial to not store the tissue specimen at RT to avoid any cell death, and tissue and/or RNA degradation.

## Procedure:

### 3 **Collection/Harvest:**

- 3.1 Animal was euthanized via cervical dislocation.
- 3.2 Animal was pinned to a necropsy tray which remained on wet ice throughout the harvest.
- 3.3 Animal was perfused with 20ml cold PBS.
- 3.4 Heart, pancreas, and kidney was collected and the specimen was fixed in 10-20 X volume of 10% NBF fixative for 24-48 hours at room temperature.
- 3.5 Samples were kept at 4 degrees Celsius until submission to histology for processing and embedding.

### 4 **Embedding:**

- 4.1 Samples to be embedded into FFPE blocks at Histology, Jax-BH.
- 4.2 Samples can be stored with Histology at Jax-BH (Rama/Juliana will coordinate further procedures)

