

Oct 06, 2022

Colony PCR

Forked from a private protocol

DOI

dx.doi.org/10.17504/protocols.io.bp2l69bnklqe/v1

An.Huang¹

¹XJTLU



An.Huang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l69bnklqe/v1>

Protocol Citation: An.Huang 2022. Colony PCR. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l69bnklqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 05, 2022



Last Modified: October 06, 2022

Protocol Integer ID: 70896

Keywords: simple colony pcr procedure, colony pcr colony pcr, pcr procedure, purified dna template, pcr, dna, procedure, whole organism of bacteria, bacteria

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Colony PCR using the whole organism of bacteria instead of purified DNA template. This simplifies PCR procedure. This protocol helps conduct a simple colony PCR procedure.

Materials

Your LB agar plate with colonies of transformed E. coli, TE solution, PCR buffer, dNTP stock solution, Forward Primer





Reverse Primer, Taq DNA polymerase

Troubleshooting



Preparation for experiment

2d

- 1 Pick several colonies of bacteria from the plate using pipette tips.
- 2 Put the selected colonies into different 15mL centrifuge tubes, each with  5 mL LB broth.
- 3 Incubate in an orbital shaker at  171 rpm, 37°C overnight.
- 4 Pipette  30 µL culture from each tube to different agar plates. Spread the culture evenly on the plate.
- 5 Incubate the plates in a biochemical incubator at  37 °C overnight.

Note

These steps help obtain adequate (and genetically pure) colonies for testing and research in the future.


- 6 Prepare several sterilized 1.5ml microcentrifuge tubes.

Note

If you have X samples to test, prepare X+2 tubes at least. You may prepare more in case you make mistakes.


Colony PCR



5m

- 7 Add  30 µL TE buffer to X+2 1.5ml microcentrifuge tubes each. Label the tubes as "1, 2, 3, ..., X, +, -".

**Note**

"+" tube means the positive control group and "-" tube means the negative control group.

8 Pick one colony from each plate  [go to step #5](#) using a sterilized pipette tip and put the colonies into different 1.5ml microcentrifuge tubes numbered "1, 2, 3, ..., X"

9 Place the tubes in a heating block, heating at  100 °C for  00:05:00 .

5m

10 Prepare Master Mix for colony PCR. The recipe for the Master Mix is as follows:

Item	For one reaction	For Y reactions (prepare one more in case pipette inaccuracies) Y+1
dH ₂ O	11.5 µl	µl
10X PCR buffer	2 µl	µl
dNTP mix	0.4 µl	µl
forward primer	0.4 µl	µl
reverse prime	0.4 µl	µl
<i>Taq</i> polymerase	0.3 µl	µl
<i>Total</i>	<i>15 µl</i>	µl


Note

Keep all PCR reagents on ice.



If you have Y reactions, prepare Master Mix for Y+1 reactions.



This means if you have X samples, you need to prepare X+3 reactions for X samples, one positive control, one negative control and another portion in case of pipette inaccuracies.

11 Label X+2 0.2mL PCR tubes as "1, 2, 3, ..., X, +, -"

12 Pipette  15 µL from the Master Mix into all X+2 0.2 mL PCR tubes.



- 13 Pipette  5 μL from the colony lysate from tube  go to step #9 "1, 2, 3, ..., X" into corresponding 0.2 mL PCR tubes.

Add  5 μL plasmid into "+" PCR tube. Add  5 μL ddH₂O into "-" PCR tube.

Note

The plasmid used here is the plasmid transduced previously into the bacteria on the original plate at Step#1.

Protocol

NAME

Plasmid transduction using competent cell

CREATED BY

An.Huang

Preview

- 14 Place the X+2 PCR tubes in Thermocycler. PCR procedure will be set as the following programme:

1 cycle of  95 °C 2 min

30 cycles of  94 °C 30 sec

 50 °C 30 sec

 72 °C 2 min

Final extension  72 °C 10 min

- 15 When the programme is finished, store the tubes at  4 °C