

Feb 13, 2025

Version 4

Colorimetric determination of urea V.4

 PLOS One

✓ Peer-reviewed method

DOI

<https://dx.doi.org/10.17504/protocols.io.14egnmzqg5d/v4>

Noah Langenfeld¹, Laurenpayne², Bruce Bugbee²

¹Utah State University; ²Crop Physiology Laboratory, Utah State University

USU Crop Physiology La...



Noah Langenfeld

Utah State University



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egnmzqg5d/v4>

External link: <https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0259760>

Protocol Citation: Noah Langenfeld, Laurenpayne , Bruce Bugbee 2025. Colorimetric determination of urea. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egnzmzqg5d/v4> Version created by **Noah Langenfeld**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working.

Created: February 13, 2025

Last Modified: February 13, 2025

Protocol Integer ID: 120244

Keywords: urea, urea assay, colorimetric urea, urea test, thiosemicarbazide, diacetyl monoxime, colorimetric determination of urea, colorimetric determination, urea in solution, urea, protocol measures the absorbance, absorbance, mm urea, complexation with diacetyl monoxime, diacetyl monoxime

Abstract

This protocol measures the absorbance of urea in solution in complexation with diacetyl monoxime at 520 nm and is linearly proportional to concentration up to 5 mM urea.

Materials

thiosemicarbazide, diacetyl monoxime, water, sulfuric acid, phosphoric acid, ferric chloride

Safety warnings

- ⚠ This protocol utilizes strong acids. Ensure proper precautions and safety equipment are used when dealing with these chemicals.

Before start

Ensure a lab coat, lab goggles, and plastic gloves are worn throughout this assay.



Mixed Acid Reagent Preparation

1

Dissolve  2.5 mg ferric chloride (97% purity) in  45 mL deionized water in a 250 mL volumetric flask.

2

Add  80 μ L of  14.8 Mass Percent concentrated phosphoric acid (85 wt. % phosphoric acid in water).

3

Safety information

Caution: Preparation of  4.7 Mass Percent sulfuric acid is highly exothermic. Slowly add acid to  100 mL water and stir frequently to avoid spattering.

Prepare  4.7 Mass Percent sulfuric acid by diluting  65.25 mL concentrated sulfuric acid ( 18 Mass Percent , 98% purity) up to  250 mL with deionized water.

4

Dilute ferric chloride, deionized water, and phosphoric acid mixture to  250 mL with  4.7 Mass Percent sulfuric acid.

5

Mix until dissolved.

Mixed Color Reagent Preparation

6

Add  20.9 mg diacetyl monoxime (99% purity) and  52.4 mg thiosemicarbazide (99% purity) to a 250 mL volumetric flask.

7

Dilute to volume ( 250 mL) with deionized water.



8 Mix until dissolved.

Urea Assay

10m 5s

9 Prepare known concentrations of urea at 0, 1, 2, 3, 4, and 5 mM for a calibration curve.

10 Fill a 600 mL beaker with  200 mL deionized water.

11 Place the 600 mL beaker on hot plate and bring to a boil.

12 To test a sample, aliquot  1 mL into a 20 mL glass test tube.

13 Add  2 mL Mixed Acid Reagent.

14 Add  2 mL Mixed Color Reagent.

15 Ensure contents of test tube are vigorously mixed by vortexing for at least

 00:00:05 .

5s

16 Repeat steps 12-15 for each sample and standard to be analyzed.

17 Ensure test tubes are labeled and place in boiling water bath for  00:20:00 .

18 Remove tubes from water bath and let cool for  00:10:00 .

10m

19 Fill plastic cuvette to line with cooled solution.



- 20 Place sample cuvette in spectrophotometer and record absorbance at 520 nm. Blank against the standard solution with both reagents, but no urea.
- 21 Use the calibration curve to determine urea concentrations of unknown solutions.