

Oct 20, 2020

Commercial automated scRNA-seq workflow

DOI

dx.doi.org/10.17504/protocols.io.bnmvmc66



DNA Pipelines R&D¹, Peter Ellis¹, Lesley Shirley¹

¹Wellcome Sanger Institute



Lesley Shirley

Wellcome Sanger Institute

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bnmvmc66

Protocol Citation: DNA Pipelines R&D, Peter Ellis, Lesley Shirley 2020. Commercial automated scRNA-seq workflow. [protocols.io https://dx.doi.org/10.17504/protocols.io.bnmvmc66](https://dx.doi.org/10.17504/protocols.io.bnmvmc66)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 20, 2020

Last Modified: October 20, 2020

Protocol Integer ID: 43413



Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This SOP describes the procedure for plate based scRNA-seq performed with a commercial available kit from New England BioLabs. Following library construction, samples are pooled in equivolume and quantified, prior to sequencing on the Illumina HiSeq 4000 platform.

Guidelines

It is vital all steps prior to cDNA amplification are performed in a designated RNase free and pre-cDNA amplification laboratory.


Note


Throughout the protocol we have indicated the liquid handling in use at Sanger for specific parts of the process. However, these steps can be performed on alternative liquid handlers.

Materials

MATERIALS

 KAPA HiFi HotStart ReadyMixPCR Kit **Kapa Biosystems Catalog #KK2602**


 NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns **New England Biolabs Catalog #E6420L**

 AccuClear® Ultra High Sensitivity dsDNA Quantitation Kit with DNA Standards **Biotium Catalog ##31028**


 Bioanalyzer chips and reagents (DNA 1000) **Agilent Technologies Catalog #5067-1504**


 AMPure XP Beads **Beckman Coulter Catalog #A63882**

STEP MATERIALS

 NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns **New England Biolabs Catalog #E6420L**

 AMPure XP Beads **Beckman Coulter Catalog #A63882**

 AccuClear® Ultra High Sensitivity dsDNA Quantitation Kit with DNA Standards **Biotium Catalog ##31028**


 NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns **New England Biolabs Catalog #E6420L**

 AMPure XP Beads **Beckman Coulter Catalog #A63882**

 AMPure XP Beads **Beckman Coulter Catalog #A63882**

 2x Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2602**

 AMPure XP Beads **Beckman Coulter Catalog #A63882**

 Bioanalyzer chips and reagents (DNA 1000) **Agilent Technologies Catalog #5067-1504**

Protocol materials

-  NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns **New England Biolabs Catalog #E6420L**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  AccuClear® Ultra High Sensitivity dsDNA Quantitation Kit with DNA Standards **Biotium Catalog ##31028**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  KAPA HiFi HotStart ReadyMixPCR Kit **Kapa Biosystems Catalog #KK2602**
-  2x Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2602**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns **New England Biolabs Catalog #E6420L**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  Bioanalyzer chips and reagents (DNA 1000) **Agilent Technologies Catalog #5067-1504**
-  AccuClear® Ultra High Sensitivity dsDNA Quantitation Kit with DNA Standards **Biotium Catalog ##31028**
-  NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns **New England Biolabs Catalog #E6420L**
-  Bioanalyzer chips and reagents (DNA 1000) **Agilent Technologies Catalog #5067-1504**
-  NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns **New England Biolabs Catalog #E6420L**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  AccuClear® Ultra High Sensitivity dsDNA Quantitation Kit with DNA Standards **Biotium Catalog ##31028**
-  NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns **New England Biolabs Catalog #E6420L**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  2x Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2602**
-  AMPure XP Beads **Beckman Coulter Catalog #A63882**
-  Bioanalyzer chips and reagents (DNA 1000) **Agilent Technologies Catalog #5067-1504**

Preparation of lysis buffer plates and FACS

- 1 **Important!** This step must be performed in a designated RNase free and pre-cDNA amplification area, keeping reagents chilled at all times

- 2



NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96

rxns

New England Biolabs Catalog #E6420L


Prepare the cell lysis buffer, which will provide sufficient volume for one 96-well plate



On ice



| Reagent | Volume (μl) |
|---------------------------------|-------------|
| NEBNext Cell Lysis Buffer (10x) | 24 |
| Murine RNase Inhibitor | 12 |
| Nuclease-Free Water | 204 |
| Total | 240 |

Mix well by pipetting.

- 3 Use the Formulatrix Mantis microfluidic liquid handler to dispense  2 μL of lysis buffer into a 96-well PCR plate.

Note

If required add a diluted stock (1/500,000) of ERCCs into the lysis buffer.


- 4 Seal dispensed plates, centrifuge immediately  1000 x g, 4°C, 00:01:00 and keep chilled on ice.
- 5 **PAUSE POINT** Lysis buffer plates can be stored at  -80 °C prior to cell sorting. Plates can be stored for < 6 months.





6

Note

When FACS sorting, take care of plate calibration/priming prior to single-cell deposition. If many plates are deposited in parallel, repeat the calibration/priming at least every 8 plates.


Defrost lysis buffer plates prior to cell sorting, centrifuge  1000 x g, 4°C, 00:01:00 and keep chilled on ice.

7 After FACS sorting, seal and centrifuge plates immediately  1000 x g, 4°C, 00:01:00 and keep chilled on ice.

8 **PAUSE POINT** Plates of sorted cells can be stored at  -80 °C for < 6 months. the quality of the data depends on the cell type and duration of storage.



Primer Annealing for first-strand synthesis


9

Prepare the primer annealing mix, which will provide sufficient volume for one 96-well plate  On ice

| Reagent | Volume (μl) |
|-----------------------------------|-------------|
| NEBNext Single Cell RT Primer Mix | 50 |
| Nuclease-Free Water | 250 |
| Total | 300 |

Mix well by pipetting.

10 The Agilent Bravo with 96 ST head will combine  1.6 μL of primer annealing mix with  2 μL of lysed cells and mix by pipetting.

11 Seal and transfer the plate to a thermocycler with the heated lid set to  100 °C and run the following program:

| Temperature | Time |
|-------------|------|
|-------------|------|



| | |
|------|-----------|
| 70°C | 5 minutes |
| 4°C | ∞ |

Reverse transcription (RT) and template switching

- 12 Prepare the RT mix, which will provide sufficient volume for one 96-well plate

On ice

| Reagent | Volume (μl) |
|----------------------------------|-------------|
| NEBNext Single Cell RT Buffer | 250 |
| NEBNext Template Switching Oligo | 50 |
| NEBNext Single Cell RT Enzyme | 150 |
| Nuclease-Free Water | 100 |
| Total | 550 |

Mix well by pipetting.

- 13 The Bravo will add 4.4 μL of RT mix to each sample and mix by pipetting.

- 14 Seal and transfer the plate to a thermocycler with the heated lid set to 100 °C and run the following program:

| Temperature | Time |
|-------------|------------|
| 42°C | 90 minutes |
| 70°C | 10 minutes |
| 4°C | ∞ |



- 15 Prepare the cDNA amplification mix, which will provide sufficient volume for one 96-well plate On ice

| Reagent | Volume (μl) |
|---|-------------|
| NEBNext Single Cell cDNA PCR Master Mix | 2500 |



| | |
|-------------------------------------|-------------|
| NEBNext Single Cell cDNA Primer Mix | 100 |
| Nuclease-Free Water | 1400 |
| Total | 4000 |


Mix well by pipetting.

- 16 The Bravo will add  32 μL of cDNA amplification mix to each sample and mix by pipetting.
- 17 Seal and transfer the plate to a thermocycler with the heated lid set to  100 °C and run the following program:

| Temperature | Time | Cycles |
|-------------|------------|------------------------------|
| 98°C | 45 seconds | 1 |
| 98°C | 10 seconds | 16-25 depending on cell type |
| 62°C | 15 seconds | |
| 72°C | 3 minutes | |
| 72°C | 5 minutes | 1 |
| 4°C | ∞ | 1 |

Purification of amplified cDNA



18

 AMPure XP Beads **Beckman Coulter Catalog #A63882**











Allow AMPure XP beads to equilibrate to room temperature (~30 minutes). Ensure solution is homogenous prior to use, mixing gently by inversion.

- 19 Centrifuged amplified cDNA plate  1000 x g, 00:01:00

- 20 Use the Agilent Bravo with a 96 LT multichannel head to perform the following steps:

- 20.1 Add 0.6 X volume of SPRI beads per sample ( 24 μL SPRI :  40 μL amplified cDNA), mix well by pipetting.



- 20.2 Incubate for  00:05:00 at  Room temperature
- 20.3 Transfer the plate to the magnet, allow  00:02:00 for the beads to settle.
- 20.4 Carefully remove and discard the supernatant without disturbing the bead pellet.
- 20.5 Wash the beads with  180 μL 80% freshly prepared ethanol for  00:00:30 remove ethanol and discard.
- 20.6 Repeat ethanol wash.
- 20.7 Allow beads to dry  00:05:00
- 20.8 Remove the plate from the magnet, add  50 μL nuclease-free water and resuspend by mixing well.
- 20.9 Incubate for  00:02:00 at  Room temperature
- 20.10 Transfer the plate to the magnet, allow  00:05:00 for the beads to settle.
- 20.11 Transfer supernatant into a new 96-well PCR plate, taking care not to disturb the pellet.

Quality control of amplified cDNA

21







Note

Purified amplified cDNA is quantified with a fluorescence based assay. We use the AccuClear Ultra High Sensitivity dsDNA Quantitation kit with 7 DNA standards (Biotium) according to manufacturer's instructions.

To streamline the workflow, we do not normalise sample input for library preparation. Instead, we calculate an average concentration and transfer a fixed volume such that 5-25 ng of each successfully amplified cDNA sample enters library preparation.



AccuClear® Ultra High Sensitivity dsDNA Quantitation Kit with DNA Standards **Biotium Catalog ##31028**

- 22 Pipette  20 μL of each DNA standard into wells A1 - G1 of a 96-well PCR plate. Add nuclease-free water to H1.
- 23 Dilute the AccuClear dye (100X) to working concentration by mixing  300 μL dye with  30 mL AccuClear buffer in a 50 ml Falcon. Mix thoroughly by vortexing and transfer to a 384-well reservoir.
- 24 Use the SPT Labtech Mosquito LV to stamp  200 nL of amplified cDNA and  1 μL of known standards in triplicate into a 384-well assay plate. Immediately proceed to the next step.
- 25 Use the Agilent Bravo with a 384ST multichannel head to add  50 μL 1 X AccuClear dye from the reservoir to the assay plate, mix thoroughly by pipetting.
- 26 Measure fluorescence values on a BMG FLUOstar Omega plate reader calibrated for use with AccuClear dye.
- 27 Confirm known standards are performing as expected.
- 28 Dilute any samples >125 ng/ μL with nuclease free water so they are in the range of 10 - 125 ng/ μL and repeat quantitation.

**Note**

We use 5X the volume of standard vs sample in our assay setup, which should allow a quantitative range of 0.15 ng/μl - 125 ng/μl. This deviates from the standard kit SOP which has a stated range of 0.03 ng/μl - 25 ng/μl.

- 29 Taking an average across the plate. Transfer ~ 10 ng of cDNA into a new 96-well PCR plate for sequencing library preparation.
- 30 **PAUSE POINT** Purified amplified cDNA can be stored at -20 °C for several weeks prior to library preparation.

Illumina sequencing library preparation

31

Note

We use the NEBNext® Single Cell/Low Input RNA Library Prep Kit for Illumina for library preparation, which we have automated on the Agilent Bravo NGS platform with some modifications. We use a custom adapter set, however any TruSeq adapters are suitable.















NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina - 96 rxns
New England Biolabs Catalog #E6420L



AMPure XP Beads **Beckman Coulter Catalog #A63882**


- 32 Allow AMPure XP beads to equilibrate to room temperature (~30 minutes). Ensure solution is homogenous prior to use, mixing gently by inversion.
- 33 Centrifuged plate containing ~ 10 ng of purified cDNA 1000 x g, 00:01:00
- 34 Use the Agilent Bravo with a 96 LT multichannel head to perform the following steps:
- 34.1 Add 0.9 X volume of SPRI beads per sample, mix well by pipetting.




- 34.2 Incubate for  00:05:00 at  Room temperature
- 34.3 Transfer the plate to the magnet, allow  00:02:00 for the beads to settle.
- 34.4 Carefully remove and discard the supernatant without disturbing the bead pellet.
- 34.5 Wash the beads with  180 μL 80% freshly prepared ethanol for  00:00:30 remove ethanol and discard.
- 34.6 Repeat ethanol wash.
- 34.7 Allow beads to dry  00:05:00
- 34.8 Remove the plate from the magnet, add  13 μL TE pH 8.0 and resuspend by mixing well.
- 34.9 Incubate for  00:02:00 at  Room temperature
- 34.10 Transfer the plate to the magnet, allow  00:05:00 for the beads to settle.
- 34.11 Transfer  12.4 μL into a new 96-well PCR plate, taking care not to disturb the pellet.
- 35 Prepare fragmentation/end prep mix, which will provide sufficient volume for one 96-well plate  On ice

| Reagent | Volume (μL) |
|-------------------------------------|--------------------------|
| NEBNext Ultra II FS Reaction Buffer | 336 |
| NEBNext Ultra II FS Reaction Enzyme | 96 |
| Total | 432 |

Mix well by pipetting.

- 36 The Bravo will add  3.6 µL of fragmentation/end prep mix to each sample and mix by pipetting.

- 37 Seal and transfer the plate to a thermocycler with the heated lid set to  100 °C and run the following program:

| Temperature | Time |
|-------------|------------|
| 72°C | 15 minutes |
| 65°C | 30 minutes |
| 4°C | ∞ |



- 38 Prepare adapter ligation mix, which will provide sufficient volume for one 96-well plate

 On ice

| Reagent | Volume (µl) |
|--------------------------------------|-------------|
| NEBNext Ultra II Ligation Master Mix | 1440 |
| NEBNext Ultra II Ligation Enhancer | 48 |
| TruSeq Duplexed Adapter (100 µM) | 12 |
| Nuclease-Free Water | 108 |
| Total | 1608 |

Mix well by pipetting.

- 39 The Bravo will add  13.4 µL of ligation mix to each sample and mix by pipetting.


- 40 The plate is incubated on deck at  20 °C for  00:15:00 , however this step may also be performed on a thermocycler.

41



Note



We use alternative TruSeq compatible adapters, which do not require the USER enzyme incubation step. If using NEBNext adapters, follow the steps in the NEB protocol to add USER enzyme to the ligation reaction.


42

 AMPure XP Beads **Beckman Coulter Catalog #A63882**



Allow AMPure XP beads to equilibrate to room temperature (~30 minutes). Ensure solution is homogenous prior to use, mixing gently by inversion.

42.1 Add 0.7 X volume of SPRI beads per sample ( 20 μL SPRI :  29.4 μL amplified cDNA), mix well by pipetting.


42.2 Incubate for  00:05:00 at  Room temperature


42.3 Transfer the plate to the magnet, allow  00:02:00 for the beads to settle.



42.4 Carefully remove and discard the supernatant without disturbing the bead pellet.


42.5 Wash the beads with  180 μL 80% freshly prepared ethanol for  00:00:30 remove ethanol and discard.

42.6 Repeat ethanol wash.

42.7 Allow beads to dry  00:05:00

42.8 Remove the plate from the magnet, add  25 μL nuclease-free water and resuspend by mixing well.

42.9 Incubate for  00:02:00 at  Room temperature

42.10 Transfer the plate to the magnet, allow  00:05:00 for the beads to settle.

42.11 Transfer supernatant into a new 96-well PCR plate, taking care not to disturb the pellet.

43

Note



We use KAPA HiFi HotStart ReadyMix and unique dual indexed (UDI) tag plates for library PCR.


Note: this deviates from the standard NEB protocol which uses NEBNext Ultra II Q5 Master Mix and different cycling conditions.

 2x Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2602**


44 Prepare PCR mix, which will provide sufficient volume for one 96-well plate  On ice

| Reagent | Volume (μl) |
|-----------------------------|-------------|
| KAPA HiFi HotStart ReadyMix | 3000 |
| Total | 3000 |












45 The Bravo will add  25 μL PCR mix and  25 μL sample into a lyophilised plate of UDIs and mix thoroughly by pipetting. The final concentration of each UDI in the PCR reaction is 2 μM.





46 Seal and transfer the plate to a thermocycler with the heated lid set to  100 °C and run the following program:

| Temperature | Time | Cycles |
|-------------|------------|--------|
| 98°C | 45 seconds | 1 |
| 98°C | 10 seconds | 8 |
| 62°C | 15 seconds | |
| 72°C | 3 minutes | |
| 72°C | 5 minutes | 1 |
| 4°C | ∞ | 1 |

- 47 **PAUSE POINT** amplified libraries can be stored at  -20 °C for several weeks prior to library purification.

Pooling and purification of amplified libraries

- 48 In a post-PCR lab, use the Hamilton STAR or Beckman NX-8 to combine  5 µL of each sample per plate to form an equivolume pool of 96 samples.
- 49  AMPure XP Beads **Beckman Coulter Catalog #A63882**
Allow AMPure XP beads to equilibrate to room temperature (~30 minutes). Ensure solution is homogenous prior to use, mixing gently by inversion.
- 49.1 Manually transfer  400 µL of the equivolume pool into a 1.5 ml Eppendorf tube
- 49.2 Add 0.95 X volume of SPRI beads ( 380 µL SPRI :  400 µL amplified libraries), mix well by pipetting.
- 49.3 Incubate for  00:05:00 at  Room temperature
- 49.4 Transfer the tube to a magnet, allow  00:05:00 for the beads to form a pellet.
- 49.5 Carefully remove and discard the supernatant without disturbing the bead pellet.
- 49.6 Wash the beads with  180 µL 80% freshly prepared ethanol for  00:00:30 remove ethanol and discard.
- 49.7 Repeat ethanol wash.
- 49.8 Allow beads to dry  00:05:00

- 49.9 Remove the tube from the magnet, add  400 µL nuclease-free water and resuspend by mixing well.
- 49.10 Incubate for  00:02:00 at  Room temperature
- 49.11 Transfer tube to magnet, allow  00:05:00 for the beads to form a pellet.
- 49.12 Transfer supernatant into a new tube, taking care not to disturb the pellet.

Quality control and normilisation of sequencing libraries

50

Note

Library sequencing pools are quantified on an Agilent Bioanalyzer. Pools are then diluted to 2.8 nM for sequencing.


Equipment

| | |
|---|----------------|
| Bioanalyzer | NAME |
| Bioanalyzer | TYPE |
| Agilent | BRAND |
| G2991AA | SKU |
| https://www.agilent.com/en/product/bioanalyzer-automated-electrophoresis/bioanalyzer-instrument/2100-bioanalyzer-instrument-228250 | LINK |
| Any bioanalyzer will suffice. | SPECIFICATIONS |





Bioanalyzer chips and reagents (DNA 1000) **Agilent Technologies Catalog #5067-1504**

- 51 Run  1 μL of the library pool in triplicate on a Bioanalyzer using the DNA 1000 kit.
- 52 Taking an average of the readings add nuclease-free water to the library pool to produce a final concentration of 2.8 nM.

Sequencing

53

Note

We sequence samples on an Illumina HiSeq 4000 instrument (paired-end, 75-bp reads) according to the manufacturer's protocol. We typically aim for an average depth of 1 million reads per cell, plexing up to 384 samples per run.