




Jan 23, 2024

🌐 CompDuplex: Accurate detection of somatic mutations by duplex-seq with comprehensive genome coverage

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx3x4og8j/v1>



Muchun Niu¹, Chenghang Chuck Zong¹

¹Department of Molecular and Human Genetics, Baylor College of Medicine

SMAHT



Muchun Niu

Baylor College of Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx3x4og8j/v1>

Protocol Citation: Muchun Niu, Chenghang Chuck Zong 2024. CompDuplex: Accurate detection of somatic mutations by duplex-seq with comprehensive genome coverage. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx3x4og8j/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 17, 2024

Last Modified: January 23, 2024

Protocol Integer ID: 93708

Keywords: Duplex-seq, Somatic mutation, Accurate mutation calling, comprehensive genome coverage somatic mutation, seq with comprehensive genome coverage somatic mutation, comprehensive genome coverage of compduplex, accurate detection of somatic mutation, comprehensive genome coverage, method with comprehensive genome coverage, whole genome characterization of somatic mosaicism, accuracy of somatic mutation, somatic mutation, limited genome coverage, whole genome characterization, sequencing method, whole genome, sequencing, human genome, genome, sequencing platform, whole genome scale, end sequencing, streamlined chemistry of compduplex assay, complementary strands of dna, compduplex assay, dna, somatic mosaicism

Funders Acknowledgements:

SMaHT Network (NIH)

Grant ID: UG3NS132132

Disclaimer

Patent application covering CompDuplex chemistry has been filed by Baylor College of Medicine.

Abstract

Somatic mutations continuously accumulate in the human genome, posing vulnerabilities towards aging and increased risk of various diseases. However, accurate detection of somatic mutations at the whole genome scale is still challenging. By tagging and independently sequencing the two complementary strands of DNA, the recent development of duplex-sequencing methods has greatly improved the detection accuracy, however, the limited genome coverage and the compromised compatibility with existing sequencing platforms have constrained the broad applications of these methods.

To overcome these technical challenges, here we developed a duplex sequencing method with comprehensive genome coverage, which we refer to as CompDuplex-seq. The streamlined chemistry of CompDuplex assay allows efficient generation of libraries readily compatible with standard Illumina 2×150 paired-end sequencing. In addition, we validated the accuracy of somatic mutation calling and comprehensive genome coverage of CompDuplex by profiling a single-cell expanded clone. To summarize, CompDuplex chemistry supports genome-wide coverage while maintaining high accuracy, which we believe will facilitate the whole genome characterization of somatic mosaicism in various biological systems.



Protocol materials

- ⊠ Ampure XP beads **Beckman Catalog #A63881**
- ⊠ Tn5 transposase, unloaded **Diagenode Catalog #C01070010-20**
- ⊠ Glycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516**
- ⊠ Tagmentation Buffer 1 **Illumina, Inc. Catalog #20015171**
- ⊠ Water **Invitrogen - Thermo Fisher Catalog #46-2224**
- ⊠ NEBNext Ultrall Q5 Master Mix **New England Biolabs Catalog #M0544X**
- ⊠ rCutSmart Buffer **New England Biolabs Catalog #B6004S**
- ⊠ BtgZI - 500 units **New England Biolabs Catalog #R0703L**
- ⊠ 0.5M EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889-100ML**
- ⊠ 10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #93443-100ML**
- ⊠ 1M Tris-HCl, pH 8.0 **Invitrogen - Thermo Fisher Catalog #15568025**
- ⊠ Protease **Qiagen Catalog #19155**
- ⊠ T4 Polynucleotide Kinase Reaction Buffer **New England Biolabs Catalog #B0201S**
- ⊠ StickTogether DNA Ligase Buffer **New England Biolabs Catalog #B0535S**
- ⊠ T7 DNA Ligase - 750,000 units **New England Biolabs Catalog #M0318L**
- ⊠ iTaq Universal SYBR Green Supermix **Bio-Rad Laboratories Catalog #172-5112**
- ⊠ 5M NaCl **Invitrogen - Thermo Fisher Catalog #AM9760G**
- ⊠ Water **Invitrogen - Thermo Fisher Catalog #46-2224**



Custom Tn5 transposase assembly

1h

1

2m

Prepare 500 μL of **2X Annealing Buffer**.

40 μL 1M Tris-HCl, pH 8.0 Invitrogen - Thermo Fisher Catalog #15568025

10 μL 5M NaCl Invitrogen - Thermo Fisher Catalog #AM9760G

450 μL Water Invitrogen - Thermo Fisher Catalog #46-2224

2 **Transposon annealing.**

27m

Prepare the following mix in a PCR tube.

20 μL **2X Annealing Buffer**

10 μL [M] 200 micromolar (μM) ME_REV

10 μL [M] 200 micromolar (μM) T_BtgZ1_ME

Oligonucleotides sequences

ME_REV: /5Phos/CTGTCTCTTATACACATCT

T_BtgZ1_ME: ATGTGTGGAGCGATG AGATGTGTATAAGAGACAG

Mix up the reaction, spin down, and perform the following reactions on a thermal cycler.

95 °C 5 min

65 °C 5 min, Ramp rate 0.1°C/s

4 °C Hold, Ramp rate 0.1°C/s

This results in a [M] 50 micromolar (μM) transposon mix .

3 **Custom Tn5 assembly.**

31m

Prepare the following mix in a 1.5 mL tube.

20 μL Tn5 transposase, unloaded Diagenode Catalog #C01070010-20

20 μL [M] 50 micromolar (μM) transposon mix



Gently pipet to mix, and incubate at 23°C for 30 min.

Add $20\ \mu\text{L}$ Glycerol Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516

to the tube, mix well, aliquot, label as "Tn5_BtgZ1", and store at -20°C .

CompDuplex Procedure

4 1. Re-purification of genomic DNA

Dilute $20\ \text{ng}$ Sample to $10\ \mu\text{L}$ nucleus-free water in a PCR tube. Add $5\ \mu\text{L}$ Ampure XP beads Beckman Catalog #A63881 (0.5X) to the tube, mix well, and incubate at Room temperature for 5 min.

Put the tube on the magnetic stand. When the aqueous phase becomes transparent, remove the supernatant, wash twice with $120\ \mu\text{L}$ 80 % volume Ethanol. Air-dry the beads.

5 2. On-bead custom Tn5 tagmentation

Prepare the following **tagmentation mix**.

$2\ \mu\text{L}$ Tagmentation Buffer 1 Illumina, Inc. Catalog #20015171

$0.48\ \mu\text{L}$ 1:30 diluted Tn5_BtgZ1

$7.52\ \mu\text{L}$ Water Invitrogen - Thermo Fisher Catalog #46-2224

Resuspend the air-dried beads with $10\ \mu\text{L}$ **tagmentation mix**.

Incubate on a thermal cycler using the following program: 20°C for 5 min, 55°C for 15 min.

Quench the reaction by adding $3\ \mu\text{L}$ of 0.2M EDTA.

Incubate on a thermal cycler at 50°C for 30 min.



6 3. Gap filling

Add  3 μL of 0.2M MgCl_2 to the tube.

Add  16 μL

 NEBNext Ultrall Q5 Master Mix **New England Biolabs Catalog #M0544X** to the tube.

Incubate on a thermal cycler at  65 $^{\circ}\text{C}$ for 3 min.

Quench the reaction by adding  2 μL of 0.5M EDTA.

Perform double size selection with 0.48X/0.75X

 Ampure XP beads **Beckman Catalog #A63881**, elute in  7.5 μL water.

7 4. Restriction enzyme digestion

Prepare the following restriction enzyme digestion mix (Total  15 μL).

 7.5 μL Gap filling product

 1.5 μL  rCutSmart Buffer **New England Biolabs Catalog #B6004S**

 0.7 μL  BtgZI - 500 units **New England Biolabs Catalog #R0703L**

 5.3 μL Water

Incubate on a thermal cycler with the following program:  25 $^{\circ}\text{C}$ for 10 min,

 37 $^{\circ}\text{C}$ for 3 h,  10 $^{\circ}\text{C}$ for 12 h.

8 5. Restriction enzyme release

Prepare the following **restriction enzyme release mix** (Total  500 μL).

 15 μL  1M Tris-HCl, pH 8.0 **Invitrogen - Thermo Fisher Catalog #15568025**

 12 μL

 0.5M EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889-100ML**

 10 μL 1M KCl

 5 μL

 10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #93443-100ML**



🧴 18 μ L 10% NP40

🧴 420 μ L Water

🧴 20 μ L [M] 20 mg/mL Protease Qiagen Catalog #19155

Add 🧴 15 μ L **restriction enzyme release mix** to the sample.

Incubate on a thermal cycler with the following program: 🌡️ 50 °C for 45 min,

🌡️ 25 °C for 2 h.

Purify with 0.85X Ampure XP beads Beckman Catalog #A63881 , elute in

🧴 10 μ L water.

9 6. Y-shape ligation adapter annealing

Prepare the following [M] 15 micromolar (μ M) T1T2_TCTT ligation adapter .

🧴 1 μ L

T4 Polynucleotide Kinase Reaction Buffer New England Biolabs Catalog #B0201S

🧴 10.5 μ L [M] 100 micromolar (μ M) **T1_BtgZ1_ME**

🧴 7.5 μ L [M] 100 micromolar (μ M) **Rev_T2_BtgZ1**

🧴 31 μ L Water.

Oligonucleotides sequences

T1_BtgZ1_ME: TCGTCGGCAGCGTC AGATGTGTAT

Rev_T2_BtgZ1: /5Phos/ TCTT ATACATCT CCGAGCCCACGAGAC

Mix up the reaction, spin down, and perform the following reactions on a thermal cycler.

🌡️ 65 °C 5 min

🌡️ 20 °C Hold, Ramp rate 0.1°C/s

10 7. Y-shape ligation adapter ligation

Prepare the following ligation mix (Total 🧴 50 μ L).



🧴 10 µL Sample from Step 8

🧴 25 µL

🧴 StickTogether DNA Ligase Buffer **New England Biolabs Catalog #B0535S**

🧴 3 µL [M] 15 micromolar (µM) T1T2_TCTT ligation adapter

🧴 3 µL 🧴 T7 DNA Ligase - 750,000 units **New England Biolabs Catalog #M0318L**

🧴 9 µL Water

Incubate at 🌡 25 °C for 30 min.

Purify with 0.8X 🧴 Ampure XP beads **Beckman Catalog #A63881**, elute in

🧴 51 µL water.

11 **8. Quantification of library complexity.**

Dilute 🧴 1 µL of ligation product with 🧴 9 µL water. Use the 1:10 diluted product to quantify library complexity using qPCR. Any Illumina Nextera libraries with known concentration can be used as a standard. An example qPCR procedure is shown below.

🧴 5 µL

🧴 iTaq Universal SYBR Green Supermix **Bio-Rad Laboratories Catalog #172-5112**

🧴 0.25 µL [M] 10 micromolar (µM) **T1ME**

🧴 0.25 µL [M] 10 micromolar (µM) **T2ME**

🧴 3.5 µL water

🧴 1 µL template

Oligonucleotides sequences

T1ME: TCGTCGGCAGCGTC AGATGTGTATAAGAGACAG

T2ME: GTCTCGTGGGCTCGG AGATGTGTATAAGAGACAG

🌡 94 °C 2 min

30 cycles of

🌡 94 °C 20s

🌡 68 °C 20s

🌡 72 °C 1 min, fluorescence scanning

Melting curve



12 **9. Library amplification**

We recommend to amplify 14 cycles for a library with 20 million DNA fragments.

🧪 25 μ L

🧬 NEBNext Ultrall Q5 Master Mix **New England Biolabs Catalog #M0544X**

🧪 2.5 μ L [M] 10 micromolar (μ M) **Illumina Nextera N5XX index primer**

🧪 2.5 μ L [M] 10 micromolar (μ M) **Illumina Nextera N7XX index primer**

🧪 20 μ L ligation product and water

Oligonucleotides sequences

Illumina Nextera N5XX index primer: AATGATACGGCGACCACCGAGATCTACAC
NNNNNNNN TCGTCGGCAGCGTC

Illumina Nextera N7XX index primer: CAAGCAGAAGACGGCATACGAGAT
NNNNNNNN AGATGTGTATAAGAGACAG

PCR cycles:

🌡️ 98 °C 30s

14 cycles of

🌡️ 98 °C 10s

🌡️ 63 °C 30s

🌡️ 72 °C 1 min

🌡️ 72 °C 3 min

Purify with 0.8X 🧬 Ampure XP beads **Beckman Catalog #A63881** , elute in

🧪 20 μ L water.