



Mar 31, 2023

Complex I activity assay



In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l27r8qg1y/v1>

Maria Jose Perez J.¹, michela.deleidi¹

¹German Center for Neurodegenerative Diseases (DZNE), Tübingen, 72076 Germany



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l27r8qg1y/v1>

Protocol Citation: Maria Jose Perez J., michela.deleidi 2023. Complex I activity assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l27r8qg1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 28, 2023

Last Modified: June 01, 2024

Protocol Integer ID: 79543

Keywords: Complex I activity assay, - Extinction coefficient, Enzyme activity, ASAPCRN, assay, complex, activity, protocol

Abstract

This protocol describes the complex I activity assay.

Attachments



404-877.docx

15KB

Guidelines

Reference for analysis: <https://www.nature.com/articles/nprot.2012.058>

Materials

KIT:

Mitocheck Complex I activity assay kit

 MitoCheck® Complex I Activity Assay Kit **Cayman Chemical Company Catalog #700930**

 Qproteome Mitochondria Isolation Kit **Qiagen Catalog #37612**

Before start

- All assays are carried out at 25 °C .
- After mitochondrial isolation (Qproteome Mitochondria Isolation Kit. QIAGEN Cat. No. / ID: 37612), resuspend the final pellet in 50 µL of storage buffer, keep isolated mitochondria On ice .
- In a ventilated hood, weigh out 6.5 mg of KCN and dissolve in 1 mL of 0.1 Mass Percent NaOH (stock solution of 100 millimolar (mM))
- Label two polystyrene tubes as A and B. For 20 reactions prepare:

A	B
Tube A (1 ml)	Tube B (675 µl)
910 µl of Complex I activity buffer	625 µl of Complex I activity buffer
20 µl of 100mM KCN (1 mM)	30 µl of NADH assay reagent
50 µl FF-BSA Assay Reagent	20 µl of Ubiquinone assay reagent
20 µl of Vehicle	



Protocol

5m

- 1 Distribute the contents of tube A and B in strips suitable for multichannel use.
- 2 In a Half Volume 96-well clear plate add  50 μL of the contents of tube A to each well. 
- 3 Add  20 μL of sample to each well. 
- 4 Place plate in plate reader and add  30 μL of B to each well. 
- 5 Immediately measure absorbance at 340 nm in kinetic read mode (30 seconds intervals for  00:05:00 at  25 $^{\circ}\text{C}$) 5m

Calculations

- 6 The specific activity of complex I is calculated as $\text{nmol min}^{-1} \text{mg}^{-1}$ of protein according to the following equation:

Enzyme activity ($\text{nmol min}^{-1} \text{mg}^{-1}$) = $(\Delta \text{ Absorbance/min} \times 1,000) / [(\text{extinction coefficient} \times \text{volume of sample used in ml}) \times (\text{sample protein concentration in mg ml}^{-1})]$.
- 7 Extinction coefficient for NADH $6.2 \text{ mM}^{-1} \text{ cm}^{-1}$.