



Mar 21, 2025

Confocal Microscopy of DFP-Treated Cells

 [Nature Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj9e9blk5/v1>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj9e9blk5/v1>

External link: <https://doi.org/10.1038/s41556-025-01712-y>

Protocol Citation: Elias Adriaenssens 2025. Confocal Microscopy of DFP-Treated Cells. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n2bvj9e9blk5/v1>

Manuscript citation:

Adriaenssens, E., Schaar, S., Cook, A.S.I. *et al.* Reconstitution of BNIP3/NIX-mitophagy initiation reveals hierarchical flexibility of the autophagy machinery. *Nat Cell Biol* **27**, 1272–1287 (2025). <https://doi.org/10.1038/s41556-025-01712-y>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 05, 2025

Last Modified: March 21, 2025

Protocol Integer ID: 120015

Keywords: ASAPCRN, confocal microscopy of dfp, confocal microscopy, cells after fixation, treated cell, treated cells this protocol detail, dfp, microscopy, cell

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

This protocol details the confocal microscopy of DFP-Treated cells after fixation and staining.

Materials

-  Bovine Serum Albumin [BSA] Fisher Scientific Catalog #9048-46-8
-  DAPI Fluoromount-G® Southern Biotech Catalog #0100-20



Procedure

1d 2h 20m

- 1 Seed the cells on glass coverslips (12 mm #1.5) at a concentration of 120,000 cells/well, and after treatment with DFP for ⌚ 24:00:00 , fix in 100% methanol for ⌚ 00:20:00 at 🌡️ -20 °C .
- 2 After washing with PBS, block the samples with 5% (v/v) BSA (9048-46-8, Sigma-Aldrich) diluted in PBS for ⌚ 01:00:00 at 🌡️ Room temperature .
- 3 Dilute the primary and secondary antibodies in 5% BSA and incubate for ⌚ 01:00:00 at 🌡️ Room temperature with three PBS washing steps in between.
- 4 Mount the cells on microscopy slides in DAPI Fluoromount-G mounting medium (0100-20, Southern Biotech), which stains the nuclei, add a cover glass, and store at 🌡️ 4 °C until use.
- 5 Perform the confocal microscopy with a Zeiss LSM700 laser scanning confocal microscopy with Plan-Apochromat 40×/1.30 Oil DIC, WD 0.21 mm objective or Zeiss LSM900 with Airyscan 2 confocal super-resolution microscopy with Plan-Apochromat 63×/1.4 Oil DIC, WD 0.19 mm objective.

1d 0h 20m



1h



1h

