



Aug 26, 2025

🌐 Construction of protein-based artificial kinetochore in mouse oocytes

📖 [Nature Cell Biology](#)

DOI

dx.doi.org/10.17504/protocols.io.e6nvw41k9lmk/v1

Yuanzhuo Zhou^{1,2}, Kohei Asai^{1,2}, Hirohisa Kyogoku^{1,3}, Tomoya S. Kitajima^{1,2}

¹Laboratory for Chromosome Segregation, RIKEN Center for Biosystems Dynamics Research (BDR), Kobe, Japan;

²Graduate School of Biostudies, Kyoto University, Kyoto, Japan;

³Graduate School of Agricultural Science, Kobe University, Kobe, Japan

Designing protein-base...



Yuanzhuo Zhou

Laboratory for Chromosome Segregation, RIKEN Center for Bios...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw41k9lmk/v1>

External link: <https://doi.org/10.1038/s41556-025-01792-w>



Protocol Citation: Yuanzhuo Zhou, Kohei Asai, Hirohisa Kyogoku, Tomoya S. Kitajima 2025. Construction of protein-based artificial kinetochore in mouse oocytes. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvw41k9lmk/v1>

Manuscript citation:

Zhou, Y., Asai, K., Kyogoku, H. *et al.* Designing protein-based artificial kinetochores as decoys to prevent meiotic errors in oocytes. *Nat Cell Biol* (2025). <https://doi.org/10.1038/s41556-025-01792-w>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 23, 2025

Last Modified: August 26, 2025

Protocol Integer ID: 225314

Keywords: artificial kinetochore in mouse oocyte, based artificial kinetochore, mouse oocyte, construction of protein, protein, mouse

Abstract

This protocol describes the workflow of construction of protein-based artificial kinetochore in mouse oocytes.

Materials

Pregnant mare serum gonadotropin (PMSG); mice (8–12 weeks old); M2 medium; IBMX (Sigma); mMESSAGE mMACHINE T7 kit (Thermo Fisher Scientific); mRNAs for NDC80-FKBP, NUF2, H2B-mCherry, hAG-Pfv-Sapphire-GS-FRB, FKBP-GS-FKBP; rapamycin (Funakoshi).

Troubleshooting



- 1 Inject mice (8–12 weeks old) with 0.2 ml (10 IU) pregnant mare serum gonadotropin (PMSG).
- 2 Collect fully grown GV-stage oocytes 48 hours after PMSG injection.
- 3 Incubate oocytes in M2 medium containing 200 μ M IBMX (Sigma) at 37 °C.
- 4 Transcribe and purify mRNAs using the mMESSAGE mMACHINE T7 kit (Thermo Fisher Scientific).
- 5 Prepare the mRNA mixture containing 0.3 pg NDC80-FKBP, 0.3 pg NUF2, 0.45 pg H2B-mCherry, 3.5 pg hAG-Pfv-Sapphire-GS-FRB, and 0.15 pg FKBP-GS-FKBP.
- 6 Inject the mRNA mixture into GV-arrested oocytes, followed by 2 hours incubation at 37 °C.
- 7 Induce meiotic resumption by washing out IBMX.
- 8 Culture oocytes in M2 medium containing 500 nM rapamycin (Funakoshi) to induce formation of the protein-based artificial kinetochore (NDC80-GEM-Linker).