


Oct 02, 2023

Conversion of Plasmid DNA to Minicircles

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4qzozvo5/v1>



Mathew Chu¹

¹Walter and Eliza Hall Institute of Medical Research



Mathew Chu

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4qzozvo5/v1>

Protocol Citation: Mathew Chu 2023. Conversion of Plasmid DNA to Minicircles. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4qzozvo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 01, 2023

Last Modified: October 02, 2023

Protocol Integer ID: 88641

Keywords: minicircles plasmid dna, conversion of plasmid dna, plasmid dna, dna into dna minicircle, dna minicircle, plasmid, specific integration into the host genome, delivering transgene, conversion to minicircle, transgenes into cell, host genome, cloning, bacteria, genome, dna, minicircle, transgene sequence, need for special producer strain, using standard molecular biology approach, special producer strain, standard molecular biology approach, molecular biology

Abstract

Plasmid DNA is conventionally the choice of vector for cloning and delivering transgenes into cells. However, for site-specific integration into the host genome, it is desirable that only the transgene sequence is transferred without elements of the backbone.

In this protocol, we describe a novel approach to converting plasmid DNA into DNA minicircles which consist exclusively of the transgene sequence. This approach takes advantage of the scalability of plasmids while enabling in vitro conversion to minicircles using standard molecular biology approaches, circumventing the need for special producer strains of bacteria.

Image Attribution

Image is under a creative commons licence.

Materials

- NEBuilder HiFi DNA assembly cloning kit (NEB)
- Type IIS restriction endonuclease (NEB)
- T4 DNA ligase (NEB)
- rCutSmart buffer (10X) (NEB)
- ATP (10X) (NEB)
- T5 exonuclease (NEB)
- NEBuffer 4 (10X) (NEB)
- Qiaquick PCR purification kit (Qiagen)

Adapting the Plasmid Backbone for Minicircle Production

- 1 Design a cloning site for the transgene which consists of outward-facing recognition sequences (RS) for a Type IIS restriction endonuclease:

```
cloning_site = insert_left + rev_comp(RS) + stuffer + RS +  
insert_right
```

The left and right sequences should be homologous to the transgene. A stuffer sequence of approximately 50 nt should be included between the recognition sequences to enable efficient double digestion of the cloning site.

- 2 Construct the following ssDNA oligo pool (oPool) to adapt the plasmid backbone for minicircle production:

```
adaptor_oligo = plasmid_left + NNNN + cloning_site + NNNN* +  
plasmid_right
```

Here, the left and right sequences should be homologous to the plasmid. Each oligo should consist of the cloning site flanked by a different pair of complementary 4 nt overhangs (NNNN, NNNN*).

- 3 Digest the backbone at the desired site for cloning the adaptor oligos and purify by gel separation and extraction.
- 4 Clone the adaptor oligos from the oPool directly into the linearized backbone via NEB HiFi assembly using the protocol for **bridging dsDNA with ssDNA oligos**.
- 5 Inoculate transformed bacteria directly into 100 mL of liquid broth with the appropriate selection antibiotic and culture for 24 hours at 37 C in a shaking incubator.
- 6 The next day, perform a midiprep to extract the library of adapted plasmids. Verify the diversity of the new cloning site by targeted PCR and Sanger sequencing with degenerate base calling.

Cloning the Transgene into the Adapted Plasmids



- 7 Insert the transgene into the adapted plasmid using HiFi assembly at the new cloning site.
- 8 As above, inoculate transformed bacteria directly into liquid culture and extract plasmids the next day.

Synthesis of DNA Minicircles

- 9 The following procedure for minicircle DNA synthesis is adapted from Oliynyk & Church (2022).

Using 12 ug of the plasmid library, set up the following digestion-ligation reaction:

- 6 uL Type IIS restriction endonuclease
- 6 uL T4 DNA ligase
- 10 uL ATP (10X)
- 10 uL rCutSmart buffer (10X)

Make the reaction up to 100 uL volume with molecular-grade biology water.

- 10 Perform the following digestion-ligation reaction to excise the transgene from the plasmids and circularize them:

1. 37 C for 2 hours
2. 65 C for 30 minutes
3. 4 C hold

During this protocol, plasmids will be digested at the Type IIS endonuclease recognition sequences and then re-ligated at the diverse 4 nt overhangs within the plasmid library.

Concatemer formation by trans-ligation of digested transgenes will be minimized by the diversity of the overhangs in this reaction (ensuring minicircles consist of 1 transgene sequence each).

Purification of DNA Minicircles

- 11 To remove linear dsDNA from the reaction, add the following reagents directly to the mix:
 - 5 uL T5 exonuclease
 - 5 uL NEBuffer 4 (10X)

Perform this final digestion at 37 C for 1 hour.



- 12 Using the Qiaquick PCR purification kit, purify the minicircles and elute in the desired volume (at least 10 uL). The expected yield from this protocol is ~ 1 ug per 6 ug of input plasmid DNA (~ 5 kb in size).

Protocol references

Oliynyk, R.T., Church, G.M. Efficient modification and preparation of circular DNA for expression in cell culture. *Communications Biology* **5**, 1393 (2022). <https://doi.org/10.1038/s42003-022-04363-z>