

Apr 21, 2020

Version 1

## COPAS wormsorter V.1

DOI

[dx.doi.org/10.17504/protocols.io.bfc9jjz6](https://dx.doi.org/10.17504/protocols.io.bfc9jjz6)

Saul Moore<sup>1</sup>, Ida Barlow<sup>1</sup>

<sup>1</sup>Imperial College London

Behavioural Genomics



Ida Barlow

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bfc9jjz6>

**Protocol Citation:** Saul Moore, Ida Barlow 2020. COPAS wormsorter. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bfc9jjz6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 21, 2020

**Last Modified:** October 22, 2020

**Protocol Integer ID:** 35969

**Keywords:** copas wormsorter protocol, dispensing adult worm, wormsorter protocol, adult worm, worm, copa


## Abstract

Protocol for dispensing adult worms using the COPAS 500 flowpilot

## Troubleshooting

## Prepare equipment

1



COPAS wormsorter indicating key components

Turn on the compressor at the wall – it should show a pressure of 40psi after switched on

- 2 Turn on COPAS machine with switch on the left hand side
- 3 Turn on the lasers (488 laser sufficient if using unmarked animals). Add in picture of lasers.
- 4 Turn on the computer
- 5 Discard waste contents that are in the recovery cup (small shallow cup on the left-hand side of the machine)
- 6 Check that there is water in the sheath. If the water is low, fill up with MQH2O (not M9).
- 7 Make sure that the recovery cup and sample cup are securely tightened so that there are no leaks in the system

## Prepare software



- 8 Open dbgview – should always be running in the background
- 9 Open FlowPilot software and a prepared experiment with a set gate for eg Adults. :
- 9.1 File → Load Experiment
- 9.2 File → Load sample

## Clean system

- 10 Maintenance → Flush Sample
- 11 Click 'Refill Sample' – the sample cup pressure should decrease. You can see this in the software on the left hand size (include screenshot).
- 12 Unscrew sample cup and replace with falcon filled with cleaning solution (pink in colour)
- 13 Once securely replaced click 'Done refill'
- 14 Check 'Sample on' and 'mixer on' – cleaning solution should now pass through the system; allow a 2-3 ml to pass through (make sure sheath is unchecked)

### Note

Sometimes the sample cup pressure doesn't decrease and in fact increases. You can still unscrew the sample cup but if this persists there may be a blockage.

### Note

You will get a warning about contaminating the flow cell, this normal and you can click 'Yes'



15 Uncheck 'Sample on' or click Abort to stop sample flow.

16 Repeat steps 11-15 with water

## Load sample

17 Repeat steps 11-13 with sample.

18 **Turn mixer ON.** If you do not do this you may lose all your worms that have settled to the bottom of the tube!!!

19 Maintenance → Prime Flow Cell; to flush sample through the system and remove air bubbles

20 Maintenance → Flush sample

21 Check 488nm (and 568nm) laser boxes

22 Check 'Use sort gate' for stored sort gate – include screenshot of software here

23 Click 'Acquire' – sample should pass through the system and number of events per second will be shown:

- Aim for 10-20 events per second
- If too few/too many events increase/decrease 'Sample cup pressure' so that it is between 1.5-2psi
- To ensure only one event per droplet go to Setup→Coincidence, select 'Pure, no double'. This increases accuracy in the number of worms dispensed but the time to dispense may increase.

## Test dispensing

24 Click on the plate icon on the top bar



- 25 Select number of objects to sort
- 26 Select the wells you would like to fill (for testing we use a spare 60mm plate and fill wells A1, A2, B1, B2)
- 27 Select which gate to use
- 28 Apply
- 29 Place 60mm plate in front left corner of left-hand stage with A1 in the left corner.
- 30 Click 'Fill plate'
- 31 Keep an eye on the number of events per second
- 32 Ensure the 'Diverter pressure' is checked
- 33 Check under microscope that the correct number of objects were dispensed per 'well'
- 34 If too many objects, decrease sample cup pressure and repeat steps 8-11 or select Pure no double to increase accuracy.

## Fill plate

- 35 Click on the plate icon on the top bar
- 36 'Clear plate'
- 37 Select number of objects per well and click 'Apply to All' or select which wells you would like to fill.



- 38 Apply
- 39 Place 96 well plate in left-hand stage
- 40 Ensure 'Diverter pressure' is checked'; if it is not then liquid comes out of the dispenser constantly and you get flooding.
- 41 'Fill plate'
- 42 Keep an eye on the number of events per second still and monitor how much sample fluid is coming through the system

## Clean system

- 43 Repeat steps 11-15
- 44 Keep sample cup with water secured so that the system is air-tight and closed
- 45 Turn off all equipment (Computer, lasers, compressor, worm sorter).