


Jun 22, 2023

Version 1

Copying Archived LTEE Samples V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l61yykvqe/v1>



Zachary Blount¹, Jeffrey Barrick²

¹Michigan State University; ²The University of Texas at Austin

the-ltee



Jeffrey E Barrick

Michigan State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l61yykvqe/v1>

External link: <https://the-ltee.org>

Protocol Citation: Zachary Blount, Jeffrey Barrick 2023. Copying Archived LTEE Samples . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bp2l61yykvqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 10, 2022

Last Modified: June 22, 2023

Protocol Integer ID: 64296

Keywords: archived ltee sample, ltee sample, copies of freezer stock, homogeneous cell population, genetic diversity, making copy, representative copy, term evolution experiment, population samples versus clonal sample, new mutation, copy, population sample, clonal sample, fossil, freezer stock

Funders Acknowledgements:

National Science Foundation

Grant ID: DEB-1951307

Abstract

This protocol describes how to replenish the frozen "fossil" record of the *E. coli* Long-Term Evolution Experiment (LTEE) by making copies of freezer stocks. Different methods are used to copy *population samples* versus *clonal samples*. For the former, the goal is to make a representative copy that maintains the genetic diversity that is present in the population. For the latter, the goal is to make a copy of a homogeneous cell population without allowing it to evolve new mutations.

The culture media and supplies needed for this protocol are described here:

Protocol



NAME

LTEE Media Recipes

CREATED BY

Jeffrey E Barrick

Preview

Copying Archived Population Samples

3d

- 1 Use this procedure to make copies of archived population samples from the LTEE that are as close as possible to exact copies of the original samples.

Note

Strictly speaking, making an exact copy of a bacterial population is impossible. Stored samples are frozen with glycerol as a cryoprotectant (which *E. coli* cells can metabolize), and are revived in **DM1000**. Thus, a frozen population that is being revived is exposed to a resource environment that contains both glycerol and a far higher glucose concentration than is found in DM25. As these conditions are very different than those the population experienced during the LTEE prior to archiving, outgrowth during revival is an unavoidable perturbation. However, this procedure has been designed to ameliorate and reduce the effects of this perturbation so as to maintain the *representative genetic diversity* within the population sample, minimize changes in the *frequencies of different genotypes*, and avoid *further evolution caused by new mutations*.

- 1.1 Disinfect the bench with [M] 70 % (v/v) ethanol or [M] 10 % (v/v) bleach and light a Bunsen burner to create an updraft, which will reduce the possibility of contaminants falling into materials during work.
- 1.2 Prepare an autoclaved 50-mL Erlenmeyer flask capped with a 20-mL beaker for each sample that will be revived, plus a blank. Fill each flask with  9.9 mL of **DM1000**.
- 1.3 Retrieve vials of frozen stock from -80°C storage, place in ice in an insulated ice bucket, and allow the stock to thaw completely.
- 1.4 Vortex each freezer vial and pipette  120 μL of the thawed stock into one of the flasks. Wipe off the micropipettor with a Kimwipe moistened with [M] 70 % (v/v) ethanol before pipetting from each thawed stock to prevent cross-contaminating the archived samples.

**Note**

Normally  100 μL of culture is transferred each day in the LTEE to a new flask. Reviving  120 μL here accounts for the 1:5 volume of  80 % (v/v) glycerol added to the culture prior to freezing so that, at least nominally, we are preserving a normal sampling of the genetic diversity that would be present in the population of *E. coli* cells after an LTEE transfer.

- 1.5 Incubate flasks  120 rpm, 37°C  Overnight
- 1.6 Disinfect bench with  70 % (v/v) ethanol or  10 % (v/v) bleach and light a Bunsen burner to create an updraft, which will reduce the possibility of contaminants falling into materials during work.
- 1.7 Prepare another set of the same number of flasks. Fill each flask with  9.9 mL of **DM25**.
- 1.8 Prepare one test tube filled with  9.75 mL of sterile saline for each sample.
- 1.9 Dilute the grown cultures 1:40 by pipetting  250 μL into a test tube filled with saline.
- Note**
- This dilution attempts to match the initial density in the new DM25 cultures to what it would be after a normal daily transfer. It may not do so perfectly (especially for Cit⁺ populations), but it should work to reasonably preserve the representative genetic diversity within the population.
- 1.10 Thoroughly vortex each test tube three times before using a micropipettor to transfer  100 μL from it into one of the new flasks filled with **DM25**.



1.11 Incubate the flasks  120 rpm, 37°C for  24:00:00 .

1d

1.12 Add  2 mL of  80 % (v/v) glycerol to each flask. Swirl or gently vortex to mix. Pipette the culture-glycerol mix into glass vials that have been labeled with the appropriate sample information. A small glass vial (0.5 dram) can be filled with  1.25 mL each, while large glass vials (2 dram) can be filled with  6 mL each. Be careful not to overfill the vials, or they can break as the contents freeze and expand.

Copying Archived Clonal Samples

2d

2 Use this procedure to make copies of samples of clonal isolates from the LTEE.

Note

Clonal isolates were obtained by picking single colonies from agar plates after multiple rounds of purification by streaking. They should consequently have minimal genetic diversity. (Ideally, they should have no genetic diversity, but mutations occur during the growth of colonies, so low frequencies of genetic variants are unavoidable.) Unlike with the above case of copying samples of mixed populations, the principal concern with copying clonal samples is minimize the potential for *further evolution caused by new mutations*.

- 2.1 Disinfect the bench with  70 % (v/v) ethanol or  10 % (v/v) bleach and light a Bunsen burner to create an updraft, which will reduce the possibility of contaminants falling into materials during work.
- 2.2 Prepare an autoclaved 50-mL Erlenmeyer flask capped with a 20-mL beaker for each sample that will be revived, plus a blank. Fill each flask with  10 mL of **DM1000**.
- 2.3 Retrieve vials of frozen stock from -80°C storage, place in ice in an insulated ice bucket, and allow the stock to thaw completely.
- 2.4 Vortex each freezer vial and pipette  12 µL of the thawed stock into one of the flasks. Wipe off the micropipettor with a Kimwipe moistened with  70 % (v/v) ethanol before pipetting from each thawed stock to prevent cross-contaminating the archived samples.

**Note**

Alternatively, do not let the frozen stock melt. Scrape  2 μL –  10 μL of frozen culture from the freezer vial using a pipette tip attached to a pipettor and use it to inoculate the flask. The important point in this step is to inoculate with enough cells that it is unlikely that a low frequency genetic variant in the stock might be overrepresented in the inoculum.

2.5 Incubate the flasks  120 rpm, 37°C  Overnight

2.6 Add  2 mL of [M] 80 % (v/v) glycerol to each flask. Swirl or gently vortex to mix. Pipette the culture-glycerol mix into glass vials that have been labeled with the appropriate sample information. A small glass vial (0.5 dram) can be filled with  1.25 mL each, while large glass vials (2 dram) can be filled with  6 mL each. Be careful not to overfill the vials, or they can break as the contents freeze and expand.