

Nov 21, 2020

# Coral tissue and skeleton DNA extraction

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.bi9bkh2n](https://dx.doi.org/10.17504/protocols.io.bi9bkh2n)

Molly A Moynihan<sup>1</sup>

<sup>1</sup>Marine Biological Laboratory



Molly A Moynihan

Marine Biological Laboratory



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bi9bkh2n](https://dx.doi.org/10.17504/protocols.io.bi9bkh2n)

**Protocol Citation:** Molly A Moynihan 2020. Coral tissue and skeleton DNA extraction. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bi9bkh2n>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 01, 2020

**Last Modified:** November 21, 2020

**Protocol Integer ID:** 39939








## Abstract

This extraction protocol is based on the Qiagen DNeasy PowerBiofilm kit with modifications from Sunagawa et al. 2010.

Sunagawa, S., Woodley, C. M. & Medina, M. Threatened Corals Provide Underexplored Microbial Habitats. PLoS ONE 5, e9554 (2010).

## Materials

### MATERIALS

-  Proteinase K (20 mg/ml)
-  nuclease free water
-  Lysozyme
-  1000uL sterile filter tips
-  DNeasy PowerBiofilm Kit **Qiagen Catalog #24000-50**
-  20ul sterile filter tips
-  200ul sterile filter tips






Equipment		
Block heater		NAME
	--	BRAND
	--	SKU

Equipment		
24 Vortex Adapter for 2mL tubes		NAME
	-	BRAND
	-	SKU





Equipment	
Centrifuge	NAME
Benchtop Centrifuge	TYPE
Eppendorf	BRAND
5405000441	SKU
<a href="https://online-shop.eppendorf.us/US-en/Centrifugation-44533/Centrifuges-44534/Centrifuge-5425-PF-243560.html">https://online-shop.eppendorf.us/US-en/Centrifugation-44533/Centrifuges-44534/Centrifuge-5425-PF-243560.html</a>	LINK
Any benchtop centrifuge will suffice	SPECIFICATIONS



## Preparation




- 1 If needed, prepare lysozyme by dissolving  12.5 mg of lysozyme powder ( $\geq 40,000$  units/mg) per  1 mL of nuclease free water. Multiple aliquots can be prepared and stored at   $-20\text{ }^{\circ}\text{C}$ .
- 2 Clean surfaces and pipettes with 10% bleach and 70% ethanol.
- 3 Thaw samples and lysozyme (if frozen).
- 4 Warm solution solution MBL at   $65\text{ }^{\circ}\text{C}$  for  00:10:00.
- 5 Pre-label tubes.

## Extraction

- 6 Add between  200  $\mu\text{L}$  to  400  $\mu\text{L}$  of airbrushed coral tissue or  200  $\text{cm}^3$  -  500  $\text{cm}^3$  of pulverized skeleton to each bead tube.

The optimal amount of starting material may vary for each sample, depending on how much PBS was used to airbrush sample, the amount of material of the skeleton's endolithic community, and the coral species.










Too much starting material can result in low yields, particularly for coral species with thick tissue.

- 7 Add  350 mL warm MBL solution to each bead tube.
- 8 Add  100  $\mu\text{L}$  solution FB to each bead tube. Vortex briefly.
- 9 Add  10  $\mu\text{L}$  lysozyme (12.5mg/ml) (see Step 1) to the bead sample mixture.



- Incuatate at Room temperature for 00:10:00 .
- 10 Add 20  $\mu$ L proteinase-K (20mg/ml) to each sample and incubate at 65  $^{\circ}$ C for 01:00:00 .
- 11 Bead beat the sample with a Vortex Adapter for 00:15:00 .
- 12 Centrifuge the tubes at 13000 x g, Room temperature, 00:01:00 .
- 13 Transfer the supernatant to a clean 2ml collection tube.
- 14 Add 200  $\mu$ L solution IRS and vortex briefly to mix. Incubate at 4  $^{\circ}$ C for 00:05:00
- 15 Centrifuge the tubes at 13000 x g, Room temperature, 00:01:00 .
- 16 Avoiding the pellet, transfer all of the supernatant to a 2ml collection tube.
- 17 Add 900  $\mu$ L solution MR and vortex briefly.  
(Note: if solution MR has precipitated, warm at 55 $^{\circ}$ C for 5-10 minutes)
- 18 Load 650  $\mu$ L supernatant onto an MB spin column and centrifuge at 13000 x g, Room temperature, 00:01:00 . Discard the flow-through and repeat until all the supernatant has been processed.
- 19 Place the MB Spin Column into a clean 2ml Collection Tube.
- 20 Shake Solution PW. Add 650  $\mu$ L solution PW and centrifuge at 13000 x g, Room temperature, 00:01:00 .



- 21 Discard the flow-through and add  650  $\mu\text{L}$  ethanol and centrifuge at  13000 x g, Room temperature, 00:01:00 .
- 22 Discard the flow-through and centrifuge again at  13000 x g, Room temperature, 00:02:00 to ensure all ethanol is removed from the filter.
- 23 Place the MB Spin Column basket into a clean 2ml collection tube.
- 24 Add  60  $\mu\text{L}$  to  100  $\mu\text{L}$  of nuclease free water to the center of the white filter membrane, depending on expected yield and desired concentration.  
  
Make sure entire membrane is wet. Incubate at room temperature for  00:05:00 .
- 25 Centrifuge at  13000 x g, Room temperature, 00:01:00 and discard the MB Spin Column.
- 26 Proceed to DNA quantification and dilute (if needed) with nuclease-free water.  
  
Store at  -20 °C for short term storage (e.g. to be used within the same week) or  -80 °C for long term storage.