

Jun 25, 2025

Version 2

CosMx for FFPE V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl41mrolo5/v2>

Joanna Bi¹

¹Stanford University

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: Jeff.spraggins@vanderbilt.edu



Joanna Bi

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl41mrolo5/v2>

Protocol Citation: Joanna Bi 2025. CosMx for FFPE. protocols.io

<https://dx.doi.org/10.17504/protocols.io.3byl41mrolo5/v2> Version created by [Joanna Bi](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 25, 2025

Last Modified: June 25, 2025

Protocol Integer ID: 220950

Keywords: cosmx for ffpe cosmx protocol, ffpe cosmx protocol, ffpe tissue, cosmx, ffpe, tissue

Abstract

CosMx protocol for FFPE tissue

- 1 FFPE blocks were sectioned into 5 μm sections using a microtome, placed onto BOND Plus Slides (Leica, SKU S21.2113.A), and dried as instructed prior to use (MAN-10159-02, NanoString). Sections were positioned within the imageable area as defined in the manual.
- 2 Dried sections were stored at 4°C until use.
- 3 Slides were baked at 60°C for 18 hours.
- 4 The FFPE sections were then deparaffinized and subjected to antigen retrieval at 100°C using Bond Epitope Retrieval Solution 2 (Leica, AR9640).
- 5 Incubation frames were assembled on the slides, and sections were permeabilized with Proteinase K at 40°C for 30 minutes, followed by a 5-minute fiducial incubation.
- 6 Post-fixation was carried out using 10% NFB, followed by Tris-Glycine treatment and blocking with NHS-acetate.
- 7 The sections were then incubated overnight with the CosMx Human Universal Cell Characterization Panel, which consists of a 950-target probe panel and 10 controls (NanoString, cat# 121500002).
- 8 After incubation, sections were washed to remove excess probes.
- 9 Cell membranes were then stained using morphology markers including CD298, PanCK, CD45, and CD68 from the CosMx Segmentation and Supplemental Markers Kit (NanoString, cat# 121500020, 121500021, 121500022).
- 10 Prepared slides were loaded onto the CosMx SMI instrument for imaging, following the NanoString manual (MAN-10161-02-1, Software Version 1.1). After pre-scanning each slide, regions of interest (ROIs) were selected as fields of view (FOVs) for further characterization.
- 11 The acquired scan data was then uploaded to NanoString's cloud-based analysis platform, AtoMx, for downstream processing and analysis.