

Dec 16, 2022

Cost-efficient Yeast genome Flongle library

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwjb2wlmk/v1>

Yutaro Hori¹

¹The University of Tokyo



Yutaro Hori

The University of Tokyo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwjb2wlmk/v1>

Protocol Citation: Yutaro Hori 2022. Cost-efficient Yeast genome Flongle library . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e6nvwjb2wlmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 16, 2022

Last Modified: December 16, 2022

Protocol Integer ID: 74086

Keywords: efficient yeast genome flongle library, flongle library, genome, yeast, flongle, library

Abstract

A cost-efficient protocol for constructing a Flongle library.

Materials

PEG/NaCl precipitation buffer (see Rocky Mountain adventure protocol)

	A	B
	PEG 8000	9% (w/v)
	NaCl	1 M
	Tris-HCl (pH 8.0)	10 mM

PEG wash buffer

Dilute PEG/NaCl precipitation buffer with the same amount of water.



- 1 Assemble the following components to end repair and add A to the genomic sample:

10m

	A	B
	Ultra II buffer	1.75
	End repair buffer	1.75
	Ultra II enzyme	1
	End repair enzyme	1
	DNA	x (500 ng)
	MQ	44.5 - x

Incubate at 20 °C for 00:05:00 , then 65 °C for 00:05:00 .

- 2 Add 50 µL of Ampure beads (or equivalent beads) and mix well.

Wash with 75 % EtOH twice and elute with 11 µL of 10 millimolar (mM) Tris-HCl (pH 8.0).

Keep the tube on the magnetic rack and do not disturb the beads while washing with 75% EtOH.

Use 1 µL of the solution to check the concentration with Qubit. The concentration needs to be at least 30 ng/µL .

- 3 Assemble the following components for ligation:

1h 21m

	A	B
	DNA	10 uL
	TLB	4 uL
	Quick ligase	1.5 uL
	AMX	0.8 uL

Incubate at Room temperature for 00:30:00 .



Add 2.3 μL of 5 Molarity (M) NaCl and incubate at Room temperature for 00:30:00 and then cfg at max speed for 00:20:00 .

Remove sup and add PEG wash buffer, cfg at max speed for 00:01:00 .

Remove sup and add 11 μL of 10mM Tris-HCl (pH 8.0) and incubate at 4 $^{\circ}\text{C}$ for Overnight .

Check the concentration with Qubit.

4 Assemble the following priming solution:

	A	B
	FLB	2.5 uL
	FB	97.5 uL

and the library mix solution:

	A	B
	DNA	x uL (100 ng)
	SB II	10 uL
	LB II	4 uL
	MQ	6 - x uL

. Load them and start sequencing.