

Feb 14, 2024

Coupling of TMEM192 antibody to MyOne™ Epoxy Dynabeads™

 [The Journal of Clinical Investigation](#)

DOI

dx.doi.org/10.17504/protocols.io.q26g7p2ykgwz/v1

Pawel Lis¹, Daniel Saarela¹, Dario R Alessi¹

¹MRC-PPU at The University of Dundee



Francesca Tonelli

MRC-PPU at The University of Dundee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7p2ykgwz/v1>

Protocol Citation: Pawel Lis, Daniel Saarela, Dario R Alessi 2024. Coupling of TMEM192 antibody to MyOne™ Epoxy Dynabeads™. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7p2ykgwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 12, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 95179

Keywords: ASAPCRN, lysosomes, TMEM192, tmem192 antibody, rabbit monoclonal tmem192 antibody, isolation of untagged lysosome, untagged lysosome, lysosome, coupled bead, bead

Funders Acknowledgements:

Aligning Science Across Parkinson's




Grant ID: ASAP-000463

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol describes coupling of  600 µg of rabbit monoclonal TMEM192 antibody (Abcam recombinant Anti-TMEM192 antibody [EPR14330-67], BSA and Azide free, ab232600) to  20 mg of MyOne™ Epoxy Dynabeads™ (Invitrogen™, 34001D) to obtain  2 mL of final suspension. The coupled beads generated using this protocol can be used for the isolation of untagged lysosomes from cells and tissues.

Attachments



230607-

Dynabeads_cou...

130KB

Guidelines

TM

Note

This protocol can be used for coupling of any monoclonal antibody to MyOne™ Epoxy Dynabeads. Antibody to be coupled should be provided in PBS at a concentration of at least 1.2 mg/ml, and must be free of sodium azide, amine-based buffers, glycerol and protein stabilisers (BSA). Maximum binding capacity of MyOne™ Epoxy Dynabeads™ is estimated to be ~ 20 µg - 30 µg of antibody per mg of beads. In this protocol 30 µg per mg of beads is used, to ensure complete saturation of beads.

Materials

Buffers:

	A	B
	C1	0.1 M Sodium Phosphate (Na ₂ HPO ₄ :NaH ₂ PO ₄) buffer pH 7.4
	C2	3 M Ammonium Sulphate ((NH ₄) ₂ SO ₄) in 0.1 M Sodium Phosphate buffer pH 7.4
	HB	100 mM Glycine pH 11.3, 0.01% Tween-20
	LB	200 mM Glycine pH 2.8, 0.01% Tween-20
	SB	50 mM Tris-HCl (NH ₂ C(CH ₂ OH) ₃ ·HCl) pH 7.4 with 140 mM NaCl and 0.1% Tween-20
	SBS	50 mM Tris-HCl (NH ₂ C(CH ₂ OH) ₃ ·HCl) pH 7.4 with 140 mM NaCl, 0.1% Tween-20 and 0.2% NaN ₃

Note

All buffers should be stored at 4 °C (up to 1 week) or at -20 °C (long-term storage) and must be brought to Room temperature before being used for coupling.

Troubleshooting



Protocol

1d 0h 11m 15s

- 1 Before opening the vial containing dried magnetic beads, equilibrate to Room temperature .
- 2 Thaw the antibody On ice and keep On ice until it is needed in step 8.
- 3 Calculate the volume of antibody needed, so that 600 μg is used – this volume should be \leq 500 μL .
- 4 Weigh 20 mg beads directly into a fresh low-binding 1.5 ml Eppendorf tube.
- 5 Resuspend beads in 1000 μL of sterile Milli-Q water, vortex for 00:00:15 , sonicate in a water bath sonicator for 00:05:00 .
- 6 Place vial on a magnetic rack for 00:01:00 , remove water using a pipette.
- 7 Repeat steps 5 and 6. After sonication there should be no bead aggregates visible.
- 8 Add the required volume of antibody to the vial containing washed beads.
- 9 Add buffer C1 up to total volume of 500 μL (C1= 500 - antibody volume). Vortex to resuspend the beads.
- 10 Add 500 μL of buffer C2 and vortex.



- 11 Incubate in a Thermomixer at 37 °C for 16–24 hours (typically 20:00:00) at 1500 rpm (make sure the beads do not settle). 20h
- 12 Place on a magnetic rack for 00:01:00 , remove liquid using a pipette. 1m
- 13 Resuspend beads in 1000 µL of buffer HB, vortex.
- 14 Place on a magnetic rack for 00:01:00 , remove liquid using a pipette. 1m
- 15 Resuspend beads in 1000 µL of buffer LB, vortex.
- 16 Place on a magnetic rack for 00:01:00 , remove liquid using a pipette. 1m
- 17 Resuspend beads in 1000 µL of buffer SB, vortex.
- 18 Repeat steps 16 and 17: 1m
- Place on a magnetic rack for 00:01:00 , remove liquid using a pipette.
 - Resuspend beads in 1000 µL of buffer SB, vortex.
- 19 Place on a magnetic rack for 00:01:00 , remove liquid using a pipette. 1m
- 20 Resuspend beads in 1000 µL of buffer SB, vortex.
- 21 Incubate in a shaker at Room temperature for 1500 rpm, 00:15:00
- 22 Place on a magnetic rack for 00:01:00 , remove liquid using a pipette. 1m



- 23 Resuspend beads in $\text{1000 } \mu\text{L}$ of buffer SBS. At this stage beads are at 20 mg/ml and should be stored in the fridge. Beads can be further diluted with buffer SBS to 10 mg/ml , which is the usual working concentration.