



May 24, 2023

## CoxII degradation assay to assess mitophagy

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygxzxr4v8j/v1>

Thanh Ngoc Nguyen<sup>1</sup>, nguyen.tha<sup>1</sup>

<sup>1</sup>Laboratory of Michael Lazarou, Walter and Eliza Hall Institute of Medical Research, Parkville, Victoria, Australia

Thanh Ngoc Nguyen: [nguyen.tha@wehi.edu.au](mailto:nguyen.tha@wehi.edu.au)

ASAP workspace



nguyen.tha

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.kxygxzxr4v8j/v1>

**Protocol Citation:** Thanh Ngoc Nguyen, nguyen.tha 2023. CoxII degradation assay to assess mitophagy. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kxygxzxr4v8j/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 07, 2022

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 66232

**Keywords:** CoxII degradation, Mitophagy, HeLa cells, ASAPCRN, analysing coxii degradation, coxii degradation via western blotting, mitophagy, coxii degradation, mitophagy this protocol

**Funders Acknowledgements:**

Aligning Science Across Parkinson's

Grant ID: ASAP-000350

## Abstract

This protocol details the procedure to assess mitophagy by analysing COXII degradation via Western blotting.

## Attachments



[ideybr7hf.docx](#)

18KB

## Materials

### Buffers and reagents:

### Growth media:

#### DMEM

|  | A                      | B       |
|--|------------------------|---------|
|  | FBS                    | 10%     |
|  | Glucose                | 4.5 g/l |
|  | GlutaMAX <sup>TM</sup> | 1x      |
|  | MEM NEAA               | 1x      |
|  | HEPES                  | 25 mM   |

- 45% D-( )-Glucose Merck MilliporeSigma (Sigma-Aldrich) Catalog #G8769
- GlutaMAX<sup>TM</sup> Supplement Thermo Fisher Scientific Catalog #35050061
- MEAA (MEM Non-Essential Amino Acids) Gibco - Thermo Fisher Scientific Catalog #11140050
- Antimycin A from Streptomyces sp. Merck MilliporeSigma (Sigma-Aldrich) Catalog #A8674 (made up in 100% Ethanol to 20 mg/ml), Oligomycin (Calbiochem, 495455; made up in DMSO to 10 mg/ml) and Q-VD-OPh MedChemExpress Catalog #HY-12305 (made up in DMSO to 10 mM).
- **Lysis buffer:** 1x LDS + 0.1 M DTT (diluted from 4x LDS (NP007; ThermoFisher); can be aliquoted and stored at -20 or -80 oC)
- 4-12% Bis-Tris NuPAGE gels (ThermoFisher)
- NuPAGETM Antioxidant (NP0005, ThermoFisher; use 0.5 ml/ 200 ml of gel running buffer)
- 20x NuPAGETM MOPS SDS running buffer (NP001, ThermoFisher)
- 20x NuPAGE transfer buffer (NP00061, ThermoFisher)
- PVDF destain: 40% Methanol, 7% Acetic Acid.
- 1x PBS
- 1x PBS/0.1% Tween20 (PBS/Tween)
- **Blocking buffer:** 5% skim milk in PBS/Tween (make fresh)
- ACTIN (Cell Signaling, 4967S), COXII (Abcam, ab110258), Parkin (Santa Cruz, sc-32282)
- Amersham ECL Prime Western Blotting Detection Reagent (RPN2232)



4h 52m 10s

## Procedures:

- 1 Seed the hela cells the day before the treatment day in 6 well plates.
  - 1.1 Each well contain  2 mL of growth media.
  - 1.2 Seed 350,000 cells for penta KO expressing BFP-Parkin and GFP-OPTN or -NDP52.
  - 1.3 Adjust the number of cells of other cell lines. So that the next day they are all in similar confluency with penta KO expressing BFP-Parkin and GFP-OPTN or -NDP52.
- 2 The next day, make sure the seeded cells are spreading out (not concentrated in the middle of the well because this can affect the results).
- 3 Aspirate off the old media and treat each well with  2 mL of growth media containing  4 micromolar ( $\mu\text{M}$ ) Antimycin A,  10 micromolar ( $\mu\text{M}$ ) Oligomycin and  10 micromolar ( $\mu\text{M}$ ) QVD for desired period.

### Note

**Note:** Make sure all drugs are vortexed well, mix the media well after adding each drug.

- 4 After the treatment, harvest the cells on ice by scraping.
  - 4.1 Pre-chill eppies and 1x PBS on ice.

### Note

**Note:** I normally put all the plates that need harvesting into a fridge and harvest one by one on ice.
  - 4.2 Aspirate the media thoroughly from the wells.



4.3 Wash the wells with  1 mL of cold 1x PBS.



Note

**Note:** Make sure swirl around after adding the PBS to wash the cells properly.

4.4 Aspirate off the PBS and add another  1 mL of cold 1x PBS.



4.5 Use a plastic cell scraper to scrape all the cells off the wells.

Note

**Note:** I use one scraper for each well. You can wash and reuse them again.

4.6 Transfer the cells-containing PBS to eppies.

4.7 Centrifuge the eppies at  3000 rpm for  00:02:00 at  4 °C .

2m



4.8 Aspirate off PBS.

4.9 Quickly centrifuge for  00:00:10 to spin down the residual PBS.

10s



4.10 Aspirate off all the PBS.

5 Lyse the cell pellets in lysis buffer and boil the samples at  99 °C with shaking for  00:07:00 .

7m



## Note

**Note:** I use the plastic clips to make sure that the lids won't pop during heating.

- 6 Let the samples cool down and spin at max speed ( Room temperature ) for 00:01:00 .

1m

- 7 Estimate the protein concentration by nanodrop

## Note

**Note:** Make sure the concentrations do not exceed 6 mg/mL . If they do, dilute with lysis buffer and reheat them for a couple of minutes at 99 °C with shaking.

- 8 Aliquot 25 µg of each sample into a new eppie and add 1x LDS to make up to 15 µL .



- 9 Set up the gel tank with MOPs buffer.

## Note

The inside chamber should be filled with 1x MOPs supplemented with antioxidants. The outside chamber doesn't need antioxidants.

- 10 Wash each well with a glass syringe.



- 11 Load markers and samples into the wells and run at 100V for 00:10:00 and 190V for 00:55:00 .

1h 5m

- 12 Subject the gels to wet transfer onto PVDF membrane using cold NuPAGE transfer buffer containing 20% Methanol for 01:00:00 at Room temperature .

1h

- 13 After transfer, incubate the PVDF membrane with PVDF destain buffer on a shaker at Room temperature for 00:02:00 .

2m





14 Wash three times with PBS/Tween ( 5 min each wash).



14.1 Wash the PVDF membrane with PBS/Tween for  00:05:00 (1/3).

5m



14.2 Wash the PVDF membrane with PBS/Tween for  00:05:00 (2/3).

5m



14.3 Wash the PVDF membrane with PBS/Tween for  00:05:00 (3/3).

5m



15 Block with blocking buffer for  00:15:00 .

15m

16 Remove blocking buffer.

17 Rinse twice with PBS/Tween and wash twice with PBS/Tween and once with 1x PBS (5 min for each wash).



17.1 Rinse the blocking buffer with PBS/Tween and wash with PBS/Tween for  00:05:00 (1/2).

5m



17.2 Rinse the blocking buffer with PBS/Tween and wash with PBS/Tween for  00:05:00 (2/2).

5m



17.3 Rinse the blocking buffer with 1x PBS and wash with 1x PBS for  00:05:00 .

5m



18 Cut the PVDF membrane and put appropriate parts into different antibodies.

#### Note

**Note:** In this case, it's ACTIN (1/5000), COXII (1/1000) and Parkin (1/1000) antibodies made up in 3% BSA in PBS/Tween.



19 Incubate on a 4 °C shaker Overnight .



#### Note

**Note:** To make sure we don't lose antibodies, I wet the tubs with PBS/Tween before putting in the antibodies.

20 The next day, recycle the antibodies back to their tubes.



21 Wash the blots three times with PBS/Tween (5 min for each wash).



21.1 Wash the blot with PBS/Tween for 00:05:00 (1/3).

5m



21.2 Wash the blot with PBS/Tween for 00:05:00 (2/3).

5m



21.3 Wash the blot with PBS/Tween for 00:05:00 (3/3).

5m



22 Incubate with appropriate HRP-conjugated secondary antibodies made up in blocking buffer for 01:00:00 .

1h



23 Wash the blots twice with PBS/Tween, once with 1x PBS (5 min for each wash).



23.1 Wash the blot with PBS/Tween for 00:05:00 (1/2).

5m



23.2 Wash the blot with PBS/Tween for 00:05:00 (2/2).

5m



23.3 Wash the blot with 1x PBS for 00:05:00 .

5m





24     Develop the blots with ECL prime.

25