



Aug 26, 2023

CRISPR/Cas9 genome editing

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzxb38gx1/v1>

wusj¹, schekman¹

¹Department of Molecular and Cell Biology, Howard Hughes Medical Institute, University of California, Berkeley, Berkeley, United States



Nancy C. Hernandez Villegas

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzxb38gx1/v1>

Protocol Citation: wusj , schekman 2023. CRISPR/Cas9 genome editing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzxb38gx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: August 03, 2023

Last Modified: August 26, 2023

Protocol Integer ID: 85898

Keywords: cas9 genome editing this protocol, cas9 genome editing, generation of dnajc5 ko, dnajc5 ko, using crispr, crispr, cas9

Funders Acknowledgements:

Randy Schekman

Grant ID: Howard Hughes Medical Institute

Abstract

This protocol describes the generation of DNAJC5 KO using CRISPR/Cas9 edition



CRISPR/Cas9 genome editing

- 1 gRNA targeting exon 4 of DNAJC5 (CACC GGAG GCCG CAGA AGAC AAAC A) was inserted into a pX330-based plasmid expressing Venus fluorescent protein
- 2 HEK293T cells were transfected with pX330-pX330-Venus-DNAJC5-Exon4-gRNA by Lipofectamine 2000 followed the commercial protocol by Thermo Fisher.
- 3 After 48 hr, we diluted the cells and single colonies were isolated from 96-well plate, expanded, and determined for DNAJC5 KO by immunoblot.