

Jul 23, 2025

CRISPR-Cas9 genome editing in *Corallochytrium limacisporum*

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygqxq6q4v8j/v1>

Patricia S. Ara¹, Elena Casacuberta¹, Claudio Scazzocchio^{2,3}, Iñaki uiz-Trillo^{1,4}, Sebastián Najle^{1,5}

¹Institut de Biologia Evolutiva (CSIC-Universitat Pompeu Fabra), Barcelona, Spain;

²Department of Life Sciences, Imperial College London, London, UK;

³Université Paris-Saclay, CEA, CNRS, Institute for Integrative Biology of the Cell (I2BC), Gif-sur-Yvette, France;

⁴ICREA, Pg. Lluís Companys 23, 08010 Barcelona, Spain;

⁵Centre for Genomic Regulation (CRG), Barcelona Institute for Science and Technology, Barcelona, Spain.



Patricia S Ara

University of Galway

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygqxq6q4v8j/v1>

Protocol Citation: Patricia S. Ara, Elena Casacuberta, Claudio Scazzocchio, Iñaki uiz-Trillo, Sebastián Najle 2025. CRISPR-Cas9 genome editing in *Corallochytrium limacisporum*. protocols.io <https://dx.doi.org/10.17504/protocols.io.kxygqxq6q4v8j/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 09, 2025

Last Modified: July 23, 2025

Protocol Integer ID: 219832

Keywords: Genome editing, *Corallochytrium limacisporum*, Unicellular holozoans, CRISPR-Cas9, cas9 genome editing in *corallochytrium limacisporum*, genome editing in *corallochytrium limacisporum*, cas9 genome editing, genome editing, cas9 rnp complex, crispr, *corallochytrium limacisporum*, unicellular marine eukaryote, resistance to rapamycin,

Funders Acknowledgements:

MICIU/ AEI /10.13039/501100011033

Grant ID: PID2020-120609GB-I00

MICIU/AEI /10.13039/501100011033

Grant ID: PID2023-153273NB-I00

Moore Foundation

Grant ID: and MMI Experimental Model Systems grant #4973.01

MICIU/AEI/10.13039/501100011033

Grant ID: PRE2018-085438

Abstract

This protocol walks you through how to achieve genome editing in *Corallochytrium limacisporum*, a unicellular marine eukaryote which is closely related to animals. Here we describe how to achieve resistance to rapamycin or carboxin. In brief, the cells are Nucleofected, to introduce a sgRNA-Cas9 RNP complex. Next, the cells are selected on agar plates and isolated for genotyping.

Materials

- Marine Broth (BD Difco #2216)
- Marine Agar (Scharlau #01-291)
- T-25 flasks (Corning #430639)
- 90 mm cell culture dishes (Falcon #353002)
- 12-well plate (Corning #3513)
- 96-well plate (Nunc #167008)
- Cell scraper (Corning #3010)
- 5 mm glass beads (Sigma-Aldrich #18406)
- P3 Primary Cell 4D-Nucleofector X Kit (Lonza #V4XP-3032)
- Amaxa® 4D-Nucleofector®
- EnGen Spy Cas9 NLS (New England Biolabs, #M0646M)
- PBS (Sigma #P5368-10PAK)
- Rapamycin Sigma-Aldrich #37094
- Carboxin Sigma-Aldrich #45371
- Ca-Mg Free Sea Water: NaCl 460 mM, NaHCO₃ 2.15 mM, KCl 10.7 mM, Na₂SO₄ 7 mM, Tris-HCl pH 8.0 10 mM



Before start

- *C. limacisporum* is prone to contamination as its growth medium is very rich, always proceed under sterile conditions (laminar flow cabinet) to avoid contamination.
- Unless otherwise stated, *C. limacisporum* is grown in Marine Broth, its usual growth medium.
- All reagents and their references can be found in the Materials section.



Seed cells for transfection

20m

- 1 Prepare a T-25 flask of *C. limacisporum* so that it is growing exponentially on the day of transfection. Usually,  100-200 μL from a confluent culture in a 25 cm^2 flask with  5 mL of Marine Broth works well.
- 2 Incubate the cells at  23 $^{\circ}\text{C}$  Overnight in an incubator without agitation.

20m



Transfection

40m

- 3 Before starting, prepare the following reagents and keep them on ice:
 - 1X Filter-sterilized PBS
 - Lonza nucleofector strip (part of P3 Primary Cell 4D-Nucleofector X Kit)
 - Lonza P3 Buffer (part of P3 Primary Cell 4D-Nucleofector X Kit)
 - Repair DNA template (if needed)

Note


Ensure that the supplement has been added to the Lonza P3 buffer. Supplemented Lonza P3 Buffer expires 3 months after preparation. To avoid wasting P3 buffer, aliquot the components to prepare 1/3 of the volume:

-  225 μL P3 solution
-  50 μL supplement

- 4 Prepare RNP complex (sgRNA+Cas9) mix. One mix should be prepared separately for each sgRNA to be used. When adding the final component, gently swirl it with the pipette tip and incubate at RT for at least 30 min.

30m

RNP mix:

-  10 μL Lonza P3 Buffer
-  2.5 μL sgRNA at  60 micromolar (μM)
-  2 μL EnGen Spy Cas9 NLS

Proceed with the next steps during the  00:30:00 incubation.

**Note**

For co-transfection experiments, prepare the two mixes independently:

Mix 1:

- 5 μ L of P3 buffer
- 2 μ L EnGen Spy Cas9 NLS
- 2.5 μ L sgRNA1

Mix 2:

- 5 μ L of P3 buffer
- 2 μ L EnGen Spy Cas9 NLS
- 2.5 μ L sgRNA2

- 5 Scrape and homogenize the cells, then determine the concentration by counting a 1:100 dilution of the homogenized cell suspension using a Neubauer chamber.

- 6 Pellet 1.5×10^6 cells per reaction at  2000 x g, 00:05:00

5m

Note

Cells for multiple transfections can be pelleted together in the same tube.

Note

During centrifugation: Prepare a 12-well plate with 1 mL of Marine Broth for every transfection

- 7 Remove supernatant from pelleted cells and wash the pellet with  100 μ L **ice-cold** sterile PBS.

- 8 Spin the cells again  2000 x g, 00:05:00

5m

**Note**

During centrifugation: Prepare the Lonza nucleofector.
Turn on Lonza nucleofector, select the vessel (strip) and select the wells that are going to be used.

- Solution: **Primary cell P3**
- Pulse Code: **EN138**

- 9 Remove supernatant from the washed cells and resuspend in  6.5 μ L ice-cold P3 buffer for every transfection that will be done.

Note

Don't leave the cells in P3 buffer for longer than necessary.

- 10 Prepare mix for nucleofection:

 14.5 μ L RNP mix  [go to step #4](#)

 6.5 μ L cells resuspended in P3  [go to step #9](#)

 2 μ L repair DNA template (if needed)

Note

For co-transfection with two guides:

- 9.5 μ L RNP mix 1
- 9.5 μ L RNP mix 2
- 6.5 μ L cells resuspended in P3
- 2 μ L repair template (if needed)

Note

When doing HDR, different length repair templates can be used at 100 uM. No repair template is necessary for NHEJ.

- 11 Add transfection mixes to the Lonza strip, place the strip in the machine and apply **EN138** pulse



- 12 Quickly remove strip from Lonza and add  80 μ L Marine Broth to allow the cells to recover
- 13 Transfer transfected cells to 12-well plate with 1 mL Marine Broth
- 14 Grow cells at  23 $^{\circ}$ C  Overnight for selection on rapamycin, or  48:00:00 for selection on carboxin

Plating onto Marine Agar

- 15 Melt  20 mL of Marine Agar (Scharlau #01-291) for every plate that you are going to prepare.

Note

If you are using rapamycin always use fresh plates and wait until it is quite cool to the touch but not solidifying yet, as rapamycin is very sensitive to heat. Carboxin plates can be prepared in advance and stored at 4 $^{\circ}$ C for up to a week.

- 16 Add selection drug to the agar and add  20 mL agar per dish. Allow them to solidify until it is dry.

Note

Rapamycin selection concentration: 125 ng/mL
Carboxin selection concentration: 30 μ g/mL

- 17 Scrape the transfected wells and add  100-200 μ L of the culture onto the agar.

Note

If pelleting a larger volume, spin the cells down first and resuspend the pellet in 100 μ L of Marine Broth first.



- 18 Add 3-4 glass beads (Sigma-Aldrich #18406) and shake the plate until the agar is completely dry. Gently remove the glass beads by pouring them from the plate. Seal the plate with parafilm.
- 19 Incubate the plates at  23 °C for 10-14 days.

Note

Monitor the plates to ensure contaminants don't appear, these will usually appear faster than *C. limacisporum* colonies. *C. limacisporum* will usually start to be detectable after around 5-7 days and appear as very small, transparent colonies. Eventually, they will become round and pink.

Colony picking and streaking

- 20 Prepare Marine Agar plates containing the selection drug as above  [go to step #15](#) and label them with the approximate number of colonies that have appeared, or however many you want to isolate and screen.
- 21 Gently pick the colonies by touching them with a sterile pipette tip and gently streaking onto a fresh agar plate labelled with clone identifiers. Seal the plate with parafilm.
- 22 Incubate the plates at  23 °C allowing the streaks to grow for 2-3 days until the streak is visible.

Transfer clones to liquid medium

- 23 Once the streak is visible, pick a part of the streak using a sterile tip and place it in  100 µL Marine Broth with selecting drug in a 96-well plate.
- 24 Allow the cells to grow by incubating the plates at  23 °C for 2-3 days before genotyping.

Genotyping

- 25 Pipette up and down to resuspend cells in the 96-well plate.



- 26 Transfer  50 μL of the culture to a tube. Fresh Marine Broth can be added to the well to maintain the culture if needed, otherwise they can also be retrieved from the streaked agar plates.

Note

50 μL of culture in a 0.2 mL tube or PCR strip is enough to genotype with, but larger volumes can be used by scaling up the amount of Ca-Mg free sea water used to wash the cells in step 28.

- 27 Pellet the cells at  2500 x g, 00:03:00 . If using PCR strips, use a Minispin with strip adaptor for 3 minutes.

Note

From this step forward, sterile environment is no longer required.

- 28 Remove the supernatant and wash the cells with  100 μL Ca-Mg free sea water. Repeat the wash once.

Note**Ca-Mg free sea water:**

NaCl 460 mM, NaHCO_3 2.15 mM, KCl 10.7 mM, Na_2SO_4 7 mM, Tris-HCl pH 8.0 10 mM

- 29 After the second wash, heat the cell suspension for  00:10:00 at  95 $^{\circ}\text{C}$ in a Thermoblock or PCR machine.

- 30 Spin down cell remains at  2500 x g, 00:05:00 and use  1-2 μL of the supernatant as template in genotyping PCR.



Note

Lysed cells can be stored at 4°C before genotyping for up to a week.