


Jul 05, 2019

CRISPR-Cas9 Genome Editing in Human Primary T Cells

DOI

dx.doi.org/10.17504/protocols.io.4vdgw26



Stemcell Technologies¹

¹STEMCELL Technologies

STEMCELL Technologies

Tech. support email: techsupport@stemcell.com



Stemcell Technologies

STEMCELL Technologies

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4vdgw26

External link: https://www.stemcell.com/genome-editing-of-human-primary-t-cells-using-the-arcitect-crispr-cas9-system.html?utm_source=protocolsio&utm_medium=referral

Protocol Citation: Stemcell Technologies 2019. CRISPR-Cas9 Genome Editing in Human Primary T Cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4vdgw26>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 27, 2019

Last Modified: July 05, 2019

Protocol Integer ID: 25221

Keywords: gene editing, CRISPR, Cas9, electroporation, cell therapy, immunotherapy, regenerative medicine, CAR T cell therapy, genetics, molecular biology, transfection

Abstract

CRISPR-Cas9, an RNA-guided genome editing technology, is revolutionizing cell biology due to the ease and efficiency by which it enables genetic manipulation of mammalian cells. Through targeted modification of specific genes or regulatory regions, researchers can now rapidly generate precise genetic models to study normal and diseased cell physiology. Beyond genetic manipulation for research purposes, CRISPR-Cas9 genome editing also holds great potential for therapeutic applications, including immunotherapy and regenerative medicine. This protocol provides instructions for CRISPR-Cas9 genome editing of human primary T cells, including the isolation and activation of primary human T cells, preparation of a CRISPR-Cas9 RNP complex, delivery of an RNP complex into activated primary T cells using electroporation, and methods to evaluate genome editing efficiency.

For a discussion of alternative strategies and expected outcomes, and protocol modifications, please see:

Technical Bulletin: Genome Editing of Human Primary T Cells.

Guidelines

ArciTect™

CRISPR-Cas9 Genome Editing Precision for Cell Biologists

The **ArciTect™** family of products for CRISPR-Cas9 genome editing provides you with a rapid, flexible and precise system to modify the genome as you see fit. The guide RNA (gRNA) complex, consisting of crRNA and tracrRNA, can be customized using the **online order tool**. A protocol for the efficient genome editing of human embryonic stem (hES) cells or human induced pluripotent stem (hiPS) cells has been optimized using the ArciTect™ product family.

Reduce Off-Target Effects

The ArciTect™ product family is a ribonucleoprotein (RNP)-based Cas9 genome editing system. Unlike previous CRISPR technologies which utilize plasmid or mRNA-based systems, the ArciTect™ system shows timely degradation of the RNP complex to minimize cleavage of off-target regions.

Why use ArciTect™?

- **CUSTOMIZABLE.** Design crRNA to target your sequence of interest.
- **FLEXIBLE.** Multiple variations of Cas9 to suit your specific genome editing needs.
- **RAPID.** No need for transcription and translation.
- **REDUCED OFF-TARGET EFFECTS.** Timely degradation of the RNP complex to minimize potential off-target cutting.



Materials

MATERIALS


- ☒ Gentamicin **Gibco - Thermo Fisher Scientific Catalog #15750037**
- ☒ Human Peripheral Blood Leuko Pak, Fresh **STEMCELL Technologies Inc. Catalog #70500**
- ☒ EasySep™ Human T Cell Isolation Kit **STEMCELL Technologies Inc. Catalog #17951**
- ☒ ImmunoCult™-XF T Cell Expansion Medium **STEMCELL Technologies Inc. Catalog #10981**
- ☒ L-Glutamine **STEMCELL Technologies Inc. Catalog #07100**
- ☒ Human Recombinant IL-2 (CHO-expressed) **STEMCELL Technologies Inc. Catalog #78036**
- ☒ ImmunoCult™ Human CD3/CD28 T Cell Activator **STEMCELL Technologies Inc. Catalog #10971**
- ☒ ArciTect™ crRNA **STEMCELL Technologies Inc.**
- ☒ ArciTect™ tracrRNA Kit **STEMCELL Technologies Inc.**
- ☒ Nuclease-Free Water **STEMCELL Technologies Inc. Catalog #79001**
- ☒ Costar® 24-Well Flat-Bottom Plate Tissue Culture-Treated **STEMCELL Technologies Inc. Catalog #38017**
- ☒ Falcon® Conical Tubes 15 mL **STEMCELL Technologies Inc. Catalog #38009**
- ☒ Costar® Microcentrifuge Tubes **STEMCELL Technologies Inc. Catalog #38038**
- ☒ Neon™ Transfection System 100 µL Kit **Thermo Fisher Scientific Catalog #MPK10025**
- ☒ ArciTect™ Cas9 Nuclease **STEMCELL Technologies Inc. Catalog #76002**
- ☒ ArciTect™ Cas9-eGFP Nuclease **STEMCELL Technologies Inc. Catalog #76006**
- ☒ ImmunoCult™ Human CD3/CD28/CD2 T Cell Activator **STEMCELL Technologies Inc. Catalog #10970**
- ☒ ArciTect™ High-Fidelity DNA Polymerase Kit **STEMCELL Technologies Inc. Catalog #76026**
- ☒ ArciTect™ T7 Endonuclease I Kit **STEMCELL Technologies Inc. Catalog #76021**

Optional: Source frozen primary T cells in Step 1.

- ☒ Human Peripheral Blood Pan-T Cells, Frozen **STEMCELL Technologies Inc. Catalog #70024**

Before start

This protocol describes the following workflow for T cell genome editing.



The guide RNA (gRNA) complex, consisting of crRNA and tracrRNA, can be designed and customized using our [online order tool](#).



T Cell Isolation and Activation

- 1 Isolate human T cells from peripheral blood using EasySep™ Human T Cell Isolation Kit. Refer to the Product Information Sheet for details.
- 2 Count cells and adjust to 1×10^6 cells/mL in ImmunoCult™-XF T Cell Expansion Medium supplemented with:
 - 2 mM L-Glutamine
 - 50 µg/mL Gentamicin
 - 10 ng/mL Human Recombinant IL-2
- 3 Activate human T cells by adding 25 µL/mL of ImmunoCult™ Human CD3/CD28 T Cell Activator. Incubate cell suspension at 37°C and 5% CO₂ for 72 hours.

3d

Note

Optional: Assess T cell activation by binding the activation marker CD25 with Anti-Human CD25 Antibody, Clone BC96 and performing flow cytometry.

Preparation of Stock Solutions

- 4 Briefly centrifuge the vials of ArciTect™ crRNA and ArciTect™ tracrRNA before opening.
- 5 Depending on the size of the ArciTect™ crRNA and ArciTect™ tracrRNA vials being used, add nuclease-free water to each vial to give a final concentration of 200 µM, as indicated in Table 1.

ArciTect™ crRNA OR ArciTect™ tracrRNA Vial Size	Catalog #	Volume of Nuclease-Free Water (µL)
2 nmol	76010 / 76016	10
10 nmol	76011 / 7601	50

	7	
20 nmol	76012 /7601 8	100

Table 1. Resuspension volumes for 200 μM * ArciTect™ crRNA or ArciTect™ tracrRNA.*200 μM is equal to 200 pmol/ μL

- Mix thoroughly. If not used immediately, aliquot and store at -80°C for up to 1 month. After thawing the aliquots, use immediately. Do not re-freeze.

Preparation of ArciTect™ CRISPR-Cas9 RNP Complex

- Prepare a 60 μM guide RNA solution by combining components in an RNase-free microcentrifuge tube as indicated below in Table 2. Volumes are sufficient for a single electroporation reaction (10 μL /electroporation reaction); adjust as required. Mix thoroughly.

Com pone nt	Volu me per Elect ropor ation (μL)
Nucle ase- free water	2
ArciT ect™ Anne aling Buffle r (5X)	2
200 μM ArciT ect™ crRN A	3
200 μM ArciT ect™ tracr RNA	3
Total volu me	10

Table 2. Preparation of 60 μ M Guide RNA solution per electroporation reaction.**Note**

Note: ArciTect™ Annealing Buffer (5X) and ArciTect™ tracrRNA are both included in the ArciTect™ tracrRNA kit.

- 8 In a thermocycler or heating block, incubate guide RNA solution at 95°C for 5 minutes followed by 60°C for 1 minute. Cool to room temperature and place on ice.

6m

Note

Note: If not used immediately, store guide RNA solution at -80°C for up to 1 month.

- 9 To prepare the ribonucleoprotein (RNP) complex, combine 60 μ M guide RNA (prepared in step 8) and **ArciTect™ Cas9 Nuclease** (or **ArciTect™ Cas9-eGFP**) in an RNase-free microcentrifuge tube as indicated below in Table 3. Volumes are sufficient for a single electroporation; adjust as required for the desired number of wells.

Com pone nt	Volu me per Elect ropor ation (μ L)
Resu spens ion Buffle r T	6.8
4 μ g/ μ L ArciT ect™ Cas9 Nucle ase*	7.2
60 μ M guide RNA	6
Total volu me	10

Table 3. Preparation of RNP Complex Mixture.

**Note**


***Note:** If using 3 µg/µL ArciTect™ Cas9-eGFP Nuclease, add 11.4 µL of Cas9-eGFP to 6 µL of 60 µM guide RNA and 2.6 µL of Resuspension Buffer T for a total volume of 20 µL. Note that 1 µg/µL ArciTect™ Cas9 should not be used for Neon® electroporation, as the total Resuspension buffer volume should represent ~90% of the total reaction volume.

Note: Resuspension Buffer T is included in the Neon™ Transfection System 100 µL Kit.

- 10 Mix the RNP Complex Mixture by gently pipetting up and down 2 times; avoid creating air bubbles or foam in the tube.
- 11 Incubate the RNP Complex Mixture at room temperature (15 - 25°C) for 10 - 15 minutes.

15m

Electroporation of RNP Complex into T Cells

- 12 For each electroporation condition (including positive and negative controls), prepare 2 mL of ImmunoCult™-XF T Cell Expansion Medium supplemented with:
 - 2 mM L-Glutamine
 - 50 µg/mL Gentamicin
 - 10 ng/mL Human Recombinant IL-2.
- 13 For each condition, add 2 ml of supplemented medium (prepared in step 12) to 1 well of a 6-well plate and place in a 37°C incubator. Store the remainder of the supplemented medium at 2 - 8°C.
- 14 Transfer 1.2×10^6 cells from the activated T cell suspension (prepared in section A) to a 15 mL conical tube. Centrifuge at  300 x g , for 5 minutes at room temperature.
- 15 Aspirate medium and resuspend cells in 100 µL of Resuspension Buffer T.
- 16 Using individual RNase-free microcentrifuge tubes for each condition, combine RNP Complex Mixture (prepared in steps 7 - 11) with activated T cells as indicated in Table 4.

Com pone nt	Volu me Per Elect ropor

	ation (μ L)
RNP Comp lex	20
Activ ated T cell suspe nsion	100
Total Volu me	120

Table 4. Mixture of RNP Complex and T Cells for Electroporation

- 17 Mix gently by pipetting up and down carefully 2 times; avoid creating air bubbles or foam in the tube.

Note

Note: If air bubbles are present in the tip when the cells are electroporated, cell viability and transfection efficiency will be significantly reduced.

Note: Refer to the manufacturer's instructions for electroporation. Electroporation conditions may require optimization for different cell types.

- 18 Use a 100 μ L Neon® pipette tip to draw up 100 μ L of the mixture and place into the electroporation chamber containing 3 mL of Electrolytic Buffer E2. Electroporate mixture according to conditions indicated in Table 5.

Electroporati on Parameter	
Electr ical poten tial	1400 V
Pulse width	30 millis econd s
Numb er of pulse s	1

Table 5. Recommended Electroporation Conditions for Human T Cells Using a Neon® Electroporation Device

19 Immediately add the electroporated cells to 1 well of the 6-well plate prepared in step 2. Incubate at 37°C and 5% CO₂ for 2 - 3 hours.

3h

20 Count cells and adjust cell density to 2.5×10^5 viable cells/ mL by adding IL-2-supplemented ImmunoCult™-XF T Cell Expansion Medium (prepared in step 12).

Note

Optional: If using Cas9-eGFP nuclease, assess electroporation efficiency 12 - 24 hours after electroporation by flow cytometry.

21 Incubate at 37°C and 5% CO₂ for 48 - 72 hours for genome editing to occur. Harvest cells for assessment of genome editing efficiency, genomic DNA can be amplified by PCR using primers flanking the target region and **ArciTect™ High-Fidelity DNA Polymerase Kit**, followed by sequencing of PCR products. Alternatively, **ArciTect™ T7 Endonuclease Kit** can be used to assess editing efficiency (% INDEL formation) following PCR amplification.

3d

Expected result




Figure 1. TRAC gRNA design and evaluation. (A) Schematic of the TRAC locus with crRNA sequences 1 - 4 aligned to the TRAC first exon. (B) Representative TCRαβ/CD3 flow cytometry results at 48 hours post RNP complex electroporation (containing the corresponding crRNA1 - 4) into activated T cells. (C) Genome editing (cleavage) efficiency was assessed at 48 hours post electroporation using the ArciTect™ T7 Endonuclease I Kit. U: Uncut; C: Cut; (D) gRNA target sequences for TRAC knockout.




Figure 2. T Cells express the CD25 activation marker after treatment with ImmunoCult™ Human T Cell Activators. Isolated human T cells were activated with either ImmunoCult™ Human CD3/ CD28 or CD3/CD28/CD2 T Cell Activator for 2 or 3 days. Activation status was assessed by CD25 flow cytometry. Each data point per condition represents an individual donor; n = 4 - 8 donors. Error bars represent standard error of the mean.




Figure 3. High efficiency TRAC knockout across activation conditions and dynamics. (A) TRAC knockout efficiency in human T cells activated with either ImmunoCult™ Human CD3/CD28 or CD3/CD28/CD2 T Cell Activator for 2 or 3 days was assessed by binding the TCRαβ and CD3 receptors with antibodies followed by flow cytometry analysis. Each data point per condition represents an individual donor; n = 4 - 8 donors. Error bars represent standard error of the mean. (B) Genome editing (cleavage) efficiency was assessed at 48 hours post electroporation in human T cells activated with ImmunoCult™ Human CD3/CD28 T Cell Activator for 3 days using the ArciTect™ T7 Endonuclease I Kit. Mock electroporated: - RNP; RNP electroporated: + RNP. (C - D) Representative dot plots of TCRαβ and CD3 flow cytometry analysis from (C) mock electroporated and (D) RNP electroporated human T cells activated with ImmunoCult™ Human CD3/CD28 T Cell Activator for 3 days. (E) Representative dot plot of CD4 and CD8 flow cytometry analysis of human T cells activated with ImmunoCult™ Human CD3/CD28 T Cell Activator for 3 days.

