



Nov 12, 2023

Version 1

CRISPR knock-in validation V.1

DOI

dx.doi.org/10.17504/protocols.io.n2bvj3qb5lk5/v1

Leonardo A Parra-Rivas¹

¹University of California, San Diego



Leonardo A Parra-Rivas

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj3qb5lk5/v1>

Protocol Citation: Leonardo A Parra-Rivas 2023. CRISPR knock-in validation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj3qb5lk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 12, 2023

Last Modified: November 12, 2023

Protocol Integer ID: 90815

Keywords: validation crispr, crispr knock, crispr

Abstract

CRISPR knock-in validation

Troubleshooting

1

Cultured hippocampal neurons tagged with oScarlet were prepared with Lucigen QuickExtract DNA extraction solution (Biosearch Technologies, cat# QE09050). Briefly, neurons were lysed by adding 50 ul of QuickExtract solution to each sample.

2

The samples were mixed by pipetting and incubated at 68 °C for 15 min followed by 95 °C for 10 min in a thermocycler before being stored at -20 °C for downstream analysis.

3

Two different PCR reactions were performed to amplify DNA products of ~500 bp corresponding to the 5' and 3' integration junctions. PCR#1 used the α -syn forward primer: TGTGCTTTCTCTTCCCTCTCTG and the reverse oScarlet primer CCGTCCTCGAAGTTCATCAC, whereas PCR#2 used the α -syn forward primer: ATAACACTTCGTGCAGCACC and the reverse oScarlet primer ACAGGATGTCCCAGGAGAAG.

4

PCR products were extracted from the agarose gel using the Monarch DNA gel extraction kit (New England Biolabs, cat# T1020S), and samples were submitted for sequencing for analysis at MCLAB.