

Mar 28, 2025

# 🌐 Cross-Linking and Strong Cation Exchange (SCX) Fractionation

📖 [Nature](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge5q2og47/v1>

Miguel A. Gonzalez-Lozano<sup>1</sup>, Harper JW<sup>2</sup>

<sup>1</sup>Harvard Medical School; <sup>2</sup>harvard university



Harper JW

harvard university

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.261ge5q2og47/v1>

**External link:** <https://doi.org/10.1038/s41586-025-09059-y>

**Protocol Citation:** Miguel A. Gonzalez-Lozano, Harper JW 2025. Cross-Linking and Strong Cation Exchange (SCX) Fractionation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ge5q2og47/v1>

**Manuscript citation:**

Gonzalez-Lozano, M.A., Schmid, E.W., Miguel Whelan, E. *et al.* EndoMAP.v1 charts the structural landscape of human early endosome complexes. *Nature* **643**, 252–261 (2025). <https://doi.org/10.1038/s41586-025-09059-y>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 18, 2024

**Last Modified:** March 28, 2025

**Protocol Integer ID:** 110306

**Keywords:** ASAPCRN, linked peptide, positive early endosome, early endosome, purified organelle sample, crosslinker, lys cross, strong cation exchange, scx

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's**

Grant ID: 000282, 025160

## Abstract

This is a method for performing Lys-Lys cross-linking on a purified organelle sample, in this case EEA1-positive early endosomes. The crosslinker used is DSSO ((Bis(2,5-dioxopyrrolidin-1-yl) 3,3'-sulfinyldipropionate)). This cross-linker has a mass spec cleavable sulfoxide, which facilitates identification of individual cross-linked peptides.

## Materials

disuccinimidyl sulfoxide (DSSO; Thermo Scientific, A33545)

2-chloroacetamide (Sigma-Aldrich, C0267)

trypsin (Promega, V511C)

Lys-C (Wako Chemicals, 129-02541)

Sep-Pak C8 50mg Cartridge (Waters, WAT054965)

EPPS (3-[4-(2-Hydroxyethyl)-1-piperazine]propanesulfonic acid) (Thermo Scientific, J61296AE)

dithiothreitol (DTT)

Empore SPE Disks C18 to generate stage tips(Sigma Millipore, 66883-U)

PolyLC PolySulfoethyl A column (3 µm particle size, 2.1 mm inner diameter, and 100 mm length) (#102SE03--)

Acetonitrile

Formic acid

Lys-C (Wako Chemicals, 129-02541)

## Troubleshooting



## Organelle preparations

- 1 Samples for cross-linking (e.g. individual organelles) are isolated using dedicated protocols such as <https://doi.org/10.17504/protocols.io.ewov14pjyvr2/v2> for isolation of EEA1-positive early endosomes

## Protein cross-linking

- 2 Resuspend freshly prepared purified organelle pellets in KPBS (25 mM KCl, 100 mM potassium phosphate, 150 mM NaCl, pH 7.2) to a protein concentration of  $\sim 1\mu\text{g}/\mu\text{L}$ .
- 3 Cross-link the proteins by adding 1 mM DSSO (Bis(2,5-dioxopyrrolidin-1-yl) 3,3'-sulfinyldipropionate) and incubate at room temperature for 40 minutes.
- 4 Quench the cross-linking reaction by adding 50 mM Tris buffer (pH 7.5) and incubate at room temperature for 30 minutes.
- 4.1 Alternatively, other cross-linker reagents could be use. For DHSO/DMTMM, 8 mM DHSO (3,3'-sulfinyldi(propanehydrazide)) and 16mM DMTMM (4-(4,6-dimethoxy-1,3,5-triazin-2-yl)-4-methylmorpholinium chloride) at 37°C for 90 min.  
Proceed directly with protein precipitation. Samples can be digested using S-trap sample preparation following mini-spin column digestion protocol as provided by the manufacturer (Protifi, C02-mini-80).

## Protein digestion

- 5 Denature the cross-linked samples by adding 8 M urea.
- 6 Reduce the proteins by adding 5 mM dithiothreitol (DTT) and incubate for 30 minutes at 37°C.
- 7 Alkylate the reduced proteins with 40 mM chloroacetamide for 30 minutes at room temperature.
- 8 Digest the cross-linked proteins with Lys-C enzyme at a ratio of 1:75 (enzyme) and incubate at 37°C overnight.



- 9 Dilute the urea concentration of the sample to 2 M using 50 mM EPPS buffer (pH 8.0).
- 10 Further digest the proteins by adding trypsin at a ratio of 1:100 (enzyme) and incubate at 37°C for 6 hours.
- 11 Desalt the peptides using Sep-Pak solid-phase extraction column.
- 12 Dry the desalted peptides in a SpeedVac concentrator. Time will depend on the vacuum strength.

### Peptide separation by SCX

- 13 Fractionate the dried peptides by strong cation exchange (SCX) chromatography. A 70-minute linear gradient at flow rate 0.18 mL/min was used as follows:
  - 0% to 8% of mobile phase (0.5 M NaCl in 20% ACN, 0.05% formic acid) over 14 minutes.
  - Increase to 20% at 28 minutes.
  - Increase to 40% at 48 minutes.
  - Increase to 90% at 68 minutes.
- 14 Collect fractions every 30 seconds starting at 35 minutes for 10 minutes, and then every minute thereafter.
- 15 Dry the collected fractions in a SpeedVac concentrator.
- 16 Desalt the dried fractions using StageTip.
- 17 Reconstitute fractions in 5% acetonitrile (ACN) and 5% formic acid.
- 18 Analyze the reconstituted fractions by liquid chromatography-tandem mass spectrometry (LC-MS/MS).