

Aug 12, 2022

Crude Membrane Fractionation of Cultured Cells

DOI

dx.doi.org/10.17504/protocols.io.yxmvmnb99g3p/v1

Asad Malik¹, Dario R Alessi¹, Suzanne R Pfeffer²

¹Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Sciences, University of Dundee, Dow Street, Dundee DD1 5EH, UK.;

²Stanford University School of Medicine, Stanford, CA 94305-5307

asap



Dario R Alessi

Dundee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmnb99g3p/v1>

External link: <https://doi.org/10.1042/BCJ20220308>

Protocol Citation: Asad Malik, Dario R Alessi, Suzanne R Pfeffer 2022. Crude Membrane Fractionation of Cultured Cells. protocols.io <https://dx.doi.org/10.17504/protocols.io.yxmvmnb99g3p/v1>

Manuscript citation:

Malik AU, Karapetsas A, Nirujogi RS, Chatterjee D, Phung TK, Wightman M, Gourlay R, Morrice N, Mathea S, Knapp S, Alessi DR, PKC isoforms activate LRRK1 kinase by phosphorylating conserved residues (Ser1064, Ser1074 and Thr1075) within the COR GTPase domain. Biochemical Journal 479(18). doi: [10.1042/BCJ20220308](https://doi.org/10.1042/BCJ20220308)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 16, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 64712

Keywords: Cultured Cells, Crude Membrane Fractionation, Immunoblotting analysis, ASAPCRN, crude membrane fractionation of cultured cell, insoluble fraction from cultured cell, crude membrane fractionation, cultured cell, crude cellular extract, membrane versus cytosolic distribution, signaling pathway, protocol for fractionating, cytosol, cell, components from Irrk1, containing insoluble fraction, membrane, enriched fraction, cytosolic distribution, fractionating

Abstract

We present here a protocol for fractionating crude cellular extracts to prepare membrane and cytosol-enriched fractions and a nuclei-containing insoluble fraction from cultured cells. We deploy this protocol for determining the membrane versus cytosolic distribution of components from LRRK1 and LRRK2 signaling pathways.

Note

We recommend analysing the products of this fractionation scheme by quantitative immunoblotting (as described [indx.doi.org/10.17504/protocols.io.6qpvr68e3vmk/v1](https://doi.org/10.17504/protocols.io.6qpvr68e3vmk/v1)).

Note

This protocol was adapted from <https://doi.org/10.15252/emj.201798099>

Attachments



[466-976.docx](#)

781KB

Materials

MATERIALS

Reagents:

Buffer A:

[M] 10 millimolar (mM) HEPES  7.4 and cOmplete™ EDTA-free Protease Inhibitor Cocktail

Note

Added fresh before use,  cOmplete™, Mini, EDTA-free (Protease Inhibitor) **Roche Catalog ##11836170001**

Buffer B:

	A	B
	HEPES pH 7.4	250 mM
	Sodium Chloride	750 mM
	Magnesium Chloride	25 mM
	DTT	2.5 mM
	GDP	500 nM
	Sodium Fluoride	250 mM
	Sodium Pyrophosphate	25 mM
	Microcystin-LR (Enzo Life Sciences, ALX-350-012)	5 µg/ml
	cOmplete™ EDTA-free Protease Inhibitor Cocktail (added fresh before use, Roche, 11836170001)	

Note

This buffer is prepared at a 5X stock to achieve a final concentration of 1X in the resuspension buffer (4 X Buffer **A** + 1 X Buffer **B**).

Buffer C:

A	B
HEPES pH 7.4	50 mM
Sodium Chloride	150 mM
Magnesium Chloride	5 mM
DTT	0.5 mM
GDP	100 nM
Sodium Fluoride	50 mM
Sodium Pyrophosphate	5 mM
Microcystin-LR (Enzo Life Sciences, ALX-350-012)	1 µg/ml
Triton X-100	1% (v/v)
cOmpleteTMEDTA-free Protease Inhibitor Cocktail (added fresh before use, Roche, 11836170001)	

-  Gibco™ PBS pH 7.4 **Thermo Fisher Scientific Catalog #10728775**
-  Pierce & Warriner; Coomassie Plus (Bradford) Assay Kit **Thermo Fisher Catalog #23236** or equivalent).
- 4X Loading buffer:

 NUPAGE LDS sample buffer (4x) **Thermo Fisher Scientific Catalog #NP0007** or 4X SDS loading buffer:

A	B
Tris-HCl, pH6.8	250mM
SDS	8% (w/v)
Glycerol	40% (v/v)
Bromophenol blue	0.02% (w/v)

Reagents and antibodies:

⊗ Rubber tipped scraper **Merck MilliporeSigma (Sigma-Aldrich) Catalog #CLS3008**

⊗ Anti-Rab7 antibody Mouse monoclonal **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R8779**

⊗ Recombinant Anti-RAB7 (phospho S72) antibody [MJF-R38-1] (ab302494) **Abcam Catalog #ab302494**

⊗ Recombinant Anti-PKC alpha antibody [Y124] (ab32376) **Abcam Catalog #ab32376**

⊗ Recombinant Anti-Sodium Potassium ATPase antibody [EP1845Y] - Plasma Membrane Loading Control **Abcam Catalog #ab76020**

⊗ Anti-Rab7 antibody Mouse monoclonal **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R8779**

⊗ α -Tubulin Antibody **Cell Signaling Technology Catalog #2144**

⊗ Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) Antibody **Cell Signaling Technology Catalog #9101**

Equipment:

Equipment

Corning® cell lifter

NAME

Corning® cell lifter

BRAND

CLS3008

SKU

<https://www.sigmaaldrich.com/IN/en/product/sigma/cls3008>

LINK

Blade L 19 mm, handle L 180 mm, sterile, case of 100

SPECIFICATIONS

or equivalent

Equipment

Eppendorf™ 5810R Centrifuge

NAME

Centrifuge

TYPE

Eppendorf

BRAND

02-262-8187

SKU

<https://www.fishersci.com/shop/products/eppendorf-5810r-centrifuge-rotor-packages-16/022628187>^{LINK}

or equivalent

Equipment

Eppendorf® microcentrifuge

NAME

Centrifuge

TYPE

Eppendorf®

BRAND

5417

SKU

https://www.sigmaaldrich.com/IN/en/product/sigma/z366013?gclid=CjwKCAjwtlavBhBkEiwAsr7-cwVnwbWrBxOUW73qAEiKOBi00UeWVtRsRmiqBWQxzgNG-e_2iC6UNRoCZKAQAvD_BwE

LI
N
K

or equivalent

Equipment

Luer Slip 1ml IV Syringes (Medicina IVS01)

NAME

Medicina Luer Slip IV Syringes can be used with any standard, filtered or safety needles.

TYPE

Medicina

BRAND

IVS01

SKU

<https://medicina.co.uk/luer-slip-iv-syringes/>

LINK

or equivalent

Equipment

25G Luer Needle (Terumo™ NN-2525R)

NAME

Terumo Hypodermic Needles

TYPE

Terumo

BRAND

NN-2525R

SKU

<https://www.terumotmp.com/products/hypodermics/terumo-hypodermic-needles.html>

LINK

or equivalent

Equipment

Thick-walled Polycarbonate Tubes (Beckman Coulter 343775)	NAME
Thick-walled Polycarbonate Tubes	TYPE
Beckman Coulter	BRAND
343775	SKU
https://www.mybeckman.in/supplies/tubes-and-bottles/open-top-tubes/343775	LINK

or equivalent

- Ultracentrifuge (Beckman Coulter Optima TLX, or equivalent)
- Ultracentrifuge rotor (Beckman Coulter TLS.55 or equivalent)
- Plate reader for Protein quantification (BioTek Epoch, or equivalent)

Troubleshooting



Crude Membrane Fractionation

1



Note

The optimal quantity of cultured cells to use to achieve an ideal yield will vary dependent on cell type. As a guideline, we use 1 x 15cm dish of HEK293 cells per replicate seeded at 1.8×10^7 cells per dish.

Pour off media from the culture dish and aspirate completely by holding plate on edge. Wash cells twice with  5 mL of ice-cold PBS.

2 Immediately transfer the dishes to ice--this is best accomplished using wet paper towel-covered steel blocks resting  On ice .

3 Add  5 mL of ice-cold PBS and scrape the cells from the dish using a cell lifter (Sigma-Aldrich CLS3008, rubber tipped scraper, or equivalent) to ensure good yield; collect in a 15 ml tube. 

4 Pellet intact cells by centrifugation at  100 x g for  00:05:00 at  4 °C and aspirate supernatant. 

5m

5 Resuspend cells in  400 μ L of Buffer **A** by gentle pipetting. 

5.1 Transfer to an 1.5ml Eppendorf tube and incubate  On ice for  00:15:00 . 

15m

Note

Note that this is a hypotonic solution and will swell the cells;  00:05:00 . is likely sufficient at this stage. 

6 Add  100 μ L of cold Buffer **B** to the cell suspension. 

7 Using a 25-gauge needle attached to a 1 ml syringe, break the cells by passing the cell suspension through the needle 25 times.

**Note**

Breakage can be monitored by transferring a few microliters of the homogenate to a glass slide, covering with a coverslip and visualizing using a low power light microscope used to visualize cultured cells; as few as 6-10 passages may be sufficient. Broken cells will lose their reflective character and small particles of cell components will be readily detected.

- 8 Centrifuge the cell suspension at  1000 x g for  00:05:00 at  4 °C and collect the supernatant in a new 1.5ml Eppendorf tube.

5m

**Note**

The pellet here will contain the nuclei and other cell debris. This can be analysed by lysing in  500 µL Buffer **C**. The supernatant represents the post-nuclear supernatant.

- 9 Load the post-nuclear supernatant into thick-walled polycarbonate tubes, appropriate for ultracentrifugation in a table top ultracentrifuge. Ultracentrifuge at  150000 x g for  00:20:00 at  4 °C .

20m

**Note**

The membrane pellet will form at the bottom of the tube.

- 10 Transfer the cytosolic fraction (supernatant) to a fresh Eppendorf tube  On ice .



- 11 Wash the membrane fraction pellet with  500 µL PBS thrice to remove any potential cytosolic contaminants.

**Note**

This may not be necessary if aspiration is complete.

- 12 Resuspend membrane pellet using  500 µL of Buffer **C** using a pipet and incubate  On ice for  00:05:00  00:20:00 to allow detergent solubilization of

25m

membrane proteins.



- 13 Centrifuge membrane protein solution at 1000 x g for 00:05:00 at 4 °C to separate solubilized membrane proteins (supernatant) from insoluble membrane proteins (pellet).

5m



- 14 Determine the protein concentration of cell lysates by Bradford assay according to the manufacturer's instructions, performing measurements in triplicate.

Note

Ensure the concentration of the samples is in the linear range for the Bradford assay. If it isn't, prepare appropriate dilutions in water of each lysate. Generally, protein concentrations of near confluent cells lysed as described above should result in protein concentrations of at least 2 µg/µL .

- 15 4×SDS–PAGE sample buffer is added to samples containing 5 µg of membrane protein or an equivalent volume of cytosolic protein, and heated at 37 °C for 00:10:00 .

10m



Analysis of fractionation products by quantitative immunoblotting analysis

1d 16h 30m

- 16 The reaction products can be analysed by quantitative immunoblotting analysis (as described in dx.doi.org/10.17504/protocols.io.6qpvr68e3vmk/v1).

A	B	C	D	E
Antibody Target	Company	Cat. number	Host species	Dilution
pS72 Rab7A	Abcam Inc.	ab302494	Rabbit	1:1000
Rab7A (Total)	Sigma	R8779	Mouse	1:2000
LRRK1 (total) (C-terminus)	MRC-PPU Reagents and Services, University of Dundee	S405C	Sheep	1 µg/ml
Tubulin	Cell Signaling Technologies	2144	Mouse	1:5,000

	A	B	C	D	E
	pT202/Y204 ERK1/2	Cell Signaling Technologies	9101	Rabbit	1:1000
	PKCα	Abcam Inc.	ab32376	Mouse	1:1000
	Na-K ATPase	Abcam Inc.	ab76020	Rabbit	1:10,000

17

1d




Figure 1: Crude membrane fractionation of HEK293 Flp-in T-REx/GFP-LRRK1 WT cells following phorbol ester stimulation.

HEK293 Flp-in T-REx/GFP-LRRK1 WT cells were induced to express GFP-LRRK1 wild type by treatment with [M] 1 mg/mL doxycycline for ⌚ 24:00:00 .

18 Serum starve the cells for ⌚ 16:00:00 and then treated ± Phorbol myristic acid (PMA) ([M] 100 ng/ml) for ⌚ 00:30:00 .

16h 30m

19 Following this, Perform the fractionation as described here and samples were subjected to immunoblot analysis with the indicated antibodies; the membranes were visualized using the Odyssey CLx scan Western Blot imaging system.





Note

Adapted from <https://doi.org/10.1101/2022.06.09.495448>.