

May 29, 2024

Cryo-EM sample preparation for RCKW:DARPin complex

 Forked from [Preparation of LRRK2 RCKW cryo-EM grids](#)

 [The Journal of Biological Chemistry](#)

DOI

dx.doi.org/10.17504/protocols.io.bp2l6224kgqe/v1

Marta Sanz Murillo¹

¹University of California, San Diego



Marta Sanz Murillo

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6224kgqe/v1>

Protocol Citation: Marta Sanz Murillo 2024. Cryo-EM sample preparation for RCKW:DARPin complex. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l6224kgqe/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: May 29, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 100824

Keywords: cryo-EM, LRRK2, structural biology, ASAPCRN, em sample preparation for rckw, em grids for rckw, darpin complex, em sample preparation, rckw, cryo, lab protocol, em grid, leschziner, sample, preparation, lab

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000519

Abstract

This is Leschziner's Lab protocol for making cryo-EM grids for RCKW:DARPin complex.

Materials

LRRK2 Buffer:

- [M] 20 millimolar (mM) HEPES pH 7.4
- [M] 150 millimolar (mM) NaCl
- [M] 0.5 millimolar (mM) TCEP
- [M] 5 % volume Glycerol
- [M] 2.5 millimolar (mM) MgCl₂
- [M] 20 micromolar (μM) GDP

Note: please change salt as needed to maintain final salt of 150 mM NaCl

Troubleshooting



Safety warnings

! For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).
Take proper precautions while freezing grids.


Before start

Decide which protein concentration to use, and create the proper LRRK2 buffers in order to obtain the right salt concentration (150 mM NaCl).



Freezing Grids

20s

- 1 Plasma clean grids.
We used UltrAuFoil Holey Gold 1.2/1.3 300 mesh grids and plasma cleaned them in the Solarus II (Gatan) using the QuantiFoil Au preset.
- 2 Dilute samples to desired concentration in the **LRRK2 buffer**. Make sure final salt is at 150 mM NaCl.
For best results, make  8 μL samples, good for freezing 2 grids. This is to minimize time spent outside of storage buffer, reducing aggregation.
- 3 Apply protein to grids and plunge freeze (3-4 μL)
We used a Vitrobot (FEI) to blot away excess sample and plunge freeze
- 4 Store grids in liquid nitrogen until ready for imaging.